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T: +44-1462-444800

F: +44-1462-444801

Caution: | United States federal law restricts this device to sale by or on the order of a physician.

MicroMaxx and SonoSite are registered trademarks or trademarks of SonoSite, Inc.

Non-SonoSite product names may be trademarks or registered trademarks of their respective owners.

Protected by one or more of the following U.S. patents: 6364839. Other patents pending.

P09042-01 12/2007

Copyright 2007 by SonoSite, Inc.

All rights reserved. Printed in the USA.

MicroMaxx Ultrasound System

User Guide Supplement 3

Introduction	1
Imaging	1
Measurements and Calculations	2
Safety	3
Troubleshooting and Maintenance	10

Introduction

This user guide supplement provides information on the following additions to the MicroMaxx[®] ultrasound system:

- C8e/8-5 MHz (C8e) transducer
- Vascular exam type on C11e/8-5 MHz (C11e) transducer

See the ultrasound system user guide for additional information on safety, on intended uses, and on preparation, use, and maintenance of the ultrasound system.

Imaging

Transducer, Exam Type, and Imaging Mode

The following table describes the transducer's exam type and imaging mode.

Table 1: Transducer, Exam Type, and Imaging Mode

Transducer	Exam Type	Imaging Mode								
		2D M Mode	THI	2D MB	2D S	CPD	Color	PW	TDI PW	CW
C8e	Gyn	X	—	X	X	X	X	X	—	—
	OB	X	—	X	X	X	X	X	—	—
	Prostate	X	—	X	X	X	X	X	—	—
C11e	Vascular	X	—	X	X	X	X	X	—	—

Measurements and Calculations

Calculations

The following tables show the calculations available by exam type.

Table 2: Calculations for C8e Transducer

Exam Type	Calculations
Gyn	Gyn Follicle
OB	OB
Prostate	Volume

Table 3: Calculations for C11e Transducer

Exam Type	Calculations
Vascular	Percent Reduction Vascular Volume Flow

Safety

Guidelines for Reducing MI and TI

The following are general guidelines for reducing MI or TI. If multiple parameters are given, then the best results may be achieved by minimizing these parameters simultaneously. In some modes changing these parameters will not affect MI or TI. Changes to other parameters may also result in MI and TI reductions. Please note the 'MI' or 'TI' read out on the right side of the LCD screen.

Table 4: MI

Transducer	Depth
C8e	↑
C11e	↑

↓ Decrease or lower setting of parameter to reduce MI or TI.
 ↑ Raise or increase setting of parameter to reduce MI or TI.

Table 5: TI (TIS, TIC, TIB)

Transducer	Color Power Doppler Settings			PW Settings
	Box Depth	PRF	Depth	
C8e	↑	↓	↑	—
C11e	↑	↓	↑	↓ (Depth)

↓ Decrease or lower setting of parameter to reduce MI or TI.
 ↑ Raise or increase setting of parameter to reduce MI or TI.
 — Data are not applicable.

Transducer Surface Temperature Rise

Table 6 lists the measured surface temperature rise from ambient ($23^{\circ}\text{C} \pm 3^{\circ}\text{C}$) of transducers used on the ultrasound system. The temperatures were measured in accordance with EN 60601-2-37 section 42 where controls and settings were positioned to give maximum temperatures.

Test 1: The transducer surface temperature test on tissue mimicking material (TMM) is based on the following standard: 42.3(a) 1, Test Method B (IEC 60601-2-37, Amendment 1). The limit is a 10°C rise from ambient, as measured on the TMM, for the C11e and a 6° rise from ambient for the C8e.

Test 2: The transducer surface temperature test in air is based on the following standard: 42.3(a) 2 (IEC 60601-2-37, Amendment 1). The limit is a 27°C rise from ambient.

Table 6: Transducer Surface Temperature Rise IEC 60601-2-37

Test	C8e (Internal Use)	C11e (External Use)
1	5.7°C	9.2°C
2	14°C	19.7°C

Output Display

Table 7: Modes in Which a Thermal or Mechanical Index Is ≥ 1.0

Transducer	Index	2D/ M Mode	CPD/ Color	PW Doppler
C8e	MI	Yes	Yes	Yes
	TIC, TIB, or TIS	No	Yes	Yes
C11e	MI	No	No	No
	TIC, TIB, or TIS	No	No	Yes

Acoustic Output Table

Table 8: Table 2 Transducer Model C8e

Operating Mode: 2D

Index Label		M.I.	TIs		TIb	Tlc	
			Scan	Non-scan			Non-scan
				$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global Maximum Index Value		1.2	(a)	—	—	(b)	
Associated Acoustic Parameter	$p_{r,3}$ (MPa)	2.66					
	W_0 (mW)		#	—	—	#	
	min of $[W_{.3}(Z_1), I_{TA.3}(Z_1)]$ (mW)				—		
	Z_1 (cm)				—		
	Z_{bp} (cm)				—		
	Z_{sp} (cm)	1.2				—	
	$d_{eq}(Z_{sp})$ (cm)					—	
	F_c (MHz)	4.92	#	—	—	—	#
	Dim of A_{aprt}	X (cm)		#	—	—	—
Y (cm)			#	—	—	—	#
Other Information	PD (μ sec)	0.413					
	PRF (Hz)	7368					
	$p_r@PII_{max}$ (MPa)	3.147					
	$d_{eq}@PII_{max}$ (cm)					—	
	Focal Length	FL_x (cm)		#	—	—	#
		FL_y (cm)		#	—	—	#
$I_{PA,3}@MI_{max}$ (W/cm ²)	256.3						
Operating Control Conditions	Control 1: Exam Type		Any				
	Control 2: Optimization		Pen				
	Control 3: Depth		4.2				
	Control 4: MB		Off				

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

Table 9: Table 3 Transducer Model C8e

Operating Mode: Color/CPD

Index Label		M.I.	TIs		TIb	TIc	
			Scan	Non-scan			Non-scan
				$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global Maximum Index Value		1.5	1.1	—	—	(b)	
Associated Acoustic Parameter	$P_{r,3}$ (MPa)	3.26					
	W_0 (mW)		51.10	—	—	#	
	min of $[W_{.3}(Z_1), I_{TA,3}(Z_1)]$ (mW)				—		
	z_1 (cm)				—		
	Z_{bp} (cm)				—		
	Z_{sp} (cm)	1.1				—	
	$d_{eq}(Z_{sp})$ (cm)					—	
	F_c (MHz)	4.79	4.81	—	—	—	#
	Dim of A_{aprt} X (cm)		0.32	—	—	—	#
Y (cm)		0.4	—	—	—	#	
Other Information	PD (μ sec)	0.766					
	PRF (Hz)	2439					
	$p_r@PII_{max}$ (MPa)	3.92					
	$d_{eq}@PII_{max}$ (cm)					—	
	Focal Length FL_x (cm)		1.24	—	—		#
	FL_y (cm)		2.5	—	—		#
$I_{PA,3}@MI_{max}$ (W/cm^2)		543.2					
Operating Control Conditions	Control 1: Exam Type	Any	Any				
	Control 2: Color Opt	Low	Med				
	Control 3: Depth	4.2	2.5				
	Control 4: PRF	≤ 839	2500				
	Control 5: Color Box Position/Size	Top/Tall and Wide	Bottom/Def width				

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

Table 10: Table 4 Transducer Model C8e

Operating Mode: M-Mode

Index Label		M.I.	TIs			TIb	TIc	
			Scan	Non-scan		Non-scan		
				$A_{aprt} \leq 1$	$A_{aprt} > 1$			
Global Maximum Index Value		1.4	—	—	—	(a)	(a) & (b)	
Associated Acoustic Parameter	$p_{r,3}$ (MPa)	3.09						
	W_0 (mW)		—	—		#	#	
	min of $[W_{.3}(Z_1), I_{TA.3}(Z_1)]$ (mW)				—			
	z_1 (cm)				—			
	Z_{bp} (cm)				—			
	Z_{sp} (cm)	1.1				#		
	$d_{eq}(Z_{sp})$ (cm)					#		
	F_c (MHz)	5.05	—	—	—	#	#	
	Dim of A_{aprt}	X (cm)		—	—	—	#	#
Y (cm)			—	—	—	#	#	
Other Information	PD (μ sec)	0.425						
	PRF (Hz)	800						
	$p_r@PII_{max}$ (MPa)	3.75						
	$d_{eq}@PII_{max}$ (cm)					#		
	Focal Length	FL_x (cm)		—	—	—		#
		FL_y (cm)		—	—	—		#
$I_{PA.3}@MI_{max}$ (W/cm^2)		432.5						
Operating Control Conditions	Control 1: Exam Type		Pro and Gyn					
	Control 2: Optimization		Pen					
	Control 3: Depth		4.2					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

Table 11: Table 5 Transducer Model C8e

Operating Mode: Pulse Wave Doppler

Index Label		M.I.	TIs		TIb	TIc		
			Scan	Non-scan			Non-scan	
				$A_{aprt} \leq 1$	$A_{aprt} > 1$			
Global Maximum Index Value		1.3	—	(a)	—	2.0	(b)	
Associated Acoustic Parameter	$p_{r,3}$ (MPa)	2.12						
	W_0 (mW)		—	#		33.80	#	
	min of $[W_{.3}(Z_1), I_{TA.3}(Z_1)]$ (mW)				—			
	z_1 (cm)				—			
	Z_{bp} (cm)				—			
	Z_{sp} (cm)	1.0				1.0		
	$d_{eq}(Z_{sp})$ (cm)					0.272		
	F_c (MHz)	4.79	—	#	—	4.80	#	
	Dim of A_{aprt}	X (cm)		—	#	—	1.12	#
Y (cm)			—	#	—	0.4	#	
Other Information	PD (μ sec)	1.13						
	PRF (Hz)	1008						
	$p_r@PII_{max}$ (MPa)	2.505						
	$d_{eq}@PII_{max}$ (cm)					0.241		
	Focal Length	FL_x (cm)		—	#	—		#
		FL_y (cm)		—	#	—		#
	$I_{PA.3}@MI_{max}$ (W/cm^2)	185.2						
Operating Control Conditions	Control 1: Exam Type	Any				Any		
	Control 2: Sample Volume	1 mm				3 mm		
	Control 3: PRF	1008 Hz				10417 Hz		
	Control 4: Sample Volume Position	Zone 4				Zone 4		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

Table 12: Table 1 Transducer Model C11e

Operating Mode: Pulse Wave Doppler

Index Label		M.I.	TIs		TIb	TIC		
			Scan	Non-scan			Non-scan	
				$A_{aprt} \leq 1$	$A_{aprt} > 1$			
Global Maximum Index Value		(a)	—	(a)	—	1.8	(b)	
Associated Acoustic Parameter	$p_{r,3}$ (MPa)	#						
	W_0 (mW)		—	#		26.29	#	
	min of $[W_{.3}(Z_1), I_{TA.3}(Z_1)]$ (mW)				—			
	z_1 (cm)				—			
	Z_{bp} (cm)				—			
	Z_{sp} (cm)	#				1.1		
	$d_{eq}(Z_{sp})$ (cm)					0.236		
	F_c (MHz)	#	—	#	—	4.36	#	
	Dim of A_{aprt}	X (cm)		—	#	—	0.28	#
Y (cm)			—	#	—	0.5	#	
Other Information	PD (μ sec)	#						
	PRF (Hz)	#						
	$p_r@PII_{max}$ (MPa)	#						
	$d_{eq}@PII_{max}$ (cm)					0.226		
	Focal Length	FL_x (cm)		—	#	—		#
		FL_y (cm)		—	#	—		#
$I_{PA.3}@MI_{max}$ (W/cm^2)	#							
Operating Control Conditions	Control 1: Exam Type					Any		
	Control 2: Sample Volume					2 mm		
	Control 3: PRF					3906 Hz		
	Control 4: Sample Volume Position					Zone 1		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

Troubleshooting and Maintenance

Cleaning and Disinfecting Transducers

To disinfect the transducer, use the immersion method or the wipe method. You can use the immersion method only if your disinfectant is compatible with it. Check the product labeling. For more information about cleaning and disinfecting the C8e transducer, refer to the ultrasound system user guide.

The following table lists disinfectants that have been tested by SonoSite. See the SonoSite website for updated cleaning and disinfectant information: www.sonosite.com/documentation.

Table 13 does not have the following regulatory information for disinfectants:

- EPA Registration
- FDA 510(k) clearance (liquid sterilant, high level disinfectant)
- CE approval

Prior to use, confirm that the regulatory status of the disinfectant is appropriate for your jurisdiction and use.

Table 13: Disinfectants Compatibility with System and Transducers

Disinfection and Cleaning Solution	Country of Origin	Type	Active Ingredient	C8e
AbcoCide 14	USA	Liquid	Gluteraldehyde	A
Accel Plus	Canada	Wipe	Hydrogen Peroxide	N
Accel TB	Canada	Wipe	Hydrogen Peroxide	N
Accel Wipes	Canada	Wipe	Hydrogen Peroxide	A
Aidal Plus	Australia	Liquid	Gluteraldehyde	A
Alkacide	France	Liquid	Gluteraldehyde	A
Alkazyme	France	Liquid	Quat. Ammonia	A
Anioxyde 1000	France	Liquid	Peracetic Acid	N
Anioxy-Twin	France	Liquid	Peracetic Acid	N
Aquatabs (1000)	Ireland	Tablet	Sodium Dichloroisocyanurate	A
Aquatabs (2000)	Ireland	Tablet	Sodium Dichloroisocyanurate	A
Aquatabs (5000)	Ireland	Tablet	Sodium Dichloroisocyanurate	N

Table 13: Disinfectants Compatibility with System and Transducers (Continued)

Disinfection and Cleaning Solution	Country of Origin	Type	Active Ingredient	C8e
Ascend	USA	Liquid	Quat Ammonia	A
Asepti-HB	USA	Liquid	Quat Ammonia	A
Asepti-Steryl	USA	Spray	Ethanol	A
Asepti-Wipes	USA	Wipe	Propanol (Isopropyl Alcohol)	A
Bacillocid rasant	Germany	Liquid	Glut./Quat. Ammonia	A
Banicide	USA	Liquid	Gluteraldehyde	A
Bleach	USA	Liquid	NaCl Hypochlorite	A
Cavicide	USA	Liquid	Isopropyl	A
CaviWipes	USA	Wipes	Isopropanol	A
Chlor-Clean	UK	Liquid	Sodium Dichloroisocyanurate	A
Cidalkan Lingettes	FRA	Wipes	Ethyl Alcohol	A
Cidex 14	USA	Liquid	Gluteraldehyde	A
Cidex OPA	USA	Liquid	Ortho-phthaldehyde	A
Cidex Plus	USA	Liquid	Gluteraldehyde	A
Cleanisept Wipes	Germany	Wipes	Alkyl Ammonium Chloride	A
Clorox Disinfecting Wipes	USA	Wipe	Isopropanol	A
Control III	USA	Liquid	Quat. Ammonia	A
Coverage Spray	USA	Spray	Quat. Ammonia	A
Denatured Alcohol	USA	Liquid	Ethanol	N
DentaSept	France	Liquid	Quat. Ammonia	N
DisCide Wipes	USA	Wipe	Isopropyl Alcohol	A
DisOPA	Japan	Liquid	Ortho-phthaldehyde	A
Dispatch	USA	Spray	NaCl Hypochlorite	A
Dynacide PA	FRA	Liquid	Peracetic Acid	A
End-Bac II	USA	Liquid	Quat. Ammonia	A

Table 13: Disinfectants Compatibility with System and Transducers (Continued)

Disinfection and Cleaning Solution	Country of Origin	Type	Active Ingredient	C8e
Endozime AW Plus	France	Liquid	Propanol	A
Envirocide	USA	Liquid	Isopropyl	A
Enzol	USA	Cleaner	Ethylene Glycol	A
Expose	USA	Liquid	Isopropyl	A
Gigasept AF	Germany	Liquid	Quat. Ammonia	A
Gigasept FF	Germany	Liquid	Bersteinsäure	N
Gluteraldehyde SDS	USA	Liquid	Gluteraldehyde	A
Hexanios	France	Liquid	Polyhexanide/Quat. Ammonia	A
Hi Tor Plus	USA	Liquid	Chloride	A
Hibiclens	USA	Cleaner	Chlorhexidine	A
Hydrogen Peroxide	USA	Liquid	Hydrogen Peroxide	A
Isopropanol alcohol	n/a	Liquid	Alcohol	N
Kodan Tücher	Germany	Spray	Propanol/Alcohol	A
Kohrsolin ff	Germany	Liquid	Gluteraldehyde	A
Korsolex basic	Germany	Liquid	Gluteraldehyde	N
Korsolex extra	Germany	Liquid	Ethanol/Propanol	A
Lem-O-Quat	USA	Liquid	Alkyl/Chloride	N
LpHse	USA	Liquid	O-phenylphenol	A
Lysol	USA	Spray	Ethanol	N
Lysol IC	USA	Liquid	O-phenylphenol	A
Madacide 1	USA	Liquid	Alkyl Ammonium Chloride	A
Matar	USA	Liquid	O-phenylphenol	A
MetriCide 14	USA	Liquid	Gluteraldehyde	A
MetriCide 28	USA	Liquid	Gluteraldehyde	A
MetriZyme	USA	Cleaner	Propylene Glycol	A

Table 13: Disinfectants Compatibility with System and Transducers (Continued)

Disinfection and Cleaning Solution	Country of Origin	Type	Active Ingredient	C8e
Mikrobak forte	Germany	Liquid	Ammonium Chloride	A
Mikrozid	Germany	Wipe	Ethanol/Propanol	A
Nuclear	France	Spray	Alcohol/Biguanide	A
Precise	USA	Spray	O-phenylphenol	N
Prevention	Canada	Liquid	Hydrogen Peroxide	N
Ruthless	USA	Spray	Quat. Ammonia	A
Sagrosept	Germany	Wipe	Proponal	A
Salvanios pH 7	France	Liquid	Quat. Ammonia	A
Sani-Cloth HB	USA	Wipe	Quat. Ammonia	A
Sani-Cloth Plus	USA	Wipe	Quat. Ammonia	A
Sekusept	GER	Liquid	Gluteraldehyde	A
Sklar (4)	USA	Liquid	Isopropanol	A
Sporicidin	USA	Wipe	Phenol	A
Sporicidin	USA	Liquid	Phenol	A
Staphene	USA	Spray	Ethanol	A
Steranios	France	Liquid	Gluteraldehyde	A
Super Sani-Cloth	USA	Wipe	Isopropyl Alcohol	N
Task 105	USA	Spray	Quat. Ammonia	A
TBQ	USA	Liquid	Alkyl	A
Theracide Plus Wipes	USA	Wipe	Quat. Ammonia	A
Thericide Plus	USA	Liquid	Quat. Ammonia	A
Tor	USA	Liquid	Quat. Ammonia	A
Transeptic	USA	Cleaner	Alcohol	N
Tristel	UK	Liquid	Chlorine Dioxide	A
Tristel Wipes	UK	Wipe	Chlorine Dioxide	N

Table 13: Disinfectants Compatibility with System and Transducers (Continued)

Disinfection and Cleaning Solution	Country of Origin	Type	Active Ingredient	C8e
T-Spray	USA	Spray	Quat. Ammonia	A
T-Spray II	USA	Spray	Alkyl/Chloride	A
Vesphene IIse	USA	Liquid	Sodium/o-Phenylphenate	A
Virex II 256	USA	Liquid	Ammonium Chloride	A
Virex TB	USA	Liquid	Quat. Ammonia	A
Virox 5	CAN	Wipe	Hydrogen Peroxide	A
Virufen	France	Liquid	Alkyl Ammonium Chloride	A
Wavicide -01	USA	Liquid	Gluteraldehyde	N
Wavicide -06	USA	Liquid	Gluteraldehyde	A
Wex-Cide	USA	Liquid	O-phenylphenol	A

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Copyright © 2007 SonoSite, Inc.