# Using CIVCO Products with FUJIFILM Sonosite Systems

---

User Guide

Manufacturer	EC Authorized Representative	Australia Sponsor
FUJIFILM SonoSite, Inc. 21919 30th Drive SE Bothell, WA 98021 USA T: 1-888-482-9449 or 1-425-951-1200 F: 1-425-951-1201	FUJIFILM SonoSite B.V. Joop Geesinkweg 140 1114 AB Amsterdam, The Netherlands	FUJIFILM SonoSite Australasia Pty Ltd 114 Old Pittwater Road BROOKVALE, NSW, 2100 Australia

**CAUTION**

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite Edge, M-Turbo, S Series, SonoSite SII, SonoSite X-Porte, MicroMaxx, NanoMaxx, TITAN, SONOSITE, and the SONOSITE logo are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation. CIVCO is a registered trademark of CIVCO Medical Solutions. Infiniti is a trademark of CIVCO Medical Solutions. DICOM is a registered trademark of the National Electrical Manufacturers Association. All other trademarks are the property of their respective owners.

Patent: [www.sonosite.com/patents](http://www.sonosite.com/patents)

P01765-31

Published Mar 9, 2023

Copyright © 2023 FUJIFILM SonoSite, Inc. All rights reserved.

# Using CIVCO Products with FUJIFILM Sonosite Systems

Introduction .....	1
Document conventions .....	2
Getting help .....	2
Needle guides .....	3
Preparing for the procedure .....	5
Verifying the needle path .....	6
Performing the needle procedure .....	11
Using needle guidelines .....	11
Turning on guidelines .....	12
Moving the depth cursor .....	12
Selecting the angle on a multi-angle bracket .....	12
Clearing the guides .....	12
Cleaning and disinfecting .....	13
Temperature and humidity limits .....	13

## Introduction

This document is used in conjunction with the CIVCO reference guides.

- This document provides information on setting up the ultrasound system and verifying the needle-guide path prior to use and is included with each starter kit.
- The CIVCO reference guides provide information on attaching the bracket, attaching the needle guide, cleaning and disinfecting procedures, sterilization procedures, and additional safety and labeling information and is included with each starter kit and replacement kit.

Please read both documents before using the bracket and needle guide on the ultrasound system. Both documents are for a reader familiar with ultrasound and vascular access techniques or invasive procedures; they do not provide training in sonography, vascular access, or clinical practices.

### **WARNING**

Before using the system and needle guide, you must have training in interventional procedures as dictated by current relevant medical practices as well as in proper operation of the ultrasound system and transducer. During vascular access, the potential exists for serious complications including, without limitation, the following: pneumothorax, arterial puncture, guidewire misplacement, and risks normally associated with local or general anesthesia and recovery.


The needle guide is for use on the following ultrasound systems:

Sonosite ST ultrasound system	FC1 ultrasound system
Sonosite ZX ultrasound system	Sonosite PX ultrasound system
Sonosite LX ultrasound system	Sonosite X-Porte ultrasound system
Sonosite Edge II ultrasound system	Sonosite SII ultrasound system
Sonosite Edge ultrasound system	S Series ultrasound system
M-Turbo ultrasound system	MicroMaxx ultrasound system
NanoMaxx ultrasound system	TITAN high-resolution ultrasound system
Sonosite 180PLUS ultrasound system	iLook 25 personal imaging tool

See the ultrasound system user guide for information about the ultrasound system, transducers, additional safety, and intended uses.

## Document conventions

The document follows these conventions:

- A  **WARNING** describes precautions necessary to prevent injury or loss of life.
- A  **CAUTION** describes precautions necessary to protect the products.
- A  **NOTE** provides supplemental information.
- Numbered and lettered steps must be performed in a specific order.
- Bulleted lists present information in list format but do not imply a sequence.

## Getting help

- FUJIFILM Sonosite Technical Support:

United States and Canada	+1 877-657-8118
Europe and Middle East	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia and Pacific	+61 2 9938 8700
Other regions	+1 425-951-1330 or call your local representative
Fax	+1 425-951-6700
Email	General: <a href="mailto:ffss-service@fujifilm.com">ffss-service@fujifilm.com</a> United Kingdom: <a href="mailto:uk-service@fujifilm.com">uk-service@fujifilm.com</a> Europe, Middle East, and Africa: <a href="mailto:eraf-service@fujifilm.com">eraf-service@fujifilm.com</a> Asia and Pacific: <a href="mailto:ffss-apacme-service@fujifilm.com">ffss-apacme-service@fujifilm.com</a>
Web	<a href="http://www.sonosite.com">www.sonosite.com</a>

## Needle guides

### WARNING

- To avoid device damage or patient injury, do not use the P10, P10x, P10xp, P17, rP19x, rP19xp, P21x, or P21xp needle-guide bracket on patients with pacemakers or medical electronic implants. The needle-guide bracket for these transducers contains a magnet that is used to ensure the bracket is correctly oriented on the transducer. The magnetic field in direct proximity to the pacemaker or medical electronic implant may have an adverse effect.
- System software biopsy guidelines do not support the CIVCO Infiniti Needle Guidance system and should not be used together.

Needle guides assist in guiding a needle. The echoes of the anatomical target and the needle appear on the system display and assist in guiding the needle to the target.

When a needle guide-capable transducer is connected to the system, and a needle guide control is available on the system, the system generates a pair of guidelines that represent the anticipated path of the needle. You can use the guidelines and depth cursor to assist in guiding the needle. See [Turning on guidelines \[12\]](#).

The guidelines displayed on the screen are intended for use with brackets that have installed needle guides only.

### CAUTION

The guidelines display is not available when the centerline graphic is on. The centerline graphic aligns with the center mark of the transducer (if present) and serves as a reference mark for the center of the displayed image. The centerline graphic is not intended for use with needle-guide brackets or with any bracket-associated guidelines. For more information about the centerline graphic, see the system user guide.

The following table lists transducers that are needle guide-capable. Available needle guides may include the following types: fixed angle, multiple-fixed angle, Infiniti (variable angle), and transverse angle.

**Table 1. System and transducer compatibility**

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<sup>a</sup>The starter kit for the L25 transducer series contains sterile needle guides in various gauges and lengths. Use the transducer with either 18-gauge or 21-gauge needles. Refer to the appropriate CIVCO reference guide for which needle length to use to ensure access to the desired target depth. See the CIVCO website at [www.civco.com](http://www.civco.com) for details.

### CAUTION

The guidelines display is not available when the centerline graphic is on. The centerline graphic aligns with the center mark of the transducer (if present) and serves as a reference mark for the center of the displayed image. The centerline graphic is not intended for use with needle-guide brackets or with any bracket-associated guidelines. For more information about the centerline graphic, see the system user guide.

## Preparing for the procedure

### WARNING

To avoid or reduce the risk of infection to the patient:

- Use only FUJIFILM Sonosite or CIVCO approved needle guides, brackets, supplies, components, and accessories. Other brands may not properly fit FUJIFILM Sonosite transducers.
- Use only needle guides compatible with the systems and transducers listed in [Needle guides \[3\]](#). Contact FUJIFILM Sonosite to determine compatible needle guides for a given system and transducer.
- Make sure that the bracket is attached correctly. Refer to the appropriate CIVCO reference guide.
- When using a multi-angle bracket, make sure that you select the same angle (A, B, or C) on both the bracket and the ultrasound system.
- Do not rely on the visibility of the needle tip on the system display. Use other tactile or visual indicators to determine whether you are at or in the vessel. (Example: indentation of anterior wall, decreased resistance as the needle enters the vessel lumen, or blood return in the needle.)
- Follow Universal Precautions when performing interventional procedures.

**WARNING**

Some transducer sheaths contain natural rubber latex and talc, which can cause allergic reactions in some individuals. FUJIFILM Sonosite recommends that you identify your latex and talc-sensitive patients and be prepared to treat allergic reactions promptly.

1. Verify that all components are not damaged.
2. Verify that you are using the correct bracket for the transducer.  
The brackets for the P17, rP19 series, and P21 series transducers have a similar size and/or shape. (The bracket for the P21 series is labeled "P21." The bracket for the rP19 series is labeled "P19.")
3. Confirm that the bracket is attached correctly.
  - Make sure that the bracket is secure. Apply pressure to snap it into place.
  - For the multi-angle bracket on the L38xi, L38xp, C60xf, C60xi, rC60xi, and C60xp transducers, make sure that the angle selected on the bracket matches the angle selected on the ultrasound system. See the ultrasound system user guide for information on system set up.
4. Verify the needle path. See [Verifying the needle path \[6\]](#).
5. In addition to the **WARNINGS** and **CAUTIONS** in this document, observe the following precautionary measures:
  - Needle guides must be attached over a sterile transducer sheath.
  - FUJIFILM Sonosite does not recommend performing an anatomical survey for intracavity biopsies with the needle guide attached.
  - Be aware that the needle may not be visible on the screen. Prior to the procedure, use the depth cursor to determine the target depth. Use this as guidance for depth of needle penetration.
  - Use a straight, new needle and needle guide for each procedure.
  - Refer to the appropriate CIVCO reference guide for information on attaching the bracket and needle guide.
  - Refer to the appropriate CIVCO reference guide for additional information to prepare for the procedure.

**WARNING**

Some transducer sheaths contain natural rubber latex and talc, which can cause allergic reactions in some individuals. FUJIFILM Sonosite recommends that you identify your latex and talc-sensitive patients and be prepared to treat allergic reactions promptly.

## Verifying the needle path

Verify the needle path before the procedure to make sure that the relationship among the system on-screen guidelines, transducer, and needle guide is correct.

### NOTE

This procedure only applies to the types of needle brackets that Sonosite on-screen guidelines support: in-plane fixed, multi-angle brackets, and transverse out-of-plane brackets.

### WARNING


To avoid injury to the patient, observe the following:

- Before performing needle path verification, check that all components are free of damage.
- Use a straight, new needle for each needle path verification procedure.
- Do not use the needle guide if the needle is not following the intended path.
- The needle used for this needle path verification must not be used for the actual procedure. Always use a new, sterile needle and needle guide for each procedure.

The following items are needed to verify the needle path:


- A needle guide-capable transducer
- Bracket (Non-disposable. The type of bracket you use depends on the transducer you are using.)
- Needle-guide kit (disposable)
- New, straight needle
- Ultrasound system
- Tank filled with test object fluid (All World Scientific code 3480) or distilled water.

**Figure 1. Transducer with bracket, needle guide, and needle attached**

1. Needle 2. Needle-guide origin 3. Depth cursor measurement 4. Transducer 5. Bracket

**Figure 2. Transducer with needle in test fluid**


1. Check that all transducer, needle guide, and bracket components are free of damage.
2. Attach the bracket to the transducer, and attach the appropriate needle guide for the needle gauge. Refer to the appropriate CIVCO reference guide for instructions.
3. Set up the tank and transducer, and then insert the needle into the needle guide.
4. Immerse the transducer no more than 1/4-1/2 inch (6.35-12.7 mm) into a tank of test fluid or distilled water.
5. Connect the transducer, and turn on the ultrasound system.
6. Turn on the guidelines. See [Turning on guidelines \[12\]](#).
7. Move the needle down into the tank of test fluid or distilled water until its ultrasound image is visible on the system display.
8. Verify that the needle, as shown on the system display, is between the guidelines.

### **WARNING**

- The guidelines are intended to provide only an indication of the expected path of the needle. Mechanical tolerances and needle flexibility are factors that may result in a needle diverging from the expected path.
- To avoid injury to the patient, if the needle enters from the unexpected side of the display, verify that the needle guide and bracket are correctly attached and that the orientation of the transducer is correct.
- If the needle is still not following the expected path, do not use the needle guide. Contact FUJIFILM Sonosite or your local representative.

**Figure 3. Multi-angle and fixed-angle in-plane guides**

**NOTE**

The depth cursor measurement shown on the ultrasound system measures the distance from the tip of the needle to the origin of the needle guide.

**Figure 4. Out-of-plane guide**

### **WARNING**

- The guidelines are intended to provide only an indication of the expected path of the needle. Mechanical tolerances and needle flexibility are factors that may result in a needle diverging from the expected path.
- To avoid injury to the patient, if the needle enters from the unexpected side of the display, verify that the needle guide and bracket are correctly attached and that the orientation of the transducer is correct.
- If the needle is still not following the expected path, do not use the needle guide. Contact FUJIFILM Sonosite or your local representative.

## Performing the needle procedure

### WARNING

To avoid injury or reduce risk of infection to the patient:

- Check that a needle path verification procedure has been performed before performing the needle procedure.
- Use a straight, new, sterile needle for each needle procedure.
- Sterilize reusable needle guides after each use.
- Thin needles can bend when entering tissue. Actual position must be verified by identifying the echoes from the needle.
- If the bracket, guide, or needle breaks during the procedure, verify that all pieces of equipment are recovered.
- If unable to visualize the needle at the target, ensure that the correct needle guide has been selected and that the bracket and needle guide are correctly attached.

### CAUTION

To avoid damage to the transducer, use only gels recommended by FUJIFILM Sonosite. Using other gels can damage the transducer and void the warranty. If you have questions about gel compatibility, contact FUJIFILM Sonosite or your local representative.

In addition to the **WARNINGS** and **CAUTIONS** in this document, observe the following precautionary measures during a procedure:

- Use acoustic coupling gel for proper acoustic coupling. Apply a liberal amount of gel between the transducer and the body. Although most gels provide suitable acoustic coupling, some gels are incompatible with some transducer materials. A sample of Aquasonic gel is provided in the needle-guide kit.
- Use sterile gel when performing a needle procedure.
- The guidelines are intended to provide only an indication of the expected path of the needle. Actual position must be verified by identifying the echoes from the needle.

## Using needle guidelines

### WARNING

Avoid using guidelines and depth indicators in Zoom or Dual mode. Viewing in Zoom or Dual mode can obscure important physical structures and reference points.

Guidelines are an optional feature to help with needle guidance. The availability of this feature depends on the type of transducer you are using.

## Turning on guidelines

1. On the ultrasound system, select the correct exam type.
2. Optimize the image.
3. In 2D imaging, do one of the following:
  - (Sonosite PX, Sonosite LX, Sonosite ZX, Sonosite ST) Tap the **Needle Guide** control on. If this control is hidden, tapping **+ More Controls** displays it.
  - (Sonosite X-Porte) Tap **Needle Guide** and then tap **On**. If this control is hidden, tapping **More Controls** displays it. Guidelines are not available in Zoom or Dual mode.
  - (M-Turbo, S Series, Sonosite Edge, Sonosite Edge II, Sonosite SII) Select **Guide** on-screen.
  - (180PLUS) Press the **TEXT/PICTO** key. If pictograph is off, select the Biopsy icon to display the guidelines. If pictograph is on, select the Forward icon and select the Biopsy icon.
  - (FC1) On the touchscreen, tap **Needle Guide**.
  - (Other FUJIFILM Sonosite systems) Select **Biopsy** or **Guide** on-screen.
4. See the ultrasound system user guide for more detailed instructions on basic operation and system settings.

## Moving the depth cursor

Do one of the following:

- (Sonosite LX) For transverse-angle brackets, use the monitor touchscreen or depth control buttons to change depth. The needle depth indicators on the clinical monitor adjust accordingly.
- (Sonosite PX, Sonosite ZX) For transverse-angle, out-of-plane brackets, use the depth control buttons to change depth. The needle depth indicators on the clinical monitor adjust accordingly.
- (Sonosite ST) For transverse-angle brackets, use the depth control buttons on the touchscreen to change depth. The needle depth indicators on the clinical monitor adjust accordingly.
- (Sonosite X-Porte) On the control panel, use your finger to drag the depth cursor.
- (S Series, Sonosite SII) Turn the **Depth** knob.
- (Sonosite Edge, Sonosite Edge II) Press the buttons to increase or decrease depth.
- (180PLUS) Use the trackball. If pictograph is on, press the back icon to return to the previous menu.

The target distance measurement value changes to reflect the distance between the needle guide origin and the depth cursor. The origin is the point at which the needle enters the needle guide or needle insert. See [Figure 1, "Transducer with bracket, needle guide, and needle attached" \[7\]](#).

## Selecting the angle on a multi-angle bracket

Turn the needle guides on and do one of the following:

- (Sonosite X-Porte) Tap **Guide A**, **Guide B**, or **Guide C**.
- (Other FUJIFILM Sonosite systems) Select **A**, **B**, or **C**.

## Clearing the guides

Do one of the following:

- (Sonosite PX, Sonosite LX, Sonosite ZX, Sonosite ST) Tap the **Needle Guide** control off.
- (Sonosite X-Porte) Do either of the following:
  - Tap **Needle Guide** and then tap **Off** (multi-angle).
  - Tap **Needle Guide** (single-angle).
- (M-Turbo, Sonosite Edge, Sonosite Edge II, S Series, Sonosite SII) Do any of the following:
  - Select the angle again: **A**, **B**, or **C** (multi-angle).
  - Exit angle selection and press **Guide** (multi-angle).
  - Press **Guide** or **Biopsy** (single-angle).

## Cleaning and disinfecting

Before cleaning, disinfecting, or sterilizing the bracket, refer to the appropriate CIVCO reference guide.

## Temperature and humidity limits

For environmental limits, including shipping and storing, refer to the appropriate CIVCO reference guide.

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