PowerPark User Guide

For Sonosite ST, Sonosite PX, and Sonosite LX

Manufacturer	EC Authorized Representative	Australia Sponsor
FUJIFILM SonoSite, Inc. 21919 30th Drive SE Bothell, WA 98021 USA T: 1-888-482-9449 or 1-425-951-1200 F: 1-425-951-1201	FUJIFILM SonoSite B.V. Joop Geesinkweg 140 1114 AB Amsterdam, The Netherlands	FUJIFILM SonoSite Australasia Pty Ltd 114 Old Pittwater Road BROOKVALE, NSW, 2100 Australia

**CAUTION**

United States federal law restricts this device to sale by or on the order of a physician.

PowerPark and SONOSITE are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation. All other trademarks are the property of their respective owners.

Patent: www.sonosite.com/patents

P23655-03

Published Mar 9, 2023

Copyright © 2023 FUJIFILM SonoSite, Inc. All rights reserved.

PowerPark User Guide

Introduction	1
Document conventions	1
Getting help	2
About the PowerPark	2
Docking and undocking the stand	3
Troubleshooting	4
Turning the dock circuit breaker on or off	4
Resetting the stand module circuit breakers	4
Cleaning and disinfecting	6
Cleaning and disinfecting the PowerPark	6
Safety	7
Electrical safety	7
Equipment safety	8
Hazardous materials	8
Electromagnetic compatibility (EMC)	8
Electrostatic discharge	9
Labeling symbols	9
Specifications	12

Introduction

The *PowerPark User Guide* explains how to use, clean, and disinfect the PowerPark. It also provides information on safety and specifications. To attach or remove the PowerPark from the stand, refer to the *PowerPark Assembly Instructions*. For information on the stand or ultrasound system, see the user guide or the ultrasound system on-board help.

WARNING

To avoid the risk of patient injury, misdiagnosis, or user injury, read all the warnings in the user guide, supplemental guides, or the ultrasound system on-board help.

Document conventions

The document follows these conventions:

- A  **WARNING** describes precautions necessary to prevent injury or loss of life.
- A  **CAUTION** describes precautions necessary to protect the products.
- A  **NOTE** provides supplemental information.
- Numbered and lettered steps must be performed in a specific order.
- Bulleted lists present information in list format but do not imply a sequence.

Symbols and terms used on the PowerPark components and on the ultrasound stand, system, and transducers are explained in the ultrasound system user guide.

Getting help

- FUJIFILM Sonosite Technical Support:


United States and Canada	+1 877-657-8118
Europe and Middle East	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia and Pacific	+61 2 9938 8700
Other regions	+1 425-951-1330 or call your local representative
Fax	+1 425-951-6700
Email	General: ffss-service@fujifilm.com United Kingdom: uk-service@fujifilm.com Europe, Middle East, and Africa: eraf-service@fujifilm.com Asia and Pacific: ffss-apacme-service@fujifilm.com
Web	www.sonosite.com

About the PowerPark

The PowerPark is a power docking station for the Sonosite ST, Sonosite PX, and Sonosite LX ultrasound systems and stand. It lets you dock the stand—without manually connecting the power cord—to charge the stand battery (if present) and the ultrasound system battery.

The PowerPark includes two parts: the dock (which remains on the floor against the wall) and the stand-mounted module (model P25003) which attaches to the underside of the stand.


Figure 1. The PowerPark docking station


PowerPark includes two main components: the stand-mounted module, and the dock. 1. Ultrasound system stand 2. Stand-mounted module 3. Dock

Docking and undocking the stand


1. Make sure that the system's power cord is disconnected from AC power.

2. Wheel the stand backward into the dock.
The stand module (on the underside of the stand) snaps into place as it connects with the dock.

3. To keep the stand in its place, lock one or more of the wheels.
4. Verify PowerPark is charging by confirming that the connection indicator on the stand base glows green.

5. To undock the stand, unlock the stand's wheels, then wheel the stand away from the dock.
The connection indicator turns off.

Troubleshooting

Table 1. Symptoms and solutions


Symptom	Solution
PowerPark does not charge the battery	<ul style="list-style-type: none"> • Make sure that the electrical outlet has power. • Make sure that the dock circuit breaker is turned on. • Reset the stand-module circuit breaker if tripped. • Make sure that the stand module is fully connected to the dock.
Stand module does not connect to the dock	Make sure that the floor is clean, flat, and rigid. Avoid thick carpeting.

Turning the dock circuit breaker on or off

The circuit breaker on the PowerPark dock is typically set to the ON position but trips if there is a power event (for example, a power surge).

The circuit breaker is a means of disconnection from mains power.


1. Undock the stand.
2. Disconnect the power cord from the wall outlet, and pull the PowerPark dock away from the wall.
The circuit breaker is now accessible.
3. Press the circuit breaker to the desired position: **ON** (|) or **OFF** (O).

Resetting the stand module circuit breakers

One or both of the circuit breakers on the stand module trip if there is a power event (for example, a power surge) and may need to be reset.


1. Locate the two breakers on the end of the power module under the stand base, near the stand height adjustment pedal.

2. Press in on each circuit breaker switch until it closes. See the following figure.
 - On a correctly closed circuit breaker, the switch is pressed in.
 - On a tripped circuit breaker, the switch is popped up and exposes a white plastic area.

**NOTE**

If you cannot press down the switch and if no white area is exposed, the circuit breaker is closed.

1. Breaker closed 2. Breaker tripped

Cleaning and disinfecting

WARNING

- To avoid electrical shock, before cleaning, always disconnect the PowerPark dock from AC power.
- To avoid injury, always use protective eyewear and gloves when performing cleaning and disinfecting procedures.
- To avoid infection:
 - Ensure that the solution expiration date has not passed.
 - The level of disinfection required for a product is determined by the type of tissue it has contacted or will contact during use. Ensure that the solution strength and duration of contact are appropriate for the equipment. For more information, see the disinfectant label instructions and the recommendations of the Association for Professionals in Infection Control and Epidemiology (APIC) and FDA.

CAUTION

Use only approved cleaners or disinfectants on surfaces.

The exterior surfaces of the PowerPark stand module and the PowerPark dock should be kept clean. Whenever you clean the system stand, be sure to also clean the attached PowerPark stand module. Clean and disinfect the PowerPark dock periodically, or whenever it appears dirty or comes into contact with biomatter. Refer to the cleaners and disinfection tool available at www.sonosite.com/support/cleaners-disinfectants for a list of approved cleaners and disinfectants.

CAUTION

Do not allow cleaning solution or disinfectant to come into contact with the PowerPark electrical connectors.

Cleaning and disinfecting the PowerPark

1. Disconnect the PowerPark components from any power source:
 - For the PowerPark stand module, turn off the ultrasound system and any accessory equipment. Disconnect the stand from the PowerPark dock or unplug the power cord.
 - For the PowerPark dock, disconnect it from its AC power.

2. Clean the stand module and PowerPark dock surfaces to remove any gel or debris. Use the following procedure:
 - a. Use either a pre-moistened wipe or a soft cloth dampened with cleaner or disinfectant.
 - b. Remove all gel and debris from the PowerPark components.
3. With a new wipe, clean all parts of the PowerPark components by wiping from the clean areas to the soiled areas. This method helps to avoid cross-contamination. Follow the chemical manufacturer's instructions regarding temperature and wet contact time. Monitor the PowerPark parts for wet appearance, and reapply with a new wipe if no longer wet.
4. Allow the PowerPark components to air dry in a clean, well-ventilated space.

Safety

For maximum safety, observe all of the warnings and cautions in this section. Please report any serious safety incident that occurs in relation to the PowerPark and the system to FUJIFILM Sonosite and to the competent authority of the country in which the user and patient are located.

Electrical safety

Table 2. Electrical safety classification

Class I equipment	The PowerPark is classified as Class I equipment.
Non AP/APG	The PowerPark is not suitable for use in the presence of flammable anesthetics.

WARNING

- To avoid the risk of a shock hazard:
 - Do not plug the power cord of the PowerPark dock or PowerPark stand module into a multiple portable socket outlet (MPSO) or use an extension cord.
 - Make sure the AC power cord is disconnected from the PowerPark stand module whenever the stand is docked.
- To avoid the risk of electrical shock and possible fire:
 - Use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when the equipment is connected to a receptacle marked "Hospital Only," "Hospital Grade," or the equivalent. The grounding wire must not be removed or circumvented.
 - Use only accessories and peripherals recommended by FUJIFILM Sonosite, including the power supply. Use the PowerPark only within the voltage and current limits specified on its ratings label.
 - Inspect the power supply, plugs, cables, and power cords on a regular basis for damage.
 - Do not place the PowerPark dock in a wet environment.
- To avoid system malfunction and the risk of electrical shock, follow the PowerPark assembly instructions.

Equipment safety

CAUTION

- Excessive bending or twisting of cables can cause a failure or intermittent operation.
- Improper cleaning or disinfecting of any part of the accessories and peripherals can cause permanent damage. Refer to the cleaners and disinfection tool available at www.sonosite.com/support/cleaners-disinfectants for a list of approved cleaners and disinfectants.

Hazardous materials

WARNING

Products and accessories may contain hazardous materials. When disposing of products and accessories, be environmentally responsible and meet federal and local regulations for disposing hazardous materials.

Electromagnetic compatibility (EMC)

See the Safety section of the ultrasound system user guide for information on electromagnetic compatibility with IEC 60601-1-2, including the Manufacturer's Declaration.

Electrostatic discharge

WARNING

- Unless following ESD precautionary procedures, do not connect to or touch (with body or hand-held tools) pins (contacts) of connectors that have the ESD Sensitive Devices symbol:

- If the symbol is on a border surrounding multiple connectors, the symbol pertains to all connectors within the border.
- ESD precautionary procedures include the following:
 - Receive training about ESD, including the following (at a minimum): an introduction to the physics of electrostatic charge, the voltage levels that can occur in normal practice, and the damage that can occur to electronic components if equipment is touched by an individual who is electrostatically charged.
 - Prevent the buildup of electrostatic charge; for example, use humidification, conductive floor coverings, non-synthetic clothing, ionizers, and minimizing insulating materials.
- Ground yourself to discharge any static charge. Use a wrist strap to bond yourself to earth ground.

Labeling symbols

The following symbols are used on the products, packaging, and containers.

Table 3. Standards labeling symbols

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Manufacturer	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.1.1	Indicates the medical device manufacturer
	Serial Number	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.1.7	Indicates the manufacturer's serial number so that a specific medical device can be identified

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Catalog Number	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.1.6	Indicates the manufacturer's catalogue number so that the medical device can be identified
	Caution	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.4.4	Indicates that caution is necessary when operating the device or control close to where the symbol is placed
	Fragile handle with care	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.1	Indicates a medical device that can be broken or damaged if not handled carefully
	Keep dry	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.4	Indicates a medical device that needs to be protected from moisture
	Temperature limit	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.7	Indicates the temperature limits to which the medical device can be safely exposed
	Atmospheric pressure limitations	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.9	Indicates the range of atmospheric pressure to which the medical device can be safely exposed

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Humidity limitation	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.8	Indicates the range of humidity to which the medical device can be safely exposed
	Stacking limit by number	ISO 7000 Graphical symbols for use on equipment	2403	Do not stack over n high, where n represents the number on the label
	Consult instructions for use	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.4.3	Indicates that the operating instructions should be considered when operating the device or control close to where the symbol is placed
	Corrugated recycle	—	—	Shipping box is made of corrugated cardboard and should be recycled accordingly
	RESY – Recycling Symbol	—	—	Paper recycle
	Recycle: Electrical Equipment	BS EN 50419 Marking of Electrical and Electronic Equipment in accordance with Directive 2012/19/EU for the Waste of Electrical and Electronic Equipment (WEEE) and Directive 2006/66/EC on Batteries and Accumulators and Waste Batteries and Accumulators	Annex IX	Do Not Throw in Trash
	CE Marking	—	—	Signifies European Technical Conformity
	European community authorized representative	ISO 15223-1 Medical devices - symbols to be used with medical device labels, labeling and in-	5.1.2	Indicates the Authorized representative in the European Community

Symbol	Title	Standards development organization (SDO)	Reference number	Description
		formation to be supplied		
	Medical Device	EU MDR	EU MDR Annex I, 23.2 (q)	Indicates the item the label is adhered to is categorized as a medical device per the MDR, Annex 1, 23.2, q.
	Alternating current	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5032	Indicates on the rating plate, that the equipment is suitable for alternating current only, in order to identify appropriate terminals
	Canadian Standard Association Certification Mark	—	—	CSA certification mark signifying that the product complies with the applicable CSA and ANSI/UL requirements and is authorized for use in Canada and the US
	ETL (Electronic Testing Laboratories) certification mark	—	—	Proof product complies to North American safety standards
	China Pollution Control (10)	Ministry of Industry and Information Technology	—	Pollution Control Logo. (Applies to all parts/products listed in the China RoHS disclosure table. May not appear on the exterior of some parts/products because of space limitations.)

Specifications

Dimensions

Approximate dimensions of the PowerPark dock are as follows:

- Width: 15.00 in (38.1 cm)
- Depth: 4.88 in. (12.4 cm)
- Length: 4.88 in. (12.4 cm)
- Weight (including cord): 11.3 lbs. (5.1 kg)

Environmental limits

- Operating: 10 to 40°C; 15 to 95% R.H.; 700 hPa to 1060 hPa
- Shipping/Storage: -35 to 65°C; 15 to 95% R.H.; 500 hPa to 1060 hPa

Electrical

- PowerPark Dock Supply Input: 100–240 VAC, 50–60 Hz, 8.0–3.3 A
- PowerPark Stand-Module Supply Input: 100–240 VAC, 50–60 Hz, 6.0–2.5 A

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