

Sonosite Edge II

P21x Transducer User Guide Supplement

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

EDGE, SONOSITE and the SONOSITE logo are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

US 9,151,832; US 8,956,296; US 8,861,822; US 8,858,436; US 8,834,372; US 8,805,047; US 8,527,033; US 8,500,647; US 8,439,840; US 8,398,408; US 8,355,554; US 8,216,146; US 8,213,467; US 8,147,408; US 8,137,278; US 8,088,071; US 8,066,642; US 8,052,606; US 7,819,807; US 7,804,970; US 7,740,586; US 7,686,766; US 7,604,596; US 7,591,786; US 7,588,541; US 7,534,211; US 7,449,640; US 7,169,108; US 6,962,566; US 6,648,826; US 6,575,908; US 6,569,101; US 6,471,651; US 6,416,475; US 6,383,139; US 6,364,839; US 6,203,498; US 6,135,961; US 5,893,363; US 5,817,024; US 5,782,769; US 5,722,412; AU: 730822; AU: 727381; CA 2,372,152; CA: 2,371,711; CN103237499; CN101231457; CN 98108973.9; CN 98106133.8; CN 97113678.5; DE 69831698.3; DE 69830539.6; DE 69730563.5; DE 602004027882.3; DE 602004023816.3; DE: 60034670.6; DE 60029777.2; EP 1589878; EP 1552792; EP 1180971; EP 0875203; EP 0815793; EP 1180970; EP 0881492; ES 2229318; ES 159878; ES 1552792; ES 0881492; FR 158978; FR 1552792; FR 1180970; FR 0881492; FR 0875203; FR 0815793; GB 158978; GB 1552792; GB 1180971; GB 1180970; GB 0881492; GB 0875203; GB 0815793; IT 1589878; IT 1552792; IT 0881492; IT 0815793; JP 5782428; JP 4696150; KR 532359; KR 528102; NO 326814; NO 326202 and pending.

All other trademarks are the property of their respective owners.

Part number: P20025-02

Publication date: April 2020

Copyright © 2020 FUJIFILM SonoSite, Inc. All rights reserved.



Sonosite Edge II P2 1x Transducer User Guide Supplement

Introduction	1
Document conventions	1
Getting Help	2
Imaging modes	3
Cleaning and disinfecting	3
Safety	4
Compatible accessories and peripherals	4
Acoustic output	4
Guidelines for reducing MI and TI	4
Output display	5
Acoustic output tables	6

Introduction

This user guide supplement provides information on the P2 1x transducer, compatible with the Sonosite Edge II ultrasound system.

Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.

For a description of labeling symbols that appear on the product, see "Labeling Symbols" in the ultrasound system user guide.

Getting Help

For technical support, please contact FUJIFILM Sonosite as follows:

Phone (U.S. or Canada)	877-657-8118
Phone (outside U.S. or Canada)	425-951-1330, or call your local representative
Fax	425-951-6700
Email	ffss-service@fujifilm.com
Web	www.sonosite.com
Europe Service Center	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia Service Center	+65 6380-5581

Printed in the United States of America.

Imaging modes

Table 1: Available exam types and imaging modes

Transducer	Exam type	Imaging mode			
		2D M Mode	Color ^a	PW Doppler ^b	CW Doppler
P21x	Abdomen	✓	CVD, CPD	✓	—
	Cardiac	✓	CVD, Var	✓	✓
	Lung	✓	CVD, CPD	✓	—

^aVar = Color Doppler Variance, which is available in the cardiac exam only. CPD = Color Power Doppler, which is available in all exams except the cardiac exam type. CVD = Color Velocity Doppler.

^bFor the cardiac exam type, PW TDI is also available.

Cleaning and disinfecting

Table 2: Approved cleaners

Product	Compatible transducer	Minimum wet contact time ^a
SaniCloth AF3 ^b	P21x	3 minutes
PI-Spray II	P21x	10 minutes

^a For maximum effectiveness, the component being cleaned must remain wet with disinfectant for a minimum period of time.

^b Qualified for use as an intermediate-level disinfectant for mycobacteria.

Refer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.

Table 3: Approved high-level compatible disinfectants

Disinfectant ^a	Compatible transducer	Temperature	Disinfectant soak duration
Cidex	P21x	25°C, 77°F	45 minutes

^aRefer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.

Table 3: Approved high-level compatible disinfectants

Disinfectant ^a	Compatible transducer	Temperature	Disinfectant soak duration
Cidex OPA	P21x	20°C, 68°F	12 minutes
^a Refer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.			

Safety

WARNING

To avoid device damage or patient injury, do not use the P21x needle-guide bracket on patients with pacemakers or medical electronic implants. The needle-guide bracket for P21x contains a magnet that is used to ensure that the bracket is correctly oriented on the transducer. The magnetic field in direct proximity to the pacemaker or medical electronic implant may have an adverse effect.

Compatible accessories and peripherals

Table 4: Accessories and peripherals

Description	Maximum Cable Length
P21x transducer ^a	6.0 ft / 1.8 m
^a For transducers, maximum cable length is measured between the strain reliefs. The stated lengths do not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, or inside the transducer connector.	

Acoustic output

Guidelines for reducing MI and TI

Table 5: Guidelines for reducing mechanical index (MI)

Transducer	Depth
P21x	↑
↓ Decrease or lower setting of parameter to reduce MI. ↑ Increase or raise setting of parameter to reduce MI.	

Table 6: Guidelines for reducing thermal index (TI)

Transducer	CPD Settings						PW Settings
	Box Width	Box Height	Box Depth	PRF	Depth	Optimize	
P21x		↓		↓	↑		↓ (PRF)
↓ Decrease or lower setting of parameter to reduce TI. ↑ Increase or raise setting of parameter to reduce TI.							

Output display

Table 7: TI or MI ≥ 1.0

Transducer	Index	2D/ M Mode	CPD/ Color	PW Doppler	CW Doppler
P21x	MI	Yes	Yes	Yes	No
	TIC, TIB, or TIS	Yes	Yes	Yes	Yes

Transducer surface temperature rise

Table 8: Transducer Surface Temperature Rise, External Use (°C)

Test	P2 1x
Still air	17.2
Simulated use	8.5

Acoustic output tables

Table 9: Transducer model: P21x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	(a)		(a)		2.3
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.03					
	P (mW)		#		#		171.7
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	3.4					
	$z_{pii,\alpha}$ (cm)	3.4					
	f_{awf} (MHz)	1.83	#		#		1.94
Other information	p_{rr} (Hz)	4444					
	s_{rr} (Hz)	34.7					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	194					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	37.1					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	48.6					
	p_r at z_{pii} (MPa)	2.53					
Operating controls	Exam type	Crd					Crd
	Optimization	Gen/Pen					Pen
	Depth (cm)	4.7					27
	MB/THI	On					Off
	Sector Width	Any					Narrow

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Not applicable for this transducer/mode.

Table 10: Transducer model: P2 1x Operating mode: M Mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	(a)		1.4		1.1
Index component value			#	#	0.2	1.4	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.1					
	P (mW)		#		40.1		79.7
	$P_{1 \times 1}$ (mW)		#		16.8		
	z_s (cm)			#			
	z_b (cm)					4.90	
	z_{MI} (cm)	3.645					
	$z_{pii,\alpha}$ (cm)	3.645					
	f_{awf} (MHz)	1.93	#		1.93		1.94
Other information	p_{rr} (Hz)	800					
	s_{rr} (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	237.4					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	171.5					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	285.6					
	p_r at z_{pii} (MPa)	2.679					
Operating controls	Exam type	Abd			Abd		Abd
	Optimization	Any			Any		Pen
	Depth (cm)	7.5			10/13		32
	THI	On			On		Off

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Not applicable for this transducer/mode.

Table 11: Transducer model: P2 1x Operating mode: CPD/Color

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	1.3		1.3		2.4
Index component value			1.3	1.3	1.3	1.3	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.03					
	P (mW)		134.4		134.4		134.4
	$P_{1 \times 1}$ (mW)		121.7		121.7		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	3.4					
	$z_{pii,\alpha}$ (cm)	3.4					
	f_{awf} (MHz)	1.83	2.16		2.16		2.16
Other information	pr_r (Hz)	1819					
	srr (Hz)	28.4					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	194					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sij,\alpha}$ (mW/cm ²)	15.2					
	I_{spta} at z_{pii} or z_{sij} (mW/cm ²)	19.9					
	p_r at z_{pii} (MPa)	2.53					
Operating controls	Mode	CVD	CVD		CVD		CVD
	Exam type	Crđ	Lung		Lung		Lung
	2D Optimization/Depth	Pen/4.7 cm	Pen/4.7 cm		Pen/4.7 cm		Pen/4.7 cm
	Color Optimization/PRF	Any/Any	Med/2500 Hz		Med/2500 Hz		Med/15 24 Hz
	THI	On	N/A		N/A		N/A
	Color box size	Any	Small		Small		Small

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Not applicable for this transducer/mode.

Table 12: Transducer model: P2 1x Operating mode: CW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	1.0		3.6		3.1
Index component value			1.0	0.9	1.0	3.6	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		108.8		108.8		108.8
	P_{1x1} (mW)		108.8		108.8		
	z_s (cm)			1.2			
	z_b (cm)					1.2	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	2.0		2.0		2.0
Other information	pr (Hz)	#					
	srr (Hz)	#					
	n_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sij,\alpha}$ (mW/cm ²)	#					
	I_{spta} at z_{pii} or z_{sij} (mW/cm ²)	#					
	p_r at z_{pii} (MPa)	#					
Operating controls	Exam type		Crd		Crd		Crd
	Sample volume position		Zone 0		Zone 0		Zone 0

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Not applicable for this transducer/mode.

Table 13: Transducer model: P2 1x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	1.3		3.7		2.8
Index component value			0.8	1.3	1.0	3.7	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.73					
	P (mW)		200.3		93.8		201.2
	$P_{1 \times 1}$ (mW)		78.4		93.8		
	z_s (cm)			3.1			
	z_b (cm)					0.60	
	z_{MI} (cm)	5.0					
	$z_{pii,\alpha}$ (cm)	5.0					
	f_{awf} (MHz)	2.15	2.22		2.17		2.12
	pr (Hz)	1562					
Other information	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	216					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sij,\alpha}$ (mW/cm ²)	400.8					
	I_{spta} at z_{pii} or z_{sij} (mW/cm ²)	830.4					
	p_r at z_{pii} (MPa)	2.5					
Operating controls	Exam type	Crd	Crd		Crd		Crd
	Sample volume size (mm)	1	3		1		1
	PRF (Hz)	1563	3906		15625		3906
	Sample volume position	Zone 2	Zone 4		Zone 0		Zone 5

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Not applicable for this transducer/mode.

FUJIFILM
SONOSITE

P28967-01

