**SonoSite** Edge

**SonoSite** S

**SonoSite** M-Turbo

***User Guide Vet Supplement***

---

**Manufacturer**

FUJIFILM SonoSite, Inc.  
21919 30th Drive SE  
Bothell, WA 98021 USA  
T: 1-888-482-9449 or 1-425-951-1200  
F: 1-425-951-1201

**EC Authorized Representative**

Emergo Europe  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands

**Australia Sponsor**

FUJIFILM SonoSite Australasia Pty Ltd  
114 Old Pittwater Road  
BROOKVALE, NSW, 2100  
Australia

**Caution**

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, the SonoSite logo, Edge, S Series and M-Turbo are trademarks and registered trademarks of FUJIFILM Sonosite, Inc. in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Patents: .

Part number: P13329-05

Publication date: January 2016

Copyright © 2016 FUJIFILM SonoSite, Inc. All rights reserved.

# SonoSite Edge, SonoSite S Series, SonoSite M-Turbo User Guide Vet Supplement

<b>Imaging</b> .....	<b>1</b>
2D controls .....	1
Imaging modes and exams available by transducer .....	1
<b>Measurements and Calculations (Vet Only)</b> .....	<b>2</b>
Cardiac calculations .....	2
Volume calculations .....	2
Percent reduction calculations .....	2
<b>Safety</b> .....	<b>3</b>
Acoustic output tables .....	3

This user guide supplement describes the C11x transducers available for veterinary use. It supplements the existing C11x-transducer information in the *Edge User Guide*, *S Series User Guide*, and *M-Turbo Ultrasound System User Guide*.

## Imaging

### 2D controls

#### Note

The Sector control is unavailable on the veterinary C11x transducer.

### Imaging modes and exams available by transducer

The Abv exam type is available on the following ultrasound systems:

- ▶ SonoSite Edge (Vet)
- ▶ SonoSite S Series (Vet)
- ▶ SonoSite M-Turbo (Vet)

**Table 1: Imaging modes and exams available**

Transducer	Exam Type <sup>1</sup>	Imaging Mode				
		2D <sup>2</sup> M Mode	CPD <sup>3</sup>	Color <sup>3</sup>	PW Doppler	CW Doppler
C11x	Crd	✓	—	✓	✓	—
	Abv	✓	✓	✓	✓	—

1. Exam type abbreviations are as follows: Crd = Cardiac, Abv = Abdomen Vet
2. The optimization settings for 2D are Res, Gen, and Pen.
3. The optimization settings for Color are low, medium, and high (flow sensitivity) with a range of PRF settings.

## Measurements and Calculations (Vet Only)

### Cardiac calculations

Transducer	Exam Type
C11x	Cardiac

### Volume calculations

Transducer	Exam Type
C11x	Abv

### Percent reduction calculations

Transducer	Exam Type
C11x	Abv

# Safety

## Acoustic output tables

**Table 2: Transducer Model: C11x Operating Mode: Color**

Index Label			M.I.	TIS			TIB	TIC	
				Scan	Non-scan		Non-scan		
					$A_{aprt} \leq 1$	$A_{aprt} > 1$			
Global Maximum Index Value			(a)	(a)	—	—	—	1.1	
Associated Acoustic Parameter	Pr.3	(MPa)	#						
	$W_0$	(mW)		#	—		—	39.4	
	min of [ $W_{.3}(z_1), I_{TA.3}(z_1)$ ]	(mW)				—			
	$z_1$	(cm)				—			
	$z_{bp}$	(cm)				—			
	$z_{sp}$	(cm)	#				—		
	$d_{eq}(z_{sp})$	(cm)					—		
	$f_c$	(MHz)	#	#	—	—	—	4.37	
	Dim of $A_{aprt}$	X (cm)		#	—	—	—	—	1.12
Y (cm)			#	—	—	—	—	0.50	
Other Information	PD	( $\mu$ sec)	#						
	PRF	(Hz)	#						
	$p_r@PII_{max}$	(MPa)	#						
	$d_{eq}@PII_{max}$	(cm)					—		
	Focal Length	$FL_x$ (cm)			#	—	—		4.29
		$FL_y$ (cm)			#	—	—		4.40
$I_{PA.3}@MI_{max}$	( $W/cm^2$ )	#							
Operating Control Condition	Control 1: Mode							Color	
	Control 2: Exam Type							Cardiac	
	Control 3: PRF							2033	
	Control 4: Optimization/Depth							Low/5.1	
	Control 5: Color Box Position/Size							Top/Short & Narrow	

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