
PowerPark User Guide

For Sonosite PX

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: +1-888-482-9449 or
+1-425-951-1200
F: +1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution United States federal law restricts this device to sale by or on the order of a physician.

PowerPark and SONOSITE are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Patents: US 8,398,408; US 7,591,786; US D625014; and US D625015.

Part Number: P23655-01

Publication date: August 2020

Copyright © 2020 FUJIFILM SonoSite, Inc. All rights reserved.

PowerPark User Guide for Sonosite PX

Introduction	1
Document conventions	1
Getting help.....	2
About the PowerPark	2
Docking the stand	3
Troubleshooting	4
Cleaning and disinfecting	6
Safety	7
Electrical safety	7
Equipment safety.....	9
Hazardous materials.....	9
Electromagnetic compatibility (EMC).....	9
Electrostatic discharge	10
Labeling symbols.....	10
Specifications	14
Dimensions.....	14
Environmental limits.....	14
Electrical.....	14

Introduction

The *PowerPark User Guide* explains how to use, clean, and disinfect the PowerPark. It also provides information on safety and specifications. To attach or remove the PowerPark from the stand, refer to the *PowerPark Assembly Instructions*. For information on the stand or ultrasound system, see the user guide or the ultrasound system on-board help.

WARNING To avoid the risk of patient injury, misdiagnosis, or user injury, read all the warnings in the user guide, supplemental guides, or the ultrasound system on-board help.

Document conventions

The user guide follows these conventions:

- ▶ A  **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A  **Caution** describes precautions necessary to protect the products.
- ▶ A  **Note** provides supplemental information.

- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ♦.

For a description of labeling symbols that appear on the product, see **“Labeling symbols”** on page 10.

Getting help

For technical support, please contact Fujifilm Sonosite as follows:

Phone (U.S. or Canada)	+1-877-657-8118
Phone (outside U.S. or Canada)	+1-425-951-1330, or call your local representative
Fax	+1-425-951-6700
Email	ffss-service@fujifilm.com
Web	www.sonosite.com
Europe Service Center	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia Service Center	+61 2 9938 8700 ffss-apacme-service@fujifilm.com

Printed in the U.S.

About the PowerPark

The PowerPark is a power docking station for the Sonosite PX ultrasound system and stand. It lets you dock the stand—without manually connecting the power cord—to charge the stand battery (if present) and the ultrasound system battery.

The PowerPark includes two parts: the dock (which remains on the floor against the wall) and the stand-mounted module (which attaches to the underside of the stand).



Figure 1 PowerPark includes two main components: the stand-mounted module, and the dock.

Docking the stand


To dock the stand

- 1 Make sure the AC power cord is disconnected from the system stand.
- 2 Wheel the stand backward into the dock.

The stand module (on the underside of the stand) snaps into place as it connects with the dock.

- 3 Verify PowerPark is charging by confirming that the connection indicator on the stand base glows green.

To undock the stand

- 1 Make sure the AC power cord is disconnected from the system stand.
- 2 Wheel the stand away from the dock. The connection indicator turns off.

Troubleshooting

Table 1: Symptoms and solutions

Symptom	Solution
PowerPark does not charge the battery	▶ Make sure that the electrical outlet has power. ▶ Make sure that the dock circuit breaker is turned on. ▶ Reset the stand-module circuit breaker if tripped. ▶ Make sure that the stand module is fully connected to the dock.
Stand module does not connect to the dock	Make sure that the floor is clean, flat, and rigid. Avoid thick carpeting.


To turn the dock circuit breaker on or off

The circuit breaker on the PowerPark dock is typically set to the ON position but trips if there is a power event (for example, a power surge).

The circuit breaker is a means of disconnection from mains power.

- 1 Undock the stand.

- 2 Disconnect the power cord from the wall outlet, and pull the PowerPark dock away from the wall.
The circuit breaker is now accessible.
- 3 Press the circuit breaker to the desired position: **ON** (|) or **OFF** (O).



To reset the stand module circuit breakers

One or both of the circuit breakers on the stand module trip if there is a power event (for example, a power surge) and may need to be reset.

- 1 Locate the two breakers on the end of the power module under the stand base, near the stand height adjustment pedal.
- 2 Press in on each circuit breaker switch until it closes. See the following figure.
 - ▶ On a correctly closed circuit breaker, the switch is pressed in.
 - ▶ On a tripped circuit breaker, the switch is popped up and exposes a white plastic area.

Note If you cannot press down the switch and if no white area is exposed, the circuit breaker is closed.

Breaker closed

Breaker tripped

Cleaning and disinfecting

WARNINGS

- ▶ To avoid electrical shock, before cleaning, always disconnect the PowerPark dock from AC power.
- ▶ To avoid injury, always use protective eyewear and gloves when performing cleaning and disinfecting procedures.
- ▶ To avoid infection:
 - ▶ Ensure that the solution expiration date has not passed.
 - ▶ The level of disinfection required for a product is dictated by the type of bodily fluids or tissue it comes into contact with during use. Ensure that the solution strength and duration of contact are appropriate for the equipment. For information, see the disinfectant label instructions and the recommendations of the Association for Professionals in Infection Control and Epidemiology (APIC) and FDA.

Caution Use only approved cleaners or disinfectants on surfaces.

The exterior surfaces of the PowerPark stand module and the PowerPark dock should be kept clean. Whenever you clean the system stand, be sure to also clean the attached PowerPark stand module. Clean and disinfect the PowerPark dock periodically, or whenever it appears dirty or comes into contact with biomatter or bodily fluids. Refer to the cleaners and disinfection tool available at

www.sonosite.com/support/cleaners-disinfectants for a list of approved cleaners and disinfectants.

Caution Do not allow cleaning solution or disinfectant to come into contact with the PowerPark electrical connectors.

To clean the PowerPark stand module and dock

- 1 Disconnect the PowerPark components from any power source:
 - ▶ For the PowerPark stand module, turn off the ultrasound system and any accessory equipment. Disconnect the stand from the PowerPark dock or unplug the power cord.
 - ▶ For the PowerPark dock, disconnect it from its AC power.

2 Clean the stand module and PowerPark dock surfaces to remove any debris or bodily fluids. Use the following procedure:

- a** Use either a pre-moistened wipe or a soft cloth dampened with cleaner or disinfectant.
- b** Remove all gel, debris, and bodily fluids from the PowerPark components.

3 With a new wipe, clean all parts of the PowerPark components by wiping from the clean areas to the soiled areas. This method helps to avoid cross-contamination.

Follow the chemical manufacturer's instructions regarding temperature and wet contact time. Monitor the PowerPark parts for wet appearance, and reapply with a new wipe if no longer wet.

4 Allow the PowerPark components to air dry in a clean, well-ventilated space.

Safety

For maximum safety, observe all of the warnings and cautions in this section.

Electrical safety

Table 2: Electrical safety classification

Class I equipment	The PowerPark is classified as Class I equipment.
Non AP/APG	The PowerPark is not suitable for use in the presence of flammable anesthetics.

WARNINGS

- ▶ To avoid the risk of a shock hazard:
 - ▶ Do not plug the power cord of the PowerPark dock or PowerPark stand module into a multiple portable socket outlet (MPSO) or use an extension cord.
 - ▶ Make sure the AC power cord is disconnected from the PowerPark stand module whenever the stand is docked.
- ▶ To avoid the risk of electrical shock and possible fire:
 - ▶ Use only properly grounded equipment. Shock hazards exist if the power supply is not properly grounded. Grounding reliability can only be achieved when the equipment is connected to a receptacle marked "Hospital Only," "Hospital Grade," or the equivalent. The grounding wire must not be removed or circumvented.
 - ▶ Use only accessories and peripherals recommended by Fujifilm Sonosite, including the power supply. Use the PowerPark only within the voltage and current limits specified on its ratings label.
 - ▶ Inspect the power supply, plugs, cables, and power cords on a regular basis for damage.
 - ▶ Do not place the PowerPark dock in a wet environment.
- ▶ To avoid system malfunction and the risk of electrical shock, follow the PowerPark assembly instructions.

Equipment safety

Cautions

- ▶ Excessive bending or twisting of cables can cause a failure or intermittent operation.
- ▶ Improper cleaning or disinfecting of any part of the accessories and peripherals can cause permanent damage. Refer to the cleaners and disinfection tool available at www.sonosite.com/support/cleaners-disinfectants for a list of approved cleaners and disinfectants

Hazardous materials

WARNING Products and accessories may contain hazardous materials. When disposing of products and accessories, be environmentally responsible and meet federal and local regulations for disposing hazardous materials.

Electromagnetic compatibility (EMC)

See the Safety section of the ultrasound system user guide for information on electromagnetic compatibility with IEC 60601-1-2, including the Manufacturer's Declaration.

Electrostatic discharge

WARNINGS

- ▶ Unless following ESD precautionary procedures, all users and staff must be instructed not to connect or to touch (with body or hand-held tools) the pins of connectors that have the ESD Sensitive Devices symbol:

- ▶ If the symbol is on a border surrounding multiple connectors, the symbol pertains to all connectors within the border.
- ▶ ESD precautionary procedures include the following:
 - ▶ Receive training about ESD, including the following (at a minimum): an introduction to the physics of electrostatic charge, the voltage levels that can occur in normal practice, and the damage that can occur to electronic components if equipment is touched by an individual who is electrostatically charged.
 - ▶ Prevent the buildup of electrostatic charge; for example, use humidification, conductive floor coverings, non-synthetic clothing, ionizers, and minimizing insulating materials.
- ▶ Ground yourself to discharge any static charge. Use a wrist strap to bond yourself to earth ground.

Labeling symbols

The following symbols are used on the products, packaging, and containers.

Table 3: Standards labeling symbols

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Manufacturer	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.1.1	Indicates the medical device manufacturer, as EU directives 90/385/EEC, 93/42/EEC and 98/79/EC

Table 3: Standards labeling symbols

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Serial Number	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.1.7	Indicates the manufacturer’s serial number so that a specific medical device can be identified
	Catalog Number	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.1.6	Indicates the manufacturer’s catalogue number so that the medical device can be identified
	Consult instructions for use	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.4.3	Indicates the need for the user to consult the instructions for use
	Caution	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot for a variety of reasons, be presented on the medical device itself
	Alternating current	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5032	Indicates on the rating plate, that the equipment is suitable for alternating current only, in order to identify appropriate terminals

English

Deutsch

Español

Français

Italiano

Nederlands

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Table 3: Standards labeling symbols

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Fragile handle with care	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.1	Indicates a medical device that can be broken or damaged if not handled carefully
	Keep dry	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.4	Indicates a medical device that needs to be protected from moisture
	Medical Device	EU MDR	EU MDR Annex I, 23.2 (q)	Indicates that the device is a medical device
	Temperature limit	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.7	Indicates the temperature limits to which the medical device can be safely exposed
	Atmospheric pressure limitations	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.9	Indicates the range of atmospheric pressure to which the medical device can be safely exposed

Table 3: Standards labeling symbols

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	Humidity limitation	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.3.8	Indicates the range of humidity to which the medical device can be safely exposed
	Stacking limit by number	ISO 7000:2004 Graphical symbols for use on equipment	2403	Do not stack over n high, where n represents the number on the label.
	Corrugated recycle	—	—	Shipping box is made of corrugated cardboard and should be recycled accordingly
	RESY – Recycling Symbol	—	—	Paper recycle
	CE Marking	Council Directive 93/42/EEC	Article 17 Annex XII, Annex VII	Signifies European Technical Conformity
	Recycle: Electronic Equipment	BS EN 50419:2016 Marking of Electrical and Electronic Equipment in accordance with Directive 2012/19/EU (WEEE)	Annex IX	Do Not Throw in Trash
	Canadian Standard Association Certification Mark	—	—	CSA certification mark signifying that the product complies with the applicable CSA and ANSI/UL requirements and is authorized for use in Canada and the US

English

Deutsch

Español

Français

Italiano

Nederlands

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Table 3: Standards labeling symbols

Symbol	Title	Standards development organization (SDO)	Reference number	Description
	China Pollution Control (10)	ISO 7000:2014 Graphical symbols for use on equipment	1135	Pollution Control Logo. (Applies to all parts/products listed in the China RoHS disclosure table. May not appear on the exterior of some parts/products because of space limitations.)
	European community authorized representative	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labelling and information to be supplied	5.1.2	Indicates the Authorized representative in the European Community

Specifications

Dimensions

Approximate dimensions of the PowerPark dock are as follows:

- ▶ Width: 15.00 in (38.1 cm)
- ▶ Depth: 4.88 in. (12.4 cm)
- ▶ Length: 4.88 in. (12.4 cm)
- ▶ Weight (including cord): 11.3 lbs. (5.1 kg)

Environmental limits

Operating: 10 to 40°C; 15 to 95% R.H.; 700 hPa to 1060 hPa

Shipping/Storage: -35 to 65°C; 15 to 95% R.H.; 500 hPa to 1060 hPa

Electrical

PowerPark Dock Supply Input: 100–240 VAC, 50–60 Hz, 8–3.3 A

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