SonoSite iViz

User Guide Errata P21949-03

FUJIFILM
Value from Innovation

SonoSite

Manufacturer
FUJIFILM SonoSite, Inc.

21919 30th Drive SE
Bothell, WA 98021
USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative
Emergo Europe

Molenstraat 15
2513 BH, The Hague
The Netherlands

Australia Sponsor
FUJIFILM SonoSite Australasia Pty Ltd

114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution: | United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, iViz, and the SonoSite logo are registered and unregistered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions.

All other trademarks are the property of their respective owners.

Part Number: P21949-03

Publication Date: February 2016

Copyright © 2016 FUJIFILM SonoSite, Inc. All Rights reserved.

CE
0086

SonoSite iViz Ultrasound System

User Guide Errata

Getting started	1
Charging the battery	1
Configuring iViz	2
Connecting to a virtual private network (VPN)	2
Measurement References	3
Measurement accuracy	3
Cardiac references	4
Safety	5
iViz labeling symbols	5

The following content was either not included or incorrect in the system user guide; the revisions will be made in the next update.

Getting started

Charging the battery

To charge the battery using the iViz battery charger

- 1 Connect the line cord of the AC power adapter to a hospital-grade electrical outlet.
- 2 Connect the AC power adapter to the receptacle on the battery charger.
- 3 Remove the battery from the iViz system, and insert it in a slot on the battery charger. You can charge up to two batteries at a time.

Note

The optional iViz dual charging station allows you to charge up to four batteries at a time.

4 To determine battery charging status, check the lights on the battery charger:

Blinking green	Battery is charging.
Green	Battery is fully charged.
Amber/light red	Battery is very low.

If the light does not change to blinking or steady green after 3 hours, then try the following:

- ▶ Remove and re-insert the battery.
- ▶ Turn the battery charger off and on.

Red	Bad battery pack; recycle the battery.
-----	--

Caution

If the light on the battery charger remains amber/light red for more than 3 hours, even after multiple charging attempts, the battery may be faulty. Do not use the battery in the iViz system. Do not try to repair the battery. Remove the battery from the charger, and then contact FUJIFILM SonoSite or your representative to get a replacement. Recycle or dispose of faulty batteries in accordance with applicable laws.

Configuring iViz

Connecting to a virtual private network (VPN)

Note:

For security reasons, your iViz system must be secured with some form of password before you can connect to a VPN. For more information, see **Activating security settings**.

Before trying to connect to a VPN, you must gather the following information:

- ▶ The name of the VPN.
- ▶ The type of VPN, including any specialized settings.
- ▶ Your VPN credentials, including user name and password.

To connect to a VPN

- 1 From the Android **Settings** screen, under **Wireless & Networks**, tap **More...**
- 2 Tap **VPN**.
- 3 On the VPN screen, tap +.
- 4 Enter your VPN profile information, including the name of the VPN, the VPN type, and any additional settings.

5 Tap **Save**.

6 Log in to the VPN using your VPN user name and password.

To enable future VPN connection, check the **Save account information** check box.

7 To add another VPN, tap +.

Measurement References

Measurement accuracy

2D Measurement Accuracy and Range

2D Measurement Accuracy and Range	System Tolerance ^a	Accuracy By	Test Method ^b	Range (cm)
Axial Distance	< ±2% plus 1% of full scale	Acquisition	Phantom	0-26 cm
Lateral Distance	< ±2% plus 1% of full scale	Acquisition	Phantom	0-35 cm
Diagonal Distance	< ±2% plus 1% of full scale	Acquisition	Phantom	0-44 cm
Area ^c	< ±4% plus (2% of full scale/smallest dimension) * 100 plus 0.5%	Acquisition	Phantom	0.01-720 cm ²
Circumference ^d	< ±3% plus (1.4% of full scale/smallest dimension) * 100 plus 0.5%	Acquisition	Phantom	0.01-96 cm

^a Full scale for distance implies the maximum depth of the image.

^b An RMI 413a model phantom with 0.7 dB/cm MHz attenuation was used.

^c The area accuracy is defined using the following equation:

$$\% \text{ tolerance} = ((1 + \text{lateral error}) * (1 + \text{axial error}) - 1) * 100 + 0.5\%$$

^d The circumference accuracy is defined as the greater of the lateral or axial accuracy and by the following equation:

$$\% \text{ tolerance} = (2 (\text{maximum of 2 errors}) * 100) + 0.5\%$$

Cardiac references

Left Atrium/Aorta (LA/Ao)

Feigenbaum, H. Echocardiography. Philadelphia: Lea and Febiger, (1994), 206, Figure 4-49.

Left Ventricular End Volumes (Teichholz) in ml

Teichholz, L.E., T. Kreulen, M.V. Herman, et. al. "Problems in echocardiographic volume determinations: echocardiographic-angiographic correlations in the presence or absence of asynergy." *American Journal of Cardiology*, (1976), 37:7.

$$LVESV = (7.0 * LVDS^3)/(2.4 + LVDS)$$

where: LVESV = Left Ventricular End Systolic Volume
LVDS = Left Ventricular Dimension at Systole

$$LVEDV = (7.0 * LVDD^3)/(2.4 + LVDD)$$

where: LVEDV = Left Ventricular End Diastolic Volume (ml)
LVDD = Left Ventricular Dimension at Diastole (cm)

Left Ventricular Volume: Single Plane Method in ml

Schiller, N.B., P.M. Shah, M. Crawford, et.al. "Recommendations for Quantitation of the Left Ventricle by Two-Dimensional Echocardiography." *Journal of American Society of Echocardiography*. September-October 1989, 2:362.

$$V = \left(\frac{\pi}{4}\right) \sum_{i=1}^{20} a_i^2 \left(\frac{L}{n}\right)$$

$i=1$ to 20

where: V = Volume in ml
a = Diameter of chamber view slice i
n = Number of segments (n=20)
L = Length of the chamber view
i = Segment

Safety

iViz labeling symbols

Symbol	Definition
	Non-ionizing electromagnetic radiation. To indicate generally elevated, potentially hazardous, levels of non-ionizing radiation, or to indicate equipment or systems (for example, in the electrical medical area) that include RF transmitters or that intentionally apply electromagnetic energy for diagnosis or treatment.

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