SonoSite Edge

SonoSite S

SonoSite M-Turbo

User Guide Vet Supplement

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, the SonoSite logo, Edge, S Series and M-Turbo are trademarks and registered trademarks of FUJIFILM Sonosite, Inc. in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Patents: .

Part number: P13329-06

Publication date: November 2017

Copyright © 2017 FUJIFILM SonoSite, Inc. All rights reserved.

SonoSite Edge, SonoSite S Series, SonoSite M-Turbo User Guide Vet Supplement

Imaging	1
2D controls	1
Imaging modes and exams available by transducer	1
Measurements and Calculations (Vet Only)	2
Cardiac calculations	2
Volume calculations	2
Percent reduction calculations	2
Safety	3
Acoustic output tables	3

This user guide supplement describes the C11x transducers available for veterinary use. It supplements the existing C11x-transducer information in the *Edge User Guide*, *S Series User Guide*, and *M-Turbo Ultrasound System User Guide*.

Imaging

2D controls

Note

The Sector control is unavailable on the veterinary C11x transducer.

Imaging modes and exams available by transducer

The Abv exam type is available on the following ultrasound systems:

- ▶ SonoSite Edge (Vet)
- ▶ SonoSite S Series (Vet)
- ▶ SonoSite M-Turbo (Vet)

Table 1: Imaging modes and exams available

Transducer	Exam Type ¹	Imaging Mode				
		2D ² M Mode	CPD ³	Color ³	PW Doppler	CW Doppler
C11x	Crd	✓	—	✓	✓	—
	Abv	✓	✓	✓	✓	—

1. Exam type abbreviations are as follows: Crd = Cardiac, Abv = Abdomen Vet
2. The optimization settings for 2D are Res, Gen, and Pen.
3. The optimization settings for Color are low, medium, and high (flow sensitivity) with a range of PRF settings.

Measurements and Calculations (Vet Only)

Cardiac calculations

Transducer	Exam Type
C11x	Cardiac

Volume calculations

Transducer	Exam Type
C11x	Abv

Percent reduction calculations

Transducer	Exam Type
C11x	Abv

Safety

Acoustic output tables

Table 2: Transducer Model: C11x Operating Mode: Color

Index Label			M.I.	TIS			TIB	TIC	
				Scan	Non-scan		Non-scan		
					$A_{aprt} \leq 1$	$A_{aprt} > 1$			
Global Maximum Index Value			(a)	(a)	—	—	—	1.1	
Associated Acoustic Parameter	Pr.3	(MPa)	#						
	W_0	(mW)		#	—		—	39.4	
	min of [$W_{.3}(z_1), I_{TA.3}(z_1)$]	(mW)				—			
	z_1	(cm)				—			
	z_{bp}	(cm)				—			
	z_{sp}	(cm)	#				—		
	$d_{eq}(z_{sp})$	(cm)					—		
	f_c	(MHz)	#	#	—	—	—	4.37	
	Dim of A_{aprt}	X (cm)		#	—	—	—	—	1.12
Y (cm)			#	—	—	—	—	0.50	
Other Information	PD	(μ sec)	#						
	PRF	(Hz)	#						
	$p_r@PII_{max}$	(MPa)	#						
	$d_{eq}@PII_{max}$	(cm)					—		
	Focal Length	FL_x (cm)			#	—	—		4.29
		FL_y (cm)			#	—	—		4.40
$I_{PA.3}@MI_{max}$	(W/cm^2)	#							
Operating Control Condition	Control 1: Mode							Color	
	Control 2: Exam Type							Cardiac	
	Control 3: PRF							2033	
	Control 4: Optimization/Depth							Low/5.1	
	Control 5: Color Box Position/Size							Top/Short & Narrow	

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