SonoSite Edge
SonoSite S Series
C60xi Transducer

User Guide Supplement

Manufacturer

FUJIFILM SonoSite, Inc.

21919 30th Drive SE

Bothell, WA 98021 USA

T: 1-888-482-9449 or 1-425-951-1200

F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.

Joop Geesinkweg 140

1114 AB Amsterdam,

The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd

114 Old Pittwater Road

BROOKVALE, NSW, 2100

Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, the SonoSite logo, Edge, and S Series are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P19888-04

Publication date: June 2019

Copyright © 2019 FUJIFILM SonoSite, Inc. All rights reserved.

SonoSite Edge and SonoSite S Series C60xi User Guide Supplement

Introduction	1
Imaging	1
Needle visualization	1
Imaging modes and exams available by transducer	2
Measurements and calculations	4
Percent reduction calculations	4
Volume calculations	4
Volume flow calculations	5
Gynecology (Gyn) calculations	5
OB calculations	6
Cleaning and disinfecting	6
Compatible accessories and peripherals	7
Acoustic Output	7
Guidelines for reducing MI and TI	7
Output display	8
Acoustic output tables	9
Specifications	12
Supported transducers	12

Introduction

This user guide supplement provides information on the C60xi transducer, compatible with the SonoSite Edge and SonoSite S Series ultrasound systems.

Imaging

Needle visualization

About MBe

The MBe control is available with the C60xi transducer.

Imaging modes and exams available by transducer

Table 1: Transducer, exam type, and imaging mode (for the Edge ultrasound system)

Transducer	Exam type ^a	Imaging mode					
		2D ^b M Mode	MBe	CPD ^c	Color ^c	PW Doppler ^d	CW Doppler
C60xi	Abd	✓	—	✓	✓	✓	—
	Gyn	✓	—	✓	✓	✓	—
	Msk	✓	✓	✓	✓	✓	—
	Nrv	✓	✓	✓	✓	✓	—
	OB	✓	—	✓	✓	✓	—

^a Exam type abbreviations are as follows: Abd = Abdomen, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, IMT = Intima Media Thickness, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Vas = Vascular, Ven = Venous, Pros = Prostate.

^b The optimization settings for 2D are Res, Gen, and Pen.

^c The optimization settings for CPD and Color are low, medium, and high (flow sensitivity) with a range of PRF settings for Color depending on the setting selected.

^d For the cardiac exam type, PW TDI is also available.

Table 2: Transducer, exam type, and imaging mode (for the S Series ultrasound system)

S Series system	Transducer	Exam type ^a	Imaging mode					
			2D ^b M-Mode	MBe	CPD	Color	PW Doppler	CW Doppler
S-Cath	C60xi ^c	Abd	✓	—	✓	✓	✓	—
S-FAST	C60xi ^c	Abd	✓	—	✓	✓	✓	—
S-GYN	C60xi ^c	Gyn	✓	—	✓	✓	✓	—
		OB	✓	—	✓	✓	✓	—
S-ICU	C60xi ^c	Abd	✓	—	✓	✓	✓	—
S-MSK	C60xi ^c	Abd	✓	—	✓	✓	✓	—
		Msk	✓	✓	✓	✓	✓	—
		Nrv	✓	✓	✓	✓	✓	—
S-Nerve	C60xi ^c	Nrv	✓	✓	✓	✓	✓	—
S-Women's Health	C60xi ^c	Gyn	✓	—	✓	✓	✓	—
		OB	✓	—	✓	✓	✓	—

^a Exam type abbreviations are as follows: Abd = Abdomen, Gyn = Gynecology, Msk = Musculoskeletal, Nrv = Nerve, OB = Obstetrical

^b The optimization settings for 2D are Res, Gen, and Pen.

^c This transducer includes Tissue Harmonic Imaging.

Measurements and calculations

Percent reduction calculations

WARNING

To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate. To avoid misdiagnosis or harming the patient outcome, start a new patient information form before starting a new patient exam and performing calculations. Starting a new patient information form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

Transducer	Exam types
C60xi	Abdomen, Msk

Volume calculations

WARNING

To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate. To avoid misdiagnosis or harming the patient outcome, start a new patient information form before starting a new patient exam and performing calculations. Starting a new patient information form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

Transducer	Exam types
C60xi	Abdomen, Gyn, Msk, Nerve

Volume flow calculations

WARNING

To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate. To avoid misdiagnosis or harming the patient outcome, start a new patient information form before starting a new patient exam and performing calculations. Starting a new patient information form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

Transducer	Exam types
C60xi	Abdomen

Consider the following factors when performing volume flow measurements:

- ▶ The factors identified in the literature that affect the accuracy are as follows:
 - ▶ Difficulty ensuring uniform insonation of the vessel. The system is limited to the following sample volume sizes:
 - ▶ C60xi transducer: 2, 3, 5, 7, 10, 12 Gate Size (mm)

Gynecology (Gyn) calculations

WARNING

To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate. To avoid misdiagnosis or harming the patient outcome, start a new patient information form before starting a new patient exam and performing calculations. Starting a new patient information form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

Transducer	Exam types
C60xi	Gyn

OB calculations

WARNING

To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate. To avoid misdiagnosis or harming the patient outcome, start a new patient information form before starting a new patient exam and performing calculations. Starting a new patient information form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

Transducer	Exam types
C60xi	OB

Cleaning and disinfecting

Table 3: Approved cleaners

Product ^{a, b}	Compatible transducer
Sani-Cloth AF3	C60xi
PI-Spray II	C60xi

^a Refer to the manufacturer's instructions for concentration, temperature, and duration.

^b For a complete list of approved cleaners and disinfectants, refer to the cleaners and disinfectants tool available at www.sonosite.com/support/cleaners-disinfectants.

Table 4: Approved high-level compatible disinfectants

Disinfectant ^{a, b}	Compatible transducer
Cidex	C60xi
Cidex OPA	C60xi

^a Refer to the manufacturer's instructions for concentration, temperature, and duration.

^b For a complete list of approved cleaners and disinfectants, refer to the cleaners and disinfectants tool available at www.sonosite.com/support/cleaners-disinfectants.

Compatible accessories and peripherals

Table 5: Accessories and peripherals

Description	Maximum cable length
C60xi transducer ^{a,b}	5.5 ft/1.7 m
<p>^a For transducers, the maximum cable length is measured between the strain reliefs. The stated lengths do not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, and inside the transducer connector.</p> <p>^b This transducer supports a multi-angle biopsy guide.</p>	

Acoustic Output

Guidelines for reducing MI and TI

Table 6: MI

Transducer	Depth
C60xi	↑
<p>↓ Decrease or lower setting of parameter to reduce MI.</p> <p>↑ Increase or raise setting of parameter to reduce MI.</p>	

Table 7: TI (TIS, TIC, TIB)

Transducer	Color Power Doppler settings						PW settings
	Box width	Box height	Box depth	PRF	Depth	Optimize	
C60xi	↓	—	↑	↓	↑	—	↓ (PRF)
<p>↓ Decrease or lower setting of parameter to reduce TI.</p> <p>↑ Increase or raise setting of parameter to reduce TI.</p>							

Output display

Table 8: Models in which a TI or MI ≥ 1.0

Transducer	Index	2D/M Mode	CPD/Color	PW Doppler	CW Doppler
C60xi/5-2	MI	Yes	No	No	—
	TIC, TIB, or TIS	No	No	Yes	—

Table 9: Transducer surface temperature rise, external use (°C)

Test	C60xi
Still air	17.1
Simulated use	9.0

Acoustic output tables

Table 10: Transducer model: C60xi Operating mode: 2D

Index label			M.I.	TIS		TIB	TIC	
				Scan	Non-scan			Non-scan
					$A_{\text{aprt}} \leq 1$	$A_{\text{aprt}} > 1$		
Global maximum index value			1.0	(a)	—	—	(b)	
Associated acoustic parameter	Pr.3	(MPa)	1.69					
	W_0	(mW)	#	—		—	#	
	min of [$W_{.3}(z_1), I_{TA.3}(z_1)$]	(mW)			—			
	z_1	(cm)			—			
	z_{bp}	(cm)			—			
	z_{sp}	(cm)	4.7			—		
	$d_{eq}(z_{sp})$	(cm)				—		
	f_c	(MHz)	2.84	#	—	—	—	#
Dim of A_{aprt}	X (cm)		#	—	—	—	#	
	Y (cm)		#	—	—	—	#	
Other information	PD	(μsec)	0.579					
	PRF	(Hz)	5440					
	$p_r@PII_{\text{max}}$	(MPa)	2.679					
	$d_{eq}@PII_{\text{max}}$	(cm)				—		
	Focal Length	FL_x (cm)		#	—	—		#
		FL_y (cm)		#	—	—		#
$I_{PA.3}@MI_{\text{max}}$	(W/cm^2)	197.7						
Operating control conditions	Control 1: Exam type		Abd/ OB					
	Control 2: Optimization		Any					
	Control 3: Depth		11/ 13 cm					
	Control 4: THI		On					
	Control 5: MB (Multi Beam)		On					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data is not applicable for this transducer/mode.

Table 11: Transducer model: C60xi Operating mode: M Mode

Index label			M.I.	TIS		TIB	TIC	
				Scan	Non-scan			Non-scan
					$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global maximum index value			1.0	—	(a)	—	(a)	(b)
Associated acoustic parameter	Pr.3	(MPa)	1.62					
	W_0	(mW)		—	#		#	#
	min of $[W_{.3}(z_1), I_{TA.3}(z_1)]$	(mW)				—		
	z_1	(cm)				—		
	z_{bp}	(cm)				—		
	z_{sp}	(cm)	4.7				#	
	$d_{eq}(z_{sp})$	(cm)					#	
	f_c	(MHz)	2.85	—	#	—	#	#
	Dim of A_{aprt}	X (cm)		—	#	—	#	#
Y (cm)			—	#	—	#	#	
Other information	PD	(μ sec)	0.577					
	PRF	(Hz)	800					
	$p_r@PII_{max}$	(MPa)	2.576					
	$d_{eq}@PII_{max}$	(cm)					#	
	Focal Length	FL_x (cm)		—	#	—		#
		FL_y (cm)		—	#	—		#
$I_{PA.3}@MI_{max}$	(W/cm^2)	184.3						
Operating control conditions	Control 1: Exam type		Any					
	Control 2: Optimization		Pen					
	Control 3: Depth		7.8 cm					
	Control 4: MB (Multi Beam)		Off or On					

(a) This index is not required for this operating mode; value is <1.
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.
 # No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)
 — Data is not applicable for this transducer/mode.

Table 12: Transducer model: C60xi Operating mode: PW Doppler

Index label		M.I.	TIS			TIB	TIC	
			Scan	Non-scan		Non-scan		
				$A_{\text{aprt}} \leq 1$	$A_{\text{aprt}} > 1$			
Global maximum index value		(a)	—	(a)	—	3.1	(b)	
Associated acoustic parameter	Pr.3	(MPa) #						
	W_0	(mW)	—	#		85.64	#	
	min of $[W_{.3}(z_1), I_{\text{TA}.3}(z_1)]$	(mW)			—			
	z_1	(cm)			—			
	z_{bp}	(cm)			—			
	z_{sp}	(cm) #				1.255		
	$d_{\text{eq}}(z_{\text{sp}})$	(cm)				0.51		
	f_c	(MHz) #	—	#	—	2.233	#	
	Dim of A_{aprt}	X (cm)		—	#	—	0.6552	#
Y (cm)			—	#	—	1.3	#	
Other information	PD	(μsec) #						
	PRF	(Hz) #						
	$p_r @ PII_{\text{max}}$	(MPa) #						
	$d_{\text{eq}} @ PII_{\text{max}}$	(cm)				0.415		
	Focal Length	FL_x (cm)		—	#	—		#
		FL_y (cm)		—	#	—		#
$I_{\text{PA}.3} @ MI_{\text{max}}$	(W/cm^2) #							
Operating control conditions	Control 1: Exam type					Abd		
	Control 2: PRF					Any		
	Control 3: Sample volume					12 mm		
	Control 4: Sample volume position					Zone 1		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data is not applicable for this transducer/mode.

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Schallkopf	Untersuchungstyp ^a	Bildgebungsmodus					
		2D ^b M-Modus	MBe	CPD ^c	Farbe ^c	PW-Doppler ^d	CW-Doppler
C60xi	Abd	✓	–	✓	✓	✓	–
	Gyn	✓	–	✓	✓	✓	–
	Msk	✓	✓	✓	✓	✓	–
	Nrv	✓	✓	✓	✓	✓	–
	OB	✓	–	✓	✓	✓	–

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S Series-System	Schallkopf	Untersuchungstyp ^a	Bildgebungsmodus					
			2D ^b M-Modus	MBe	CPD	Farbe	PW-Doppler	CW-Doppler
S-Cath	C60xi ^c	Abd	✓	–	✓	✓	✓	–
S-FAST	C60xi ^c	Abd	✓	–	✓	✓	✓	–
S-GYN	C60xi ^c	Gyn	✓	–	✓	✓	✓	–
		OB	✓	–	✓	✓	✓	–
S-ICU	C60xi ^c	Abd	✓	–	✓	✓	✓	–
S-MSK	C60xi ^c	Abd	✓	–	✓	✓	✓	–
		Msk	✓	✓	✓	✓	✓	–
		Nrv	✓	✓	✓	✓	✓	–
S-Nerve	C60xi ^c	Nrv	✓	✓	✓	✓	✓	–
S-Women's Health	C60xi ^c	Gyn	✓	–	✓	✓	✓	–
		OB	✓	–	✓	✓	✓	–

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Transducer	Undersøgelsestype ^a	Billeddannelsestype					
		2D ^b M-Mode	MBe	CPD ^c	Farve ^c	PW-doppler ^d	CW-doppler
C60xi	Abd	✓	—	✓	✓	✓	—
	Gyn	✓	—	✓	✓	✓	—
	Msk	✓	✓	✓	✓	✓	—
	Nrv	✓	✓	✓	✓	✓	—
	OB	✓	—	✓	✓	✓	—

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S Series-system	Transducer	Under-søgelses-type ^a	Billeddannelsestype					
			2D ^b M-type	MBe	CPD	Color	PW-doppler	CW-doppler
S-Cath	C60xi ^c	Abd	✓	—	✓	✓	✓	—
S-HURTIG	C60xi ^c	Abd	✓	—	✓	✓	✓	—
S-GYN	C60xi ^c	Gyn	✓	—	✓	✓	✓	—
		OB	✓	—	✓	✓	✓	—
S-ICU	C60xi ^c	Abd	✓	—	✓	✓	✓	—
S-MSK	C60xi ^c	Abd	✓	—	✓	✓	✓	—
		Msk	✓	✓	✓	✓	✓	—
		Nrv	✓	✓	✓	✓	✓	—
S-Nerve	C60xi ^c	Nrv	✓	✓	✓	✓	✓	—
S-kvindens helbred	C60xi ^c	Gyn	✓	—	✓	✓	✓	—
		OB	✓	—	✓	✓	✓	—

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