SonoSite Edge
SonoSite S Series

User Guide Vet Supplement

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, the SonoSite logo, Edge, and S Series are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P13329-08

Publication date: July 2019

Copyright © 2019 FUJIFILM SonoSite, Inc. All rights reserved.

CE
2797

SonoSite Edge and SonoSite S Series User Guide Vet Supplement

Imaging	1
2D controls	1
Imaging modes and exams available by transducer	1
Measurements and calculations (Vet only)	2
Cardiac calculations	2
Volume calculations	2
Percent reduction calculations	2
Safety	3
Acoustic output tables	3

This user guide supplement describes the C11x transducers available for veterinary use. It supplements the existing C11x-transducer information in the *Edge Ultrasound System User Guide* and the *S Series Ultrasound System User Guide*.

Imaging

2D controls

Note

The Sector control is unavailable on the veterinary C11x transducer.

Imaging modes and exams available by transducer

The Abv exam type is available on the following ultrasound systems:

- ▶ SonoSite Edge (Vet)
- ▶ SonoSite S Series (Vet)

Table 1: Imaging modes and exams available

Transducer	Exam type ^a	Imaging mode				
		2D ^b M Mode	CPD ^c	Color ^c	PW Doppler	CW Doppler
C11x	Crd	✓	—	✓	✓	—
	Abv	✓	✓	✓	✓	—

^aExam type abbreviations are as follows: Crd = Cardiac, Abv = Abdomen Vet

^bThe optimization settings for 2D are Res, Gen, and Pen.

^cThe optimization settings for Color are low, medium, and high (flow sensitivity) with a range of PRF settings.

Measurements and calculations (Vet only)

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