SonoSite Edge
SonoSite S Series

User Guide Supplement

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, Edge, and S Series are trademarks and registered trademarks of FUJIFILM Sonosite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P15653-08

Publication date: June 2019

Copyright © 2019 FUJIFILM SonoSite, Inc. All rights reserved.

SonoSite Edge and SonoSite S Series User Guide Supplement

Introduction	1
Getting started	1
Touchpad and cursor	1
Using AC power	1
Imaging	2
Imaging modes and exams available by transducer	2
Measurements and calculations	2
Volume calculations	2
Measurement publications and terminology (Cardiac)	2
Safety	3
Electrical safety	3
Labeling symbols	4
Acoustic output	4
Guidelines for reducing MI and TI	4
Output display	5
Acoustic output tables	6

Introduction

This user guide supplement provides information on the C8x transducer, compatible with the SonoSite Edge and SonoSite S Series ultrasound systems. It also provides supplemental information regarding safety, touchpad care, and calculation references.

Getting started

Touchpad and cursor

Make sure to keep the touchpad dry while in use. Moisture on the touchpad can cause the cursor to respond erratically.

Using AC power

When using AC power, position the system to allow easy access to disconnect it.

To separate the system (and any connected equipment) from a supply mains

Note

Disconnecting only the DC power cable from the system or (for Edge) dock does not separate the system from the supply mains.

- ❖ Disconnect the AC power cord from the power supply or (alternatively, if using a stand) from the AC adapter on the stand base.

Imaging

Imaging modes and exams available by transducer

Table 1: Transducer, exam type, and imaging mode

Transducer	Exam type	Imaging mode				
		2D ^a M Mode	CPD ^b	Color ^b	PW Doppler	CW Doppler
C8x	Prostate	✓	✓	✓	✓	—

^a The optimization settings for 2D are Res, Gen, and Pen.

^b The optimization settings for Color are low, medium, and high (flow sensitivity) with a range of PRF settings.

Measurements and calculations

Volume calculations

Transducer	Exam type
C8x	Prostate

Measurement publications and terminology (Cardiac)

The following are additional references for the Cardiac exam type.

Left Ventricular Mass in gm for 2D

Schiller, N.B., P.M. Shah, M. Crawford, et. al. "Recommendations for Quantification of the Left Ventricle by Two-Dimensional Echocardiography." Journal of American Society of Echocardiography. September-October 1998, 2:364.

$$\text{LV Mass} = 1.05 * \{[(5/6) * A1 * (a + d + t)] - [(5/6) * A2 * (a + d)]\}$$

where:

- A1 = Short axis area, diastole (Epi)
- A2 = Short axis area, diastole (Endo)
- a = Long or semi major axis
- d = Truncated semi major axis from the widest short axis diameter to mitral annulus plane
- t = Myocardial thickness

Left Ventricular Mass in gm for M Mode

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 2nd Edition, Boston: Little, Brown and Company, (1999), 39.

$$\text{LV Mass} = 1.04 [(LVID + PWT + IVST)^3 - LVID^3] * 0.8 + 0.6$$

where:

- LVID = Internal Dimension
- PWT = Posterior Wall Thickness
- IVST = Interventricular Septal Thickness
- 1.04 = Specific gravity of the myocardium
- 0.8 = Correction factor

Safety

Electrical safety

WARNING

To avoid the risk of electrical shock, do not allow any part of the system (including the bar code scanner, external mouse, power supply, power supply connector, external keyboard, and so on), except for the transducer or ECG leads, to touch the patient.

Labeling symbols

Symbol	Title	Standards development organization	Reference number	Description
	Consult instructions for use	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.4.3	Indicates the need for the user to consult the instructions for use
	Manufacturer	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.1.1	Indicates the medical device manufacturer, as EU directives 90/385/EEC, 93/42/EEC and 98/79/EC
	European community authorized representative	ISO 15223-1 Medical devices – symbols to be used with medical device labels, labelling and information to be supplied.	5.1.2	Indicates the Authorized representative in the European Community

Acoustic output

Guidelines for reducing MI and TI

Table 2: MI

Transducer	Depth
C8x	↑
↓ Decrease or lower setting of parameter to reduce MI. ↑ Increase or raise setting of parameter to reduce MI.	

Table 3: TI (TIS, TIC, TIB)

Transducer	CPD settings						PW settings
	Box width	Box height	Box depth	PRF	Depth	Optimize	
C8x	↓				↑		↓ (Depth)

↓ Decrease or lower setting of parameter to reduce TI.
 ↑ Increase or raise setting of parameter to reduce TI.

Output display

Table 4: TI or MI ≥ 1.0

Transducer	Index	2D/M Mode	CPD/Color	PW Doppler	CW Doppler
C8x	MI	Yes	Yes	Yes	—
	TIC, TIB, or TIS	No	No	Yes	—

Table 5: Transducer surface temperature rise, internal use (°C)

Test	C8x
Still air	11.3
Simulated use	5.5

Acoustic output tables

Table 6: Transducer model: C8x Operating mode: 2D

Index label			M.I.	TIS			TIB	TIC
				Scan	Non-scan		Non-scan	
					$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global maximum index value			1.1	(a)	—	—	—	(b)
Associated acoustic parameter	$P_{r.3}$	(MPa)	2.48					
	W_0	(mW)		#	—		—	#
	min of [$W_{.3}(z_1), I_{TA.3}(z_1)$]	(mW)				—		
	z_1	(cm)				—		
	z_{bp}	(cm)				—		
	z_{sp}	(cm)	1.2				—	
	$d_{eq}(z_{sp})$	(cm)					—	
	f_c	(MHz)	5.53	#	—	—	—	#
	Dim of A_{aprt}	X (cm)		#	—	—	—	#
Y (cm)			#	—	—	—	#	
Other information	PD	(μ sec)	0.239					
	PRF	(Hz)	9524					
	$p_r@PII_{max}$	(MPa)	3.11					
	$d_{eq}@PII_{max}$	(cm)					—	
	Focal Length	FL_x (cm)		#	—	—		#
		FL_y (cm)		#	—	—		#
$I_{PA.3}@MI_{max}$	(W/cm^2)	264						
Operating control conditions	Control 1: Exam type		Pro					
	Control 2: Optimization		Pen					
	Control 3: Depth		2.5 - 3.2					
	Control 4: MB		Off					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data is not applicable for this transducer/mode.

Table 7: Transducer model: C8x Operating mode: M Mode

Index label			M.I.	TIS		TIB	TIC	
				Scan	Non-scan			Non-scan
					$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global Mmaximum index value			1.4	—	(a)	—	(a)	(b)
Associated acoustic parameter	$P_{r,3}$	(MPa)	3.16					
	W_0	(mW)		—	#		#	#
	min of [$W_{,3}(z_1), I_{TA,3}(z_1)$]	(mW)				—		
	z_1	(cm)				—		
	z_{bp}	(cm)				—		
	z_{sp}	(cm)	1.1				#	
	$d_{eq}(z_{sp})$	(cm)					#	
	f_c	(MHz)	5.07	—	#	—	#	#
	Dim of A_{aprt}	X (cm)		—	#	—	#	#
Y (cm)			—	#	—	#	#	
Other Information	PD	(μ sec)	0.427					
	PRF	(Hz)	800					
	$p_r@PII_{max}$	(MPa)	3.83					
	$d_{eq}@PII_{max}$	(cm)					#	
	Focal Length	FL_x (cm)		—	#	—		#
		FL_y (cm)		—	#	—		#
$I_{PA,3}@MI_{max}$	(W/cm^2)	482						
Operating control conditions	Control 1: Exam type		Pro					
	Control 2: Optimization		Pen					
	Control 3: Depth		4.2					
<p>(a) This index is not required for this operating mode; value is <1.</p> <p>(b) This transducer is not intended for transcranial or neonatal cephalic uses.</p> <p># No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)</p> <p>— Data is not applicable for this transducer/mode.</p>								

Table 8: Transducer model: C8x Operating mode: Color

Index label			M.I.	TIS			TIB	TIC
				Scan	Non-scan		Non-scan	
					$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global maximum index value			1.4	(a)	—	—	—	(b)
Associated acoustic parameter	$P_{r.3}$	(MPa)	3.18					
	W_0	(mW)		#	—		—	#
	min of [$W_{.3}(z_1), I_{TA.3}(z_1)$]	(mW)				—		
	z_1	(cm)				—		
	z_{bp}	(cm)				—		
	z_{sp}	(cm)	0.8				—	
	$d_{eq}(z_{sp})$	(cm)					—	
	f_c	(MHz)	4.82	#	—	—	—	#
	Dim of A_{aprt}	X (cm)		#	—	—	—	#
Y (cm)			#	—	—	—	#	
Other information	PD	(μ sec)	0.694					
	PRF	(Hz)	2548					
	$p_r@PII_{max}$	(MPa)	3.63					
	$d_{eq}@PII_{max}$	(cm)					—	
	Focal Length	FL_x (cm)		#	—	—		#
		FL_y (cm)		#	—	—		#
$I_{PA.3}@MI_{max}$	(W/cm^2)	555						
Operating control conditions	Control 1: Exam type		Pro					
	Control 2: Mode		CVD					
	Control 3: 2D optimization/depth		Pen / 1.5 - 1.9					
	Control 4: Color optimization/PRF		High / Any					
	Control 5: Color box position/size		Short & narrow / Any					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data is not applicable for this transducer/mode.

Table 9: Transducer model: C8x Operating mode: PW Doppler

Index label			M.I.	TIS			TIB	TIC
				Scan	Non-scan		Non-scan	
					$A_{aprt} \leq 1$	$A_{aprt} > 1$		
Global maximum index value			1.2	—	(a)	—	2.0	(b)
Associated acoustic parameter	$P_{r,3}$	(MPa)	2.59					
	W_0	(mW)		—	#		36.0	#
	min of [$W_{.3}(z_1), I_{TA,3}(z_1)$]	(mW)				—		
	z_1	(cm)				—		
	z_{bp}	(cm)				—		
	z_{sp}	(cm)	1.1				1.10	
	$d_{eq}(z_{sp})$	(cm)					0.28	
	f_c	(MHz)	4.79	—	#	—	4.79	#
	Dim of A_{aprt}	X (cm)		—	#	—	1.12	#
Y (cm)			—	#	—	0.40	#	
Other information	PD	(μ sec)	1.131					
	PRF	(Hz)	1008					
	$p_r@PII_{max}$	(MPa)	3.10					
	$d_{eq}@PII_{max}$	(cm)					0.28	
	Focal Length	FL_x (cm)		—	#	—		#
		FL_y (cm)		—	#	—		#
$I_{PA,3}@MI_{max}$	(W/cm^2)	296						
Operating control conditions	Control 1: Exam type		Pro				Pro	
	Control 2: Sample volume size		1 mm				1 mm	
	Control 3: Sample volume position		Zone 5				Zone 5	
	Control 4: PRF		1008				3125	

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data is reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data is not applicable for this transducer/mode.

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