TEEx, TEEp, and TEExi Transducer

FUJIFILM
Value from Innovation

SonoSite

Bite Hole Indicator User Guide

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or 1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution | Federal (United States) law restricts this device to sale by or on the order of a physician.

SonoSite, and the SonoSite logo are registered and unregistered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Patents: US 6,371,918, CA 2,373,065, DE 60021552.0, FR 1175173 and GB 1175173.

Part Number: P20995-03

Publication date: November 2017

Copyright © 2017 FUJIFILM SonoSite, Inc. All rights reserved.

TEEx, TEEExp, and TEExi

Bite Hole Indicator User Guide

Introduction	1
Document conventions	1
Getting help	1
Overview	2
Transducer components	2
Testing for bite holes	3
Preparing for bite-hole testing	3
Step 1 — Testing the transducer array (TEEx, TEExi, and TEEExp)	4
Step 2 — Testing the endoscopic shaft (TEExi and TEEExp only)	7
Step 3 — Cleaning the transducer (TEEx, TEExi, and TEEExp)	9
If no bite holes are detected	10
If the transducer fails the test	10

Introduction

It is important that you test the TEEx, TEExi, or TEEExp transducer for punctures and tears in the waterproof coating before every procedure. If any damage, irregularity, substandard functioning, or unsafe condition is observed or suspected, do not use the transducer. Call FUJIFILM SonoSite or your local representative.

Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.

Getting help

In addition to this document, you can get help with this product by contacting FUJIFILM SonoSite Technical Support:

Phone
(U.S. or Canada)


877-657-8118

Phone (outside U.S. or Canada)	425-951-1330, or call your local representative
Fax	425-951-6700
Email	service@sonosite.com
Web	www.sonosite.com
Europe Service Center	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia Service Center	+65 6380-5581

Overview

Transducer components

A basic understanding of the TEE transducer components is important to successfully performing a bite hole detection test. In this guide, you will find references to the following components:

The endoscope shaft has no electrically conducting surfaces and is covered with a layer of material that permits neither fluids nor electricity to pass through it. Electrical safety is maintained for the transducer by keeping this material intact. Punctures in this material, such as those resulting from bites or improper handling, can result in fluids entering the endoscopic shaft and the patient being exposed to electrical current. Each transducer is tested for electrical isolation and leakage current before it is shipped to a customer.

WARNING

To avoid injury to the patient, do not use the transducer if the insulating material has been punctured or otherwise compromised.

Testing for bite holes

FUJIFILM SonoSite recommends testing the endoscopic shaft and the transducer array for bite holes after every procedure. You should perform bite hole detection after cleaning the transducer, but before disinfecting it, so you can identify punctures before exposing them to disinfectants which might leak inside and damage the transducer.

WARNING

Take all necessary precautions to avoid contact with non-insulated parts that have applied voltage.

Preparing for bite-hole testing

Test both the scan head and the endoscopic shaft separately for bite holes. You can use the same equipment for both tests; however, the tester connection points are different for each test.

Assemble the following items for this 2-step test:

- ▶ Water bath with a 5% saline solution (50g NaCl/1 liter water)
- ▶ A conductor plate made of copper or aluminum with an area of at least 25 cm² (one is included with bite hole indicator kit)
- ▶ Bite hole indicator with leads
- ▶ Puncture test tool (included only with the TEExi and TEEexp transducers)

Step 1 – Testing the transducer array (TEEx, TEExi, and TEEexp)

To test the transducer array for bite holes or current leakage


Figure 1 Transducer array test setup

1 After securing the controller so it does not fall into the water bath, immerse the scan head and endoscopic shaft in liquid to above the 40 cm mark, but below the 90 cm mark.

WARNING | Do not immerse the controller, cable, or connector in any fluid.

Caution | Do not allow the endoscopic shaft to contact the conductor plate. Doing so could result in an inaccurate test result.

2 Connect the Bite hole indicator to the transducer connector and the conductor plate:

- a** Connect the black lead to the conductor plate in the salt-water bath as shown in **Figure 2**.

Note

Use the lead attached to the conductor plate. Do not submerge the indicator clip.



Figure 2 Conductor plate attachment

- b** Identify the model of TEE transducer you are testing by checking the part number label on the connector. The top line of the label will include the model name.

This information will be important in the next step.

- c Connect the red leads to the knurled socket cap screws on the transducer connector as shown in **Figure 3**. The configuration of the leads depends on the model of transducer you have. See step 2b to determine your transducer model.


Figure 3 Transducer connector attachment

3 On the Bite hole indicator, press 

4 Read the test result:

- ▶ **Leads** - The leads are not connected correctly. Check the connections, and test again.
- ▶ **Fail** - There was a bite hole detected in the transducer array. Stop the bite hole test procedure. Do not use the transducer. For remediation steps, see **"If the transducer fails the test"** on page 10.
- ▶ **Pass** - No bite holes were found.
- ▶ **All** - If all of the lights are illuminated, the battery is low. Replace the battery.

If you are testing the TEEx transducer, proceed directly to **"Step 3 – Cleaning the transducer (TEEx, TEExi, and TEExp)"** on page 9.

If you are testing the TEExi or TEExp transducer, continue with **"Step 2 – Testing the endoscopic shaft (TEExi and TEExp only)"** on page 7.

Step 2 – Testing the endoscopic shaft (TEExi and TEEp only)

To test the endoscopic shaft for bite holes or current leakage


Figure 4 Endoscopic shaft test setup

- 1 Insert the puncture test tool under the scanplane controller as shown in [Figure 5](#). The puncture test tool is included with the TEExi and TEEp transducer.


Figure 5 Insertion of the puncture test tool.

- 2 Immerse the scan head and endoscopic shaft in liquid to above the 40 cm mark, but below the 90 cm mark.

WARNING

Do not immerse the controller, the cord, or the connector in any fluid.

Caution

Do not allow the endoscopic shaft to contact the conductor plate. Doing so could result in an inaccurate test result.

3 Connect the Bite hole indicator to the puncture test tool and the conductor plate:

a Connect the black lead to the conductor plate in the salt-water bath as shown in **Figure 6**.

Note

Use the lead attached to the conductor plate. Do not submerge the indicator clip.


Figure 6 Conductor plate attachment

b Connect the red leads to the puncture test tool as shown in **Figure 7**.

Note

Each clip should be attached separately to the puncture test tool.


Figure 7 Puncture test tool attachment

4 Press 

5 Read the test result:

- ▶ **Leads** - The leads are not connected correctly. Check the connections, and test again.
- ▶ **Fail** - There was a bite hole detected in the endoscopic shaft. Stop the bite hole test procedure. Do not use the transducer. For remediation steps, see **"If the transducer fails the test"** on page 10.
- ▶ **Pass** - No bite holes were found.
- ▶ **All** - If all of the lights are illuminated, the battery is low. Replace the battery.

Step 3 – Cleaning the transducer (TEEx, TEExi, and TEEp)

To clean the transducer after the test

- 1 Rinse the endoscopic shaft for one minute with a large volume of clean, lukewarm water to remove residual saline solution (for example, eight liters).
- 2 Dry the endoscopic shaft with a clean, non-linting towel.

If no bite holes are detected

To identify the transducer as safe, you should include a sticker or certificate that travels with the transducer that includes the date of the test, the name or other identification of the tester, and the outcome of the test. If the bite hole test was performed as part of the cleaning process, continue to clean and disinfect the transducer.

If the transducer fails the test

Do not use the transducer. Do not connect the transducer to an ultrasound system. Contact FUJIFILM SonoSite for repair.

To identify the transducer as unsafe to use, you should include a sticker or certificate that travels with the transducer that includes the date of the test, the name or other identification of the tester, and the outcome of the test.

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