SonoSite Edge II

SonoSite SII

User Guide Errata

Manufacturer

FUJIFILM SonoSite, Inc.
21919 30th Drive SE
Bothell, WA 98021 USA
T: 1-888-482-9449 or
1-425-951-1200
F: 1-425-951-1201

EC Authorized Representative

FUJIFILM SonoSite B.V.
Joop Geesinkweg 140
1114 AB Amsterdam,
The Netherlands

Australia Sponsor

FUJIFILM SonoSite Australasia Pty Ltd
114 Old Pittwater Road
BROOKVALE, NSW, 2100
Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, the SonoSite logo, SonoSite Edge, and SonoSite SII are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P27224-02

Publication date: June 2019

Copyright © 2019 FUJIFILM SonoSite, Inc. All rights reserved.

SonoSite Edge II and SonoSite SII User Guide

Errata

Introduction	2
Document conventions	2
Getting Help	2
Getting Help (SonoSite Edge II)	3
Getting Started	3
Intended uses	3
System Setup	3
Connectivity settings (SonoSite SII)	3
Connectivity setup (SonoSite Edge II)	3
Network Status settings	4
Imaging	4
Imaging modes and exams available by transducer (SonoSite SII)	4
Measurements and Calculations (SonoSite SII)	4
General calculations	4
Measurement References (SonoSite SII)	5
Measurement accuracy	6
Measurement publications and terminology	6
Cleaning and Disinfecting	7
Safety	7
Clinical safety	7
Electromagnetic compatibility	7
Wireless transmission	9
Wireless USB Dongle (Panda)	9
Wireless and Security Module (Laird)	9
Labeling symbols	20
Specifications	27
Supported transducers (SonoSite SII)	27
Standards	28
Electromechanical safety standards	28
Acoustic output	29
ALARA principle	29
Applying the ALARA principle	29
Direct controls	29
Output display	30
Related guidance documents	30
Acoustic output measurement	30
Terms used in the acoustic output tables	89
Glossary (SonoSite Edge II)	90

Introduction

Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ❖.

For a description of labeling symbols that appear on the product, see "Labeling Symbols" in the user guide.

Getting Help

For technical support, please contact FUJIFILM SonoSite as follows:

Phone (U.S. or Canada)	+1-877-657-8118
Phone (outside U.S. or Canada)	+1-425-951-1330, or call your local representative
Fax	+1-425-951-6700
Email	ffss-service@fujifilm.com
Web	www.sonosite.com
Europe Service Center	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia Service Center	+65 6380-5581

Printed in the U.S.

Getting Help (SonoSite Edge II)

The following has been corrected in the SonoSite Edge II system user guide; the revision will be made in the next update.

Email

ffss-service@fujifilm.com

Getting Started

The following content was either not included or incorrect in the SonoSite Edge II and SonoSite SII system user guides; the revisions will be made in the next update.

Intended uses

Prostate imaging applications

You can assess the prostate and surrounding anatomical structures for the presence or absence of pathology.

Superficial imaging applications

You can assess the breast, thyroid, testicle, lymph nodes, hernias, musculoskeletal structures, soft tissue structures, spine, ophthalmic structures, and surrounding anatomical structures for the presence or absence of pathology. You can use the system to provide ultrasound guidance for biopsy and drainage procedures, vascular line placement, and peripheral nerve blocks.

System Setup

Connectivity settings (SonoSite SII)

All references to PDAS should be changed to SiteLink in the SonoSite SII system user guide; the revision will be made in the next update.

Connectivity setup (SonoSite Edge II)

The following reference has been updated in the SonoSite Edge II system user guide; the revision will be made in the next update.

To enable wireless connection

- ❖ See *Setting Up a Network Connection*.

Network Status settings

If your Network Status screen displays a failed wireless device message, your network password may be expired. Make sure that you have an updated network password before connecting your wireless device.

Imaging

The C8x transducer is needle guide-capable on both the SonoSite Edge II and SonoSite SII systems.

Imaging modes and exams available by transducer (SonoSite SII)

The following footnotes were missing from **Table 4-5. Imaging modes and exams available by transducer** in the SonoSite SII system user guide; the revision will be made in the next update.

^aExam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Pro = Prostate, SmP = Small Parts, Spn = Spine, Sup = Superficial, Ven = Venous.

^bThe optimization settings for 2D are Res, Gen, and Pen.

^cThe optimization settings for CPD and Color are low, medium, and high (flow sensitivity) with a range of PRF settings for Color depending on the setting selected.

Measurements and Calculations (SonoSite SII)

General calculations

Volume calculation

- WARNINGS**
- ▶ To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate.
 - ▶ To avoid misdiagnosis or harming the patient outcome, start a new patient form before starting a new patient exam and performing calculations. Starting a new patient form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

The volume calculation involves three 2D distance measurements: D1, D2, and D3. After all measurements are saved, the result appears on-screen and in the patient report.

The volume calculation is available in the following exam types: Abdomen, Arterial, Breast, Gynecological, Musculoskeletal, Nerve, Small Parts, Venous, and Superficial.

To calculate volume

Do the following for each image you need to measure:

- 1 On a frozen 2D image, tap **Calcs**.
- 2 Do the following for each measurement you need to take:
 - a From the calculations menu, under **Volume**, select the measurement name.
If **Volume** is not available in a Gyn exam, select **Gyn** and then select **Volume**.
 - b Using the touchpad or the touchscreen, position the calipers
 - c Tap **Save Calc** to save the calculation.
A check mark appears next to the saved measurement.
- 3 To save a picture of the finished calculation, tap .
- 4 Tap **Back** to exit the calculation.

Measurement References (SonoSite SII)

The following information was not included in the SonoSite SII system user guide; the revision will be made in the next update.

Measurement accuracy

Table 1: M Mode Measurement and Calculation Accuracy and Range

M Mode Measurement Accuracy and Range	System Tolerance	Accuracy By	Test Method	Range
Distance	< +/- 2% plus 1% of full scale ^a	Acquisition	Phantom ^b	0–26 cm
Time	< +/- 2% plus 1% of full scale ^c	Acquisition	Phantom ^d	0.01–10 sec
Heart Rate	< +/- 2% plus (Full Scale ^c * Heart Rate/100)%	Acquisition	Phantom ^d	5–923 bpm

^aFull scale for distance implies the maximum depth of the image.

^bAn RMI 413a model phantom with 0.7 dB/cm MHz attenuation was used.

^cFull scale for time implies the total time displayed on the scrolling graphic image.

^dFUJIFILM SonoSite special test equipment was used.

Measurement publications and terminology

General references

Hip Angle/d:D Ratio

Morin, C., Harcke, H., MacEwen, G. "The Infant Hip: Real-Time US Assessment of Acetabular Development." *Radiology* 177: p.673-677, December 1985.

Percent Area Reduction

Zwiebel W.J., J.A. Zagzebski, A.B. Crummy, et al. "Correlation of peak Doppler frequency with lumen narrowing in carotid stenosis." *Stroke*, 3: (1982), p.386-391.

$$\% \text{ Area Reduction} = [1 - A2(\text{cm}^2)/A1(\text{cm}^2)] * 100$$

where: A1 = original area of the vessel in square cm
A2 = reduced area of the vessel in square cm

Percent Diameter Reduction

Handa, Nobuo et al., "Echo-Doppler Velocimeter in the Diagnosis of Hypertensive Patients: The Renal Artery Doppler Technique," *Ultrasound in Medicine and Biology*, 12:12 (1986), p.945-952.

$$\% \text{ Diameter Reduction} = [1 - D2(\text{cm})/D1(\text{cm})] * 100$$

where: D1 = original diameter of the vessel in cm
D2 = reduced diameter of the vessel in cm

Cleaning and Disinfecting

The following website was corrected in the SonoSite Edge II and SonoSite SII user guides; the revision will be made in the next update.

www.sonosite.com/products/transducers

The following footnote was updated in the cleaning and disinfecting tables.

Refer to the cleaners and disinfection tool available at www.sonosite.com/support/cleaners-disinfectants for a more complete list of approved cleaners and disinfectants.

Safety

Clinical safety

The following warning has been updated in the SonoSite Edge II and SonoSite SII user guides; the revision will be made in the next update.

WARNING

FUJIFILM SonoSite does not recommend the use of high-frequency (HF) electromedical devices in proximity to its systems. FUJIFILM SonoSite equipment has not been validated for use with HF electrosurgical devices or procedures. Use of HF electrosurgical devices in proximity to its systems may lead to abnormal system behavior or shutdown of the system.

To avoid the risk of a burn hazard, do not use the transducer with HF surgical equipment. Such a hazard may occur in the event of a defect in the HF surgical neutral electrode connection.

Electromagnetic compatibility

The ultrasound system has been tested and found to comply with the electromagnetic compatibility (EMC) limits for medical devices to IEC 60601-1-2:2007 and IEC 60601-1-2:2014. The ultrasound system is suitable for use in the professional healthcare facility environment. Active HF surgical equipment causes high electromagnetic disturbances which may interfere with the ultrasound system operation. The ultrasound system should not be operated inside an RF-shielded room where magnetic resonance imagery is present,

because it produces high electromagnetic disturbances which may interfere with the ultrasound system operation. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

Cautions

- ▶ Medical electrical equipment requires special precautions regarding EMC and must be installed and operated according to these instructions. Portable RF communications equipment (including peripherals, such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the ultrasound system, including cables specified by FUJIFILM SonoSite. Portable and mobile RF communications equipment can affect the ultrasound system. Electromagnetic interference (EMI) from other equipment or interference sources could result in performance disruption of the ultrasound system. Evidence of disruption may include image degradation or distortion, erratic readings, equipment ceasing to operate, or other incorrect functioning. If this occurs, survey the site to determine the source of disruption, and take the following actions to eliminate the source(s).
 - ▶ Turn equipment in the vicinity off and on to isolate disruptive equipment.
 - ▶ Relocate or re-orient interfering equipment.
 - ▶ Increase distance between interfering equipment and your ultrasound system.
 - ▶ Manage use of frequencies close to ultrasound system frequencies.
 - ▶ Remove devices that are highly susceptible to EMI.
 - ▶ Lower power from internal sources within facility control (such as paging systems).
 - ▶ Label devices susceptible to EMI.
 - ▶ Educate clinical staff to recognize potential EMI-related problems.
 - ▶ Eliminate or reduce EMI with technical solutions (such as shielding).
 - ▶ Restrict use of personal communicators (cell phones, computers) in areas with devices susceptible to EMI.
 - ▶ Share relevant EMI information with others, particularly when evaluating new equipment purchases which may generate EMI.
 - ▶ Purchase medical devices that comply with IEC 60601-1-2 EMC Standards.
- ▶ Do not stack other equipment on the ultrasound system or use other equipment in close proximity and adjacent to the ultrasound system. If stacking or using other equipment in close proximity is unavoidable, then you must observe the system to verify normal operation.

Note

The emissions characteristics of the SonoSite Edge II and SII ultrasound systems make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), the ultrasound system might not offer adequate protection to radio-frequency communication services. It may be necessary to take mitigation measures, such as relocating or reorienting the equipment.

Wireless transmission

The SonoSite Edge II and SII ultrasound systems implement two wireless solutions.

- ▶ Wireless USB Dongle (Panda) is a small wireless adapter that plugs into the USB port on the:
 - ▶ Right side of the Edge II ultrasound system.
 - ▶ Upper back of the SII ultrasound system.
- ▶ Wireless and Security Module (Laird) is a module that mounts on the:
 - ▶ Lid of the Edge II ultrasound system and then plugs into the system with a right-angled USB cord
 - ▶ Transducer holder arm of the SII ultrasound system and then plugs into the system with a 12 inch USB cord

Refer to the information below for the transmission information for each one.

Wireless USB Dongle (Panda)

The Wireless USB Dongle uses the Industrial, Scientific, and Medical (ISM) frequency bands from 2.412 to 2.4835 GHz, depending on your country's regulation. The dongle implements the following methods of transmission:

- ▶ IEEE 802.11b with Direct Sequence Spread Spectrum (DSSS) at 19 dBm: Peak rate 54 Mbps, Peak throughput: 27 Mbps
- ▶ IEEE 802.11g with Orthogonal Frequency Division Multiplexing (OFDM) at 16 dBm: Peak rate 54 Mbps, Peak throughput: 27 Mbps
- ▶ IEEE 802.11n with Orthogonal Frequency Division Multiplexing (OFDM) at 15 dBm:
 - ▶ 1T1R. Peak rate: 150 Mbps, Peak throughput: 90 Mbps
 - ▶ 1T2R. Peak rate: 300 Mbps, Peak throughput: Rx 160 Mbps
 - ▶ 2T2R. Peak rate: 300 Mbps, Peak throughput: Rx 260 Mbps

Wireless and Security Module (Laird)

The Wireless and Security Module uses the Industrial, Scientific, and Medical (ISM) frequency bands from 1.400 to 2.4835 GHz, and from 5.100 to 5.800 GHz. The module implements four different methods of transmission:

- ▶ IEEE 802.11a with Orthogonal Frequency Division Multiplexing (OFDM) at 11 dBm \pm 2 dBm @ 54 Mbps
- ▶ IEEE 802.11b with Direct Sequence Spread Spectrum (DSSS) at 16 dBm \pm 2.0 dBm @ 11 Mbps
- ▶ IEEE 802.11g with Orthogonal Frequency Division Multiplexing (OFDM) at 13 dBm \pm 2.0 dBm @ 54 Mbps
- ▶ IEEE 802.11n with Orthogonal Frequency Division Multiplexing (OFDM) at 12 dBm \pm 2.0 dBm (802.11gn) @ MCS7

Compatible accessories and peripherals (SonoSite Edge II)

FUJIFILM SonoSite has tested the SonoSite Edge II ultrasound system with the following accessories and peripherals and has demonstrated compliance to the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014.

You may use these FUJIFILM SonoSite accessories and third-party peripherals with the SonoSite Edge II ultrasound system.

WARNINGS

- ▶ Use of the accessories with medical systems other than the Edge ultrasound system may result in increased emissions or decreased immunity of the medical system.
- ▶ Use of accessories other than those specified may result in increased emissions or decreased immunity of the ultrasound system.
- ▶ The ultrasound system should not be used in a domestic establishment or connected to the public mains network.

Table 2: Accessories and peripherals compatible with Edge II ultrasound system

Description	Maximum Cable Length
C8x transducer	6.0 ft/1.8 m
C11x transducer	6.0 ft/1.8 m
C35x transducer	5.5 ft/1.7 m
rC60xi transducer standard/armored	5.5 ft/1.7 m
HFL38xi transducer standard/armored	5.5 ft/1.7 m
HFL50x transducer	5.7 ft/1.7 m
HSL25x transducer	7.5 ft/2.3 m
ICTx transducer	5.5 ft/1.7 m
L25x transducer standard/armored	7.5 ft/2.3 m
L38xi transducer standard/armored	5.5 ft/1.7 m

Table 2: Accessories and peripherals compatible with Edge II ultrasound system (continued)

Description	Maximum Cable Length
L52x transducer ^a	7.9 ft/2.4 m
P10x transducer	6.0 ft/1.8 m
rP19x transducer standard/armored	6.0 ft/1.8 m
TEExi transducer	7.2 ft/2.2 m
Bar code scanner	4.8 ft/1.5 m
Battery for PowerPack	—
Battery Pack	—
Battery PowerPack	—
Black & white printer	—
Black & white printer power cable	3.3 ft/1 m
Color printer	—
Color printer power cable	3.3 ft/1 m
Color printer video cable	6 ft/1.8 m
ECG lead wires	24 in/0.6 m
ECG module	5.8 ft/1.8 m
ECG slave cable	8 ft/2.4 m
SonoSite Edge II Dock	—
SonoSite Edge II Stand	—
Footswitch	9.8 ft/3 m
Petite mouse	6 ft /1.8 m
Power cord (system)	10 ft/3 m
Power supply with DC cable	6.8 ft/2 m
Power supply AC cable	39 in/1 m
PowerPark	—
Triple Transducer Connect	—

Table 2: Accessories and peripherals compatible with Edge II ultrasound system (continued)

Description	Maximum Cable Length
USB wireless adapter	—
<p>For transducers, the maximum cable length is measured between the strain reliefs. The stated length do not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, and inside the transducer connector.</p> <p>^aThe L52x transducer is for veterinary use only.</p>	

Compatible accessories and peripherals (SonoSite SII)

FUJIFILM SonoSite has tested the SonoSite SII Ultrasound System with the following accessories and peripherals and has demonstrated compliance to the requirements of IEC 60601-1-2:2007 and IEC 60601-1-2:2014.

You may use these FUJIFILM SonoSite accessories and third-party peripherals with the SonoSite SII ultrasound system.

WARNINGS

- ▶ Use of the accessories with medical systems other than the SonoSite SII ultrasound system may result in increased emissions or decreased immunity of the medical system.
- ▶ Use of accessories other than those specified may result in increased emissions or decreased immunity of the ultrasound system.

Table 3: Accessories and peripherals compatible with SonoSite SII ultrasound system

Description	Maximum Cable Length
C8x transducer	6.0 ft/1.8 m
C11x transducer	6.0 ft/1.8 m
C35x transducer	5.5 ft/1.7 m
rC60xi transducer standard/armored	5.5 ft/1.7 m
HFL38xi transducer standard/armored	5.5 ft/1.7 m
HFL50x transducer	5.7 ft/1.7 m
HSL25x transducer	7.5 ft/2.3 m
ICTx transducer	5.5 ft/1.7 m

Table 3: Accessories and peripherals compatible with SonoSite SII ultrasound system

Description	Maximum Cable Length
L25x transducer standard/armored	7.5 ft/2.3 m
L38xi transducer standard/armored	5.5 ft/1.7 m
L52x transducer	7.9 ft/2.4 m
P10x transducer	6.0 ft/1.8 m
rP19x transducer standard/armored	6.0 ft/1.8 m
Bar code scanner	4.8 ft/1.5 m
Battery for PowerPack	—
Battery Pack	—
Battery PowerPack	—
Black & white printer	—
Black & white printer power cable	3.3 ft/1 m
Black & white printer control cable	6.0 ft/1.8 m
Black & white printer video cable	6.2 ft/1.9 m
Footswitch	9.8 ft/3 m
Footswitch USB extension cable	6.5 ft/2 m
SonoSite SII Stand	—
Power cord (system)	10 ft/3 m
Power supply with DC cable	6.8 ft/2 m
Power supply AC cable	39 in/1 m
PowerPark	—
USB wireless adapter	—
USB memory drive	—

For transducers, the maximum cable length is measured between the strain reliefs. The stated length do not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, and inside the transducer connector.

Manufacturer's declaration

The tables in this section document the intended use environment and EMC compliance levels of the system. For maximum performance, ensure that the system is used in the environments described in this table.

The system is intended for use in the electromagnetic environment specified below.

Table 4: Manufacturer's Declaration - Electromagnetic Emissions per IEC 60601-1-2:2007 and IEC 60601-1-2:2014

Emissions Test	Compliance	Electromagnetic Environment
RF emissions CISPR 11	Group 1	The Edge II and SII ultrasound systems use RF energy only for their internal functions. Therefore, their RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class A	The Edge II and SII ultrasound systems are suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network which supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

The system is intended for use in the electromagnetic environment specified below.

Table 5: Manufacturer's Declaration - Electromagnetic Immunity per IEC 60601-1-2:2007

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Electrostatic Discharge (ESD) IEC 61000-4-2	± 2.0KV, ± 4.0KV, ± 6.0KV contact ± 2.0KV, ± 4.0KV, ± 8.0KV air	± 2.0KV, ± 4.0KV, ± 6.0KV contact ± 2.0KV, ± 4.0KV, ± 8.0KV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient burst IEC 61000-4-4	± 2KV on the mains ± 1KV on signal lines	± 2KV on the mains ± 1KV on signal lines	Mains power quality should be that of a typical commercial or hospital environment.

Table 5: Manufacturer's Declaration - Electromagnetic Immunity per IEC 60601-1-2:2007

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Surge IEC 61000-4-5	± 1KV line(s) to line(s) ± 2KV line(s) to earth	± 1KV line(s) to line(s) ± 2KV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5s	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FUJIFILM SonoSite ultrasound system requires continued operation during power mains interruptions, it is recommended that the FUJIFILM SonoSite ultrasound system be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Field IEC 61000-4-8	3 A/m	3 A/m	If image distortion occurs, it may be necessary to position the FUJIFILM SonoSite ultrasound system further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the Intended installation location to assure that it is sufficiently low.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the FUJIFILM SonoSite ultrasound system including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation Distance $d = 1.2 \sqrt{P}$

Table 5: Manufacturer’s Declaration – Electromagnetic Immunity per IEC 60601-1-2:2007

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	<p>$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (10) according to the transmitter manufacturer and it is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b. Interference may occur in the vicinity of equipment marked with the following symbol:</p>  <p>(IEC 60417 No. 417-IEC-5140: “Source of non-ionizing radiation”)</p>
Note	<p>U_T is the AC mains voltage prior to application of the test level. At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>		
<p>a. Field strengths from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FUJIFILM SonoSite ultrasound system is used exceeds the applicable RF compliance level above, the FUJIFILM SonoSite ultrasound system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the FUJIFILM SonoSite ultrasound system.</p> <p>b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 6: Manufacturer's Declaration - Electromagnetic Immunity per IEC 60601-1-2:2014

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Electrostatic Discharge (ESD) IEC 61000-4-2	$\pm 8.0\text{KV}$, contact $\pm 2.0\text{KV}$, $\pm 4.0\text{KV}$, $\pm 8.0\text{KV}$ air, $\pm 15\text{KV}$	$\pm 8.0\text{KV}$, $\pm 4.0\text{KV}$, $\pm 8.0\text{KV}$ air, $\pm 15\text{KV}$	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient burst IEC 61000-4-4	$\pm 2\text{KV}$ on the mains $\pm 1\text{KV}$ on signal lines	$\pm 2\text{KV}$ on the mains $\pm 1\text{KV}$ on signal lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	$\pm 1\text{KV}$ line(s) to line(s) $\pm 2\text{KV}$ line(s) to earth	$\pm 1\text{KV}$ line(s) to line(s) $\pm 2\text{KV}$ line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U_T for 0.5 cycle 0% U_T for 5 cycles 70% U_T (30% dip in U_T) for 500 msec <5% U_T (>95% dip in U_T) for 5s	0% U_T for 0.5 cycle 0% U_T for 5 cycles 70% U_T (30% dip in U_T) for 500 msec <5% U_T (>95% dip in U_T) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the FUJIFILM SonoSite ultrasound system requires continued operation during power mains interruptions, it is recommended that the FUJIFILM SonoSite ultrasound system be powered from an uninterruptible power supply or a battery.
Power Frequency Magnetic Field IEC 61000-4-8	30 A/m	30 A/m	If image distortion occurs, it may be necessary to position the FUJIFILM SonoSite ultrasound system further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the Intended installation location to assure that it is sufficiently low.

Table 6: Manufacturer’s Declaration – Electromagnetic Immunity per IEC 60601-1-2:2014

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands	3 Vrms 6 Vrms in ISM bands	Portable and mobile RF communications equipment should be used no closer to any part of the FUJIFILM SonoSite ultrasound system including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended Separation Distance $d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	3 V/m 80 MHz to 2.7 GHz	$d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (10) according to the transmitter manufacturer and it is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol:  (IEC 60417 No. 417-IEC-5140: “Source of non-ionizing radiation”)
Proximity fields from wireless communications equipment	Per 60601-1-2:2014 Table 9	Per 60601-1-2:2014 Table 9	

Table 6: Manufacturer's Declaration - Electromagnetic Immunity per IEC 60601-1-2:2014

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment
Note	<p>U_T is the AC mains voltage prior to application of the test level. At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>		
	<p>a. Field strengths from fixed transmitters such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FUJIFILM SonoSite ultrasound system is used exceeds the applicable RF compliance level above, the FUJIFILM SonoSite ultrasound system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the FUJIFILM SonoSite ultrasound system.</p> <p>b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>		

FCC Caution: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Compatible accessories and peripherals

The following warning was added to the SonoSite Edge II and SonoSite SII user guides; the revision will be made in the next update.

WARNING

If peripherals are connected to the system, ensure that the system and peripherals are connected to the same AC Mains branch circuit.

Labeling symbols

Table 7: Standards labeling symbols



Symbol	Title	Standards development organization	Reference number	Description
	Manufacturer	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.1.1	Indicates the medical device manufacturer, as EU directives 90/385/EEC, 93/42/EEC and 98/79/EC
	Non-ionizing electromagnetic radiation	IEC 60601-1-2:2007 Medical Electrical Equipment Part 1-2: General requirements for basic safety and essential performance – collateral standard: electromagnetic compatibility	5.1.1	Indicates generally elevated, potentially hazardous levels of non-ionizing radiation, or to indicate equipment or systems e.g., in the medical electrical area that include RF transmitters or that intentionally apply RF electromagnetic energy for diagnosis or treatment
	European community authorized representative	ISO 15223-1 Medical devices - symbols to be used with medical device labels, labelling and information to be supplied.	5.1.2	Indicates the Authorized representative in the European Community
	Serial number	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.1.7	Indicates the manufacturer's serial number so that a specific medical device can be identified
	Catalog number	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.1.6	Indicates the manufacturer's catalog number so that the medical device can be identified

Table 7: Standards labeling symbols (continued)


Symbol	Title	Standards development organization	Reference number	Description
	Caution	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements	5.4.4	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions, that cannot for a variety of reasons, be presented on the medical device itself
	Fragile handle with care	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.3.1	Indicates a medical device that can be broken or damaged if not handled carefully
	Keep dry	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.3.4	.Indicates a medical device that needs to be protected from moisture
	Temperature limit	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.3.7	Indicates the temperature limits to which the medical device can be safely exposed
	Atmospheric pressure limitations	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.3.9	Indicates the range of atmospheric pressure to which the medical device can be safely exposed
	Humidity limitation	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.3.8	Indicates the range of humidity to which the medical device can be safely exposed

Table 7: Standards labeling symbols (continued)

Symbol	Title	Standards development organization	Reference number	Description
IPX7	Degree of ingress protection provided by enclosure	IEC 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance	D.3	Protected against the effects of temporary immersion
	Refer to instruction manual/booklet	IEC 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance	D.2-10	Follow instructions for use (used in accordance with IEC 60601-1)
	Consult instructions for use	ISO 15223-1:2016 Medical devices – symbols to be used with medical device labels, labelling, and information to be supplied – Part 1: General requirements	5.4.3	Indicates the need for the user to consult the instructions for use
	Alternating current	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5032	Indicates on the rating plate, that the equipment is suitable for alternating current only, in order to identify appropriate terminals
	CE marking	Council Directive 93/42/EEC	Article 17 Annex XII	Signifies European Technical Conformity
	Conformité Européenne Notified Body Reference No.: 2797	Council Directive 93/42/EEC	Article 17 Annex XII	Indicates European technical conformity and identification of notified body responsible for implementation of the procedures set out in Annexes II, IV, V, and VI.
	Dangerous voltage	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5036	Indicates hazards arising from dangerous voltage

Table 7: Standards labeling symbols (continued)

Symbol	Title	Standards development organization	Reference number	Description
	Stacking limit by number	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	2403	Indicates that the items are not to be vertically stacked higher than the specified number of items
	Caution hot	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5041	Indicates that the marked item can be hot and should not be touched without taking care
	Caution, static magnetic field hazard	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	6204	Identifies areas with potentially hazardous static magnetic fields and forces in an installation
	Type BF applied parts	IEC 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance	D.2-10	Identifies type BF applied part complying with IEC 60601-1
	Defibrillation-proof type CF applied part	IEC 60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance	D.1-27	Identifies a defibrillation-proof type CF applied part complying with IEC 60601-1
	Electrostatic sensitive device	IEC 60417:2002 Graphical Symbols For Use On Equipment	5134	Indicates packages containing electrostatic sensitive devices, or identifies a device or a connector that has not been tested for immunity to electrostatic discharge
	Regulatory Compliance Mark (RCM)	AS/NZS3820	—	Indicates C-Tick-Regulatory Compliance Mark for Australia and New Zealand Device complies with relevant Australian and New Zealand regulations for electronic devices.

Table 7: Standards labeling symbols (continued)

Symbol	Title	Standards development organization	Reference number	Description
	Batch code, date code, or lot code type of control number	ISO 15223-1 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General Requirements	5.1.5	Indicates manufacturer's batch code so that the batch or lot can be identified
	Biological risk	ISO 7010 - Graphical symbols -- Safety colors and safety signs	W009	To warn of biological Hazard
	INMETRO Safety Symbols	—	—	Indicates Brazil - Accredited certification body by the National Institute of Metrology Standardization and Industrial Quality (INMETRO)
	Canadian Standard Association Certification Mark	—	—	CSA certification mark signifying that the product complies with the applicable CSA and ANSI/UL requirements and is authorized for use in Canada and the US.
	Recycle: Electronic Equipment	BS EN 50419:2016 Marking of Electrical and Electronic Equipment in accordance with Directive 2012/19/EU (WEEE).	Annex IX	Do Not Throw in Trash
	Corrugated recycle	—	—	Shipping box is made of corrugated cardboard and should be recycled accordingly
	Date of manufacture	ISO 7000- Graphical symbols for Use on Equipment	5.1.3	To indicate the date on which a product was manufactured
	Direct current (DC)	—	—	—

Table 7: Standards labeling symbols (continued)



Symbol	Title	Standards development organization	Reference number	Description
	21 Part 15	Federal Communications Commission (FCC) Declaration of conformity	Federal Communications Commission	FCC-Tested to Federal Communications Commission requirements Device complies with relevant FCC regulations for electronic devices
GEL	Gel	—	—	—
	Resy- Recycling Symbol	—	—	Paper recycle
	Degree of Ingress Protection Provided by Enclosure	IEC 60601-1 Medical Electrical Equipment Part 1: General requirements for basic safety and essential performance	D.3	Protected against the effects of temporary immersion in water. Submersible Protected against the effects of temporary immersion.
	Degree of Ingress Protection Provided by Enclosure	IEC 60601-1 Medical Electrical Equipment Part 1: General requirements for basic safety and essential performance	D.3	Protected against the effects of temporary immersion in water. Water-Tight Equipment Protected against the effects of extended immersion
	—	—	—	Indicates handle with care
	—	—	—	Indicates follow manufacturer's instructions for disinfecting time
	—	—	—	Indicates disinfect transducer
	Maximum weight load	IEC 60601-1 Medical Electrical Equipment Part 1: General requirements for basic safety and essential performance	7.2.21	Indicates total weight of the equipment, including the safe working load

Table 7: Standards labeling symbols (continued)

Symbol	Title	Standards development organization	Reference number	Description
	Underwriters Laboratories Certification Mark	—	—	Certification mark for electrical shock, fire, and mechanical hazards only
	UL Product Certification.	—	—	The product or company has successfully met stringent standards for product safety.
	China Pollution Control (10)	ISO 7000:2014 Graphical symbols for use on equipment	1135	Pollution Control Logo. (Applies to all parts/products listed in the China RoHS disclosure table. May not appear on the exterior of some parts/products because of space limitations.)
	China Compulsory Certificate mark ("CCC Mark"). Compulsory safety mark for compliance to Chinese national standards for many products sold in the People's Republic of China.	—	—	—
	Sterilized using ethylene oxide	ISO 15223-1 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General Requirements	5.2.3	Indicates a medical device that has been sterilized using ethylene oxide

Table 7: Standards labeling symbols (continued)

Symbol	Title	Standards development organization	Reference number	Description
	Sterilized using irradiation	ISO 15223-1 Medical devices - Symbols to be used with medical device labels, labeling and information to be supplied - Part 1: General Requirements	5.2.4	Indicates a medical device that has been sterilized using irradiation
	(SII only) Catalog number	ISO 15223-1:2016 Medical Devices - Symbols To Be Used With Medical Device Labels, Labeling, And Information To Be Supplied - Part 1: General Requirements	5.1.6	Indicates the manufacturer's Catalogue number so that the medical device can be identified.
	(SII only) Dangerous voltage	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5036	Indicates hazards arising from dangerous voltage
	(SII only) For indoor use only	ISO 7000 / IEC 60417 Graphical symbols for use on equipment	5957	Identifies electrical equipment designed primarily for indoor use

Specifications

Supported transducers (SonoSite SII)

The following redundant section was removed from the SonoSite SII user guide. The same information is present in Table 9-2 of the user guide; the revision will be made in the next update.

Table 8: Supported transducers

Description	Maximum Cable Length
C11x transducer	6.0 ft/1.8 m
rC60xi transducer standard/armored	5.5 ft/1.7 m
HFL38xi transducer standard/armored	5.5 ft/1.7 m

For transducers, the maximum cable length is measured between the strain reliefs. The stated length does not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, and inside the transducer connector.

Description	Maximum Cable Length
HFL50x transducer	5.5 ft/1.7 m
HSL25x transducer	8.0 ft/2.4 m
ICTx transducer	5.5 ft/1.7 m
L25x transducer	7.5 ft/2.3 m
L38xi transducer standard/armored	5.5 ft/1.7 m
L52x transducer	7.9 ft/2.4 m
P10x transducer	6.0 ft/1.8 m
rP19x transducer standard/armored	6.0 ft/1.8 m

For transducers, the maximum cable length is measured between the strain reliefs. The stated length does not include the lengths of cable in the following locations: underneath the strain reliefs, inside the transducer enclosure, and inside the transducer connector.

Standards

Electromechanical safety standards

Table 9: Electromechanical safety standards

Standard	Description
ANSI/AAMI ES60601-1:2005/(R) 2012, and A1:2012	Medical electrical equipment, Part 1: General requirements for basic safety and essential performance (Consolidated Edition 3.1)
CAN/CSA C22.2 No. 60601-1:2014 (Edition 3.1)	Medical electrical equipment - Part 1: General Requirements for Basic Safety and Essential Performance
IEC 60601-1:2012 (Edition 3.1)	Medical electrical equipment - Part 1: General Requirements for Basic Safety and Essential Performance
IEC 60601-2-37:2015	Medical Electrical Equipment - Part 2-37: Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment
IEC 60601-1-6:2013	Medical Electrical Equipment, Part 1-6: General requirements for basic safety and essential performance - Collateral Standard: Usability
JIS T0601-1:2012 (3rd Edition)	Japanese Industrial Standard, Medical electrical equipment - Part 1: General requirements for basic safety and essential performance

Acoustic output

ALARA principle

Applying the ALARA principle

The system imaging mode selected by the qualified ultrasound user is determined by the diagnostic information required. 2D imaging provides anatomical information; CPD imaging provides information about the energy or amplitude strength of the Doppler signal over time at a given anatomical location and is used for detecting the presence of blood flow; Color imaging provides information about the energy or amplitude strength of the Doppler signal over time at a given anatomical location and is used for detecting the presence, velocity, and direction of blood flow; Tissue Harmonic Imaging (THI) uses higher received frequencies to reduce clutter, artifact, and improve resolution on the 2D image. Understanding the nature of the imaging mode used allows the qualified ultrasound user to apply the ALARA principle.

Prudent use of ultrasound requires that patient exposure to ultrasound be limited to the lowest ultrasound output for the shortest time necessary to achieve acceptable diagnostic results. Decisions that support prudent use are based on the type of patient, exam type, patient history, ease or difficulty of obtaining diagnostically useful information, and potential localized heating of the patient due to transducer surface temperature.

The system has been designed to ensure that temperature at the face of the transducer will not exceed the limits established in IEC 60601-2-37: Particular requirement for the safety of ultrasound medical diagnostic and monitoring equipment. See **“Transducer surface temperature rise”** on page 10-9. In the event of a device malfunction, there are redundant controls that limit transducer power. This is accomplished by an electrical design that limits both power supply current and voltage to the transducer.

The sonographer uses the system controls to adjust image quality and limit ultrasound output. The system controls are divided into three categories relative to output: controls that directly affect output, controls that indirectly affect output, and receiver controls.

Direct controls

The system does not exceed a spatial peak temporal average intensity (ISPTA) of 720 mW/cm^2 for all imaging modes. (For either the Ophthalmic or Orbital exam, the acoustic output is limited to the following values: ISPTA does not exceed 50 mW/cm^2 ; TI does not exceed 1.0, and MI does not exceed 0.23.) The mechanical index (MI) and thermal index (TI) may exceed values greater than 1.0 on some transducers in some imaging modes. One may monitor the MI and TI values and adjust the controls to reduce these values. See **“Guidelines for reducing MI and TI”** on page 10-3. Additionally, one means for meeting the ALARA principle is to set the MI or TI values to a low index value and then modifying this level until a satisfactory image or Doppler mode is obtained. For more information on MI and TI, refer to Medical Ultrasound Safety, AIUM (a copy is included with each system) and IEC 60601-2-37 Annex “Guidance on the interpretation of TI and MI to be used to inform the operator.”

Output display

Related guidance documents

Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers, FDA, 2008.

Medical Ultrasound Safety, American Institute of Ultrasound in Medicine (AIUM), 2014. (A copy is included with each system.)

Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment, NEMA UD2-2004.

IEC 60601-2-37: 2015, Particular requirements for the basic safety and essential performance of ultrasonic diagnostic and monitoring equipment.

Transducer surface temperature rise

Table 10-4 and Table 10-5 list the measured surface temperature rise from ambient ($23^{\circ}\text{C} \pm 3^{\circ}\text{C}$) of transducers used on the ultrasound system. The temperatures were measured in accordance with IEC 60601-2-37 with controls and settings positioned to give maximum temperatures.

Acoustic output measurement

Since the initial use of diagnostic ultrasound, the possible human biological effects (bioeffects) from ultrasound exposure have been studied by various scientific and medical institutions. In October 1987, the American Institute of Ultrasound in Medicine (AIUM) ratified a report from its Bioeffects Committee (Bioeffects Considerations for the Safety of Diagnostic Ultrasound, *J Ultrasound Med.*, Sept. 1988: Vol. 7, No. 9 Supplement). The report, sometimes referred to as *the Stowe Report*, reviewed available data on possible effects of ultrasound exposure. Another report, *Bioeffects and Safety of Diagnostic Ultrasound*, dated January 28, 1993, provides more current information.

The acoustic output for this ultrasound system has been measured and calculated in accordance with “Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment” (NEMA UD2-2004) and IEC 60601-2-37: 2015, Particular requirements for the basic safety and essential performance of ultrasonic diagnostic and monitoring equipment.

Acoustic output tables

The format of the acoustic output tables has been updated.

Transducer model: C8x Operating mode: 2D	33
Transducer model: C8x Operating mode: M mode	34
Transducer model: C8x Operating mode: Color/CPD	35
Transducer model: C8x Operating mode: PW Doppler	36

Transducer model: C11x Operating mode: PW Doppler	37
Transducer model: C35x Operating mode: 2D	38
Transducer model: C35x Operating mode: PW Doppler	39
Transducer model: HFL38xi (Ophthalmic use) Operating mode: 2D	40
Transducer model: HFL38xi (Ophthalmic use) Operating mode: M mode	41
Transducer model: HFL38xi (Ophthalmic use) Operating mode: Color/CPD	42
Transducer model: HFL38xi (Ophthalmic use) Operating mode: PW Doppler	43
Transducer model: HFL38xi Operating mode: 2D	44
Transducer model: HFL38xi Operating mode: M mode	45
Transducer model: HFL38xi Operating mode: Color/CPD	46
Transducer model: HFL38xi Operating mode: PW Doppler	47
Transducer model: HFL50x Operating mode: 2D	48
Transducer model: HFL50x Operating mode: M mode	49
Transducer model: HFL50x Operating mode: Color	50
Transducer model: HFL50x Operating mode: PW Doppler	51
Transducer model: HSL25x (Ophthalmic use) Operating mode: 2D	52
Transducer model: HSL25x (Ophthalmic use) Operating mode: M mode	53
Transducer model: HSL25x (Ophthalmic use) Operating mode: Color/CPD	54
Transducer model: HSL25x (Ophthalmic use) Operating mode: PW Doppler	55
Transducer model: HSL25x Operating mode: 2D	56
Transducer model: HSL25x Operating mode: Color/CPD	57
Transducer model: HSL25x Operating mode: PW Doppler	58
Transducer model: ICTx Operating mode: PW Doppler	59
Transducer model: L25x (Ophthalmic use) Operating mode: 2D	60
Transducer model: L25x (Ophthalmic use) Operating mode: M mode	61
Transducer model: L25x (Ophthalmic use) Operating mode: Color/CPD	62
Transducer model: L25x (Ophthalmic use) Operating mode: PW Doppler	63
Transducer model: L25x Operating mode: 2D	64
Transducer model: L25x Operating mode: Color/CPD	65
Transducer model: L25x Operating mode: PW Doppler	66
Transducer model: L38xi Operating mode: 2D	67
Transducer model: L38xi Operating mode: M mode	68
Transducer model: L38xi Operating mode: Color/CPD	69
Transducer model: L38xi Operating mode: PW Doppler	70
Transducer model: P10x Operating mode: Color	71
Transducer model: P10x Operating mode: CW Doppler	72
Transducer model: P10x Operating mode: PW Doppler	73
Transducer model: rC60xi Operating mode: 2D	74
Transducer model: rC60xi Operating mode: M mode	75
Transducer model: rC60xi Operating mode: Color/CPD	76
Transducer model: rC60xi Operating mode: PW Doppler	77
Transducer model: rP19x (Orbital use) Operating mode: 2D	78
Transducer model: rP19x (Orbital use) Operating mode: M mode	79
Transducer model: rP19x (Orbital use) Operating mode: Color/CPD	80
Transducer model: rP19x (Orbital use) Operating mode: PW Doppler	81
Transducer model: rP19x Operating mode: 2D	82
Transducer model: rP19x Operating mode: M mode	83

Transducer model: rP19x Operating mode: Color/CPD	84
Transducer model: rP19x Operating mode: CW Doppler	85
Transducer model: rP19x Operating mode: PW Doppler	86
Transducer model: TEExi Operating mode: CW Doppler	87
Transducer model: TEExi Operating mode: PW Doppler	88

Table 10: Transducer model: C8x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.1	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.48					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
	f_{awf} (MHz)	5.53	#		#		#
Other information	pr (Hz)	9524					
	srr (Hz)	18.6					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	264					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	18.3					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	25.6					
	p_r at z_{pii} (MPa)	3.11					
Operating controls	Exam type	Pro					
	Optimization	Pen					
	Depth (cm)	2.5–3.2					
	MB	Off					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 11: Transducer model: C8x Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.91					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			#			
	z_b (cm)					#	
	z_{MI} (cm)	1.1					
	$z_{pii,\alpha}$ (cm)	1.1					
	f_{awf} (MHz)	5.07	#		#		#
Other information	pr (Hz)	800					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	433					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	149					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	226					
Operating controls	p_r at z_{pii} (MPa)	3.57					
	Exam type	Pro					
	Optimization	Pen					
	Depth (cm)	4.2					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 12: Transducer model: C8x Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.68					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
	f_{awf} (MHz)	4.82	#		#		#
Other information	p_{rr} (Hz)	2548					
	s_{rr} (Hz)	26					
	n_{pps}	12					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	381					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sij,\alpha}$ (mW/cm^2)	132					
	I_{spta} at z_{pii} or z_{sij} (mW/cm^2)	176					
	p_r at z_{pii} (MPa)	3.1					
Operating controls	Exam type	Pro					
	Mode	CVD					
	2D optimization/depth (cm)	Pen/1.5–1.9					
	Color optimization/PRF (Hz)	High/any					
	Color box position/size	Narrow/any					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 13: Transducer model: C8x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.0	(a)		1.4		(b)
Index component value			#	#	0.5	1.4	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.28					
	P (mW)		#		23.1		#
	P_{1x1} (mW)		#		23.1		
	z_s (cm)			#			
	z_b (cm)					1.0	
	z_{MI} (cm)	1.8					
	$z_{pii,\alpha}$ (cm)	1.8					
	f_{awf} (MHz)	4.80		#	4.80		#
Other information	pr (Hz)	1008					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	263					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	334					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	616					
p_r at z_{pii} (MPa)	3.1						
Operating controls	Exam type	Pro			Pro		
	Sample volume size (mm)	1			1		
	Sample volume position	Zone 4			Zone 4		
	PRF (Hz)	1008			1008		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 14: Transducer model: C11x Operating mode: PW Doppler

Index label	MI	TIS		TIB		TIC
		At surface	Below surface	At surface	Below surface	At surface
Maximum index value	(a)	(a)		1.5		1.1
Index component value		#	#	0.5	1.5	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#				
	P (mW)		#	24.6		21.7
	P_{1x1} (mW)		#	24.6		
	z_s (cm)		#			
	z_b (cm)				1.7	
	z_{MI} (cm)	#				
	$z_{pii,\alpha}$ (cm)	#				
Other information	f_{awf} (MHz)	#	#	4.37		4.36
	pr_r (Hz)	#				
	srr (Hz)	#				
	η_{pps}	#				
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#				
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#				
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#				
Operating controls	p_r at z_{pii} (MPa)	#				
	Exam type			Nrv		Nrv
	Sample volume size (mm)			1		7
	Sample volume position			Zone 1		Zone 0
	PRF (Hz)			10417		6250

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 15: Transducer model: C35x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.0	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.8					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	3.3					
	$z_{pii,\alpha}$ (cm)	3.3					
	f_{awf} (MHz)	3.45	#		#		#
Other information	pr (Hz)	1021					
	srr (Hz)	7.98					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	250					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	8.6					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	16.5					
Operating controls	p_r at z_{pii} (MPa)	2.61					
	Exam type	Msk					
	Optimization	Res					
	Depth (cm)	8.3					
	MB	N/A					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 16: Transducer model: C35x Operating mode: PW Doppler

Index label	MI	TIS		TIB		TIC
		At surface	Below surface	At surface	Below surface	At surface
Maximum index value	(a)	1.5		2.6		(b)
Index component value		1.5	1.0	1.0	2.6	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#				
	P (mW)		72.8	47.1		#
	P_{1x1} (mW)		71.1	47.1		
	z_s (cm)			1.4		
	z_b (cm)				0.50	
	z_{MI} (cm)	#				
	$z_{pii,\alpha}$ (cm)	#				
Other information	f_{awf} (MHz)	#	4.35	4.37		#
	pr_r (Hz)	#				
	srr (Hz)	#				
	η_{pps}	#				
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#				
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#				
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#				
Operating controls	p_r at z_{pii} (MPa)	#				
	Exam type		Spine		Spine	
	Sample volume size (mm)		2		1	
	Sample volume position		Zone 5		Zone 0	
			6250		15625	

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 17: Transducer model: HFL38xi (Ophthalmic use) Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.007		0.007		(b)
Index component value			0.007	0.007	0.007	0.007	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.43					
	P (mW)		0.77		0.77		#
	$P_{1 \times 1}$ (mW)		0.21		0.21		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	2.1					
	$z_{pii,\alpha}$ (cm)	2.1					
	f_{awf} (MHz)	6.59	6.75		6.75		#
Other information	pr (Hz)	11339					
	srr (Hz)	19.7					
	n_{pps}	3					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	11.4					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	0.8					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	1.3					
p_r at z_{pii} (MPa)	0.7						
Operating controls	Exam type	Oph	Oph		Oph		
	Optimization	Pen	Res		Res		
	Depth (cm)	4.9	4.9		4.9		
	MB	On	On		On		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 18: Transducer model: HFL38xi (Ophthalmic use) Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.003		0.004		(b)
Index component value			0.003	0.002	0.002	0.004	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.44					
	P (mW)		0.087		0.064		#
	P_{1x1} (mW)		0.087		0.064		
	z_s (cm)			0.9			
	z_b (cm)					1.10	
	z_{MI} (cm)	1.0					
	$z_{pii,\alpha}$ (cm)	1.0					
Other information	f_{awf} (MHz)	6.58	6.86		6.78		#
	pr (Hz)	800					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	10.3					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	1.0					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	1.7					
Operating controls	p_r at z_{pii} (MPa)	0.55					
	Exam type	Oph	Oph		Oph		
	Optimization	Pen	Res		Res		
	Depth (cm)	1.5	6.0		4.0		

(a) This index is not required for this operating mode; value is <1 .

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 19: Transducer model: HFL38xi (Ophthalmic use) Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.02		0.02		(b)
Index component value			0.02	0.02	0.02	0.02	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.39					
	P (mW)		1.11		1.11		#
	P_{1x1} (mW)		0.75		0.75		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.9					
	$z_{pii,\alpha}$ (cm)	0.9					
	f_{awf} (MHz)	5.34	5.37		5.37		#
Other information	prr (Hz)	4537					
	srr (Hz)	13.5					
	n_{pps}	13					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	5.5					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	1.3					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	2.1					
p_r at z_{pii} (MPa)	0.46						
Operating controls	Exam type	Oph	Oph		Oph		
	Mode	CVD	CVD		CVD		
	2D optimization/depth (cm)	Pen/1.5	Pen/4.9		Pen/4.9		
	Color optimization/PRF (Hz)	High/7813	High/6944		High/6944		
	Color box position/size	Bottom/small	Def/narrow		Def/narrow		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 20: Transducer model: HFL38xi (Ophthalmic use) Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.18	0.09		0.17		(b)
Index component value			0.09	0.06	0.09	0.17	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.41					
	P (mW)		3.56		3.56		#
	P_{1x1} (mW)		3.56		3.56		
	z_s (cm)			1.1			
	z_b (cm)					1.64	
	z_{MI} (cm)	0.9					
	$z_{pii,\alpha}$ (cm)	0.9					
Other information	f_{awf} (MHz)	5.34	5.33		5.33		#
	pr (Hz)	1302					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	6.6					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	10.9					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	15.0					
p_r at z_{pii} (MPa)	0.48						
Operating controls	Exam type	Oph	Oph		Oph		
	Sample volume size (mm)	1	10		10		
	Sample volume position	Zone 1	Zone 7		Zone 7		
	PRF (Hz)	1302	10417		10417		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 21: Transducer model: HFL38xi Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.05					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
	f_{awf} (MHz)	5.36	#		#		#
Other information	pr (Hz)	2127					
	srr (Hz)	11.1					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	494					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	13.3					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	19.4					
	p_r at z_{pii} (MPa)	3.81					
Operating controls	Exam type	Ven					
	Optimization	Res					
	Depth (cm)	3.3					
	MB	N/A					
	Needle vision	On					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 22: Transducer model: HFL38xi Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.12	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.14					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			#			
	z_b (cm)					#	
	z_{MI} (cm)	1.4					
	$z_{pii,\alpha}$ (cm)	1.4					
Other information	f_{awf} (MHz)	6.75	#		#		#
	pr_r (Hz)	1600					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	388					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	163.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	333.3					
Operating controls	p_r at z_{pii} (MPa)	4.35					
	Exam type	Nrv					
	Optimization	Pen					
	Depth (cm)	4.0					

(a) This index is not required for this operating mode; value is <1 .

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 23: Transducer model: HFL38xi Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.05					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
	f_{awf} (MHz)	5.36	#		#		#
Other information	pr_r (Hz)	2223					
	srr (Hz)	3.3					
	n_{pps}	14					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	494					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	27.4					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	40.1					
	p_r at z_{pii} (MPa)	3.81					
Operating controls	Exam type	SmP					
	Mode	CVD					
	2D optimization/depth (cm)	Res/3.3					
	Color optimization/PRF (Hz)	Low/401					
	Color box position/size	Def/def					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 24: Transducer model: HFL38xi Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	1.1		2.2		(b)
Index component value			1.1	0.8	1.1	2.2	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.69					
	P (mW)		47.7		47.7		#
	P_{1x1} (mW)		47.7		47.7		
	z_s (cm)			1.1			
	z_b (cm)					1.10	
	z_{MI} (cm)	1.0					
	$z_{pii,\alpha}$ (cm)	1.0					
Other information	f_{awf} (MHz)	5.34	4.86		4.86		#
	pr_r (Hz)	1008					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	308					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	102.8					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	210.0					
Operating controls	p_r at z_{pii} (MPa)	3.23					
	Exam type	Nrv	Art		Art		
	Sample volume size (mm)	1	1		1		
	Sample volume position	Zone 3	Zone 7		Zone 7		
	PRF (Hz)	1008	3125		3125		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 25: Transducer model: HFL50x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.051					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
	f_{awf} (MHz)	5.36	#		#		#
Other information	pr (Hz)	2733					
	srr (Hz)	7.2					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	493					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	8.6					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	12.6					
p_r at z_{pii} (MPa)	3.81						
Operating controls	Exam type	Any					
	Optimization	Any					
	Depth (cm)	3.3					
	Mbe	On					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 26: Transducer model: HFL50x Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.14					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			#			
	z_b (cm)					#	
	z_{MI} (cm)	1.4					
	$z_{pii,\alpha}$ (cm)	1.4					
Other information	f_{awf} (MHz)	6.75	#		#		#
	pr_r (Hz)	1600					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	388					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	163.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	333.3					
Operating controls	p_r at z_{pii} (MPa)	4.35					
	Exam type	Any					
	Optimization	Pen					
	Depth (cm)	4					

(a) This index is not required for this operating mode; value is <1 .

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 27: Transducer model: HFL50x Operating mode: Color

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.05					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
	f_{awf} (MHz)	5.36	#		#		#
Other information	pr_r (Hz)	8233					
	srr (Hz)	3.2					
	n_{pps}	14					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	494					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	26.8					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	39.2					
	p_r at z_{pii} (MPa)	3.81					
Operating controls	Exam type	Any					
	Mode	Any					
	Optimization/depth (cm)	Low/3.3					
	PRF (Hz)	Any					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 28: Transducer model: HFL50x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	1.1		1.9		(b)
Index component value			1.1	0.7	1.1	1.9	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.69					
	P (mW)		42.6		42.6		#
	P_{1x1} (mW)		42.6		42.6		
	z_s (cm)			1.1			
	z_b (cm)					1.10	
	z_{MI} (cm)	1.0					
	$z_{pii,\alpha}$ (cm)	1.0					
Other information	f_{awf} (MHz)	5.34	5.34		5.34		#
	pr (Hz)	1008					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	308					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	399.6					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	599.8					
Operating controls	p_r at z_{pii} (MPa)	3.23					
	Exam type	Any	Any		Any		
	Sample volume size (mm)	1	1		1		
	Sample volume position	Zone 3	Zone 7		Zone 7		
	PRF (Hz)	1008	1563–3125		1563–3125		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 29: Transducer model: HSL25x (Ophthalmic use) Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.02		0.02		(b)
Index component value			0.02	0.02	0.02	0.02	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.47					
	P (mW)		1.62		1.62		#
	P_{1x1} (mW)		0.70		0.70		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
	f_{awf} (MHz)	7.65	6.97		6.97		#
Other information	pr (Hz)	12580					
	srr (Hz)	12.3					
	n_{pps}	4					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	13.4					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	0.6					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	1.0					
	p_r at z_{pii} (MPa)	0.58					
Operating controls	Exam type	Oph	Oph		Oph		
	Optimization	Res	Pen		Pen		
	Depth (cm)	1.9	4.3		4.3		
	MB	On	On		On		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 30: Transducer model: HSL25x (Ophthalmic use) Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.01		0.02		(b)
Index component value			0.010	0.009	0.013	0.020	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.47					
	P (mW)		0.45		0.45		#
	P_{1x1} (mW)		0.45		0.45		
	z_s (cm)			0.9			
	z_b (cm)					0.85	
	z_{MI} (cm)	1.0					
	$z_{pii,\alpha}$ (cm)	1.0					
	f_{awf} (MHz)	7.59	6.25		6.25		#
Other information	pr (Hz)	1600					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	14.9					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	2.3					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	4.0					
	p_r at z_{pii} (MPa)	0.61					
Operating controls	Exam type	Oph	Oph		Oph		
	Optimization	Res	Pen		Pen		
	Depth (cm)	1.9	4.3		4.3		

(a) This index is not required for this operating mode; value is <1 .

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 31: Transducer model: HSL25x (Ophthalmic use) Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.06		0.06		(b)
Index component value			0.06	0.06	0.06	0.06	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.42					
	P (mW)		2.9		2.9		#
	$P_{1 \times 1}$ (mW)		1.9		1.9		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.7					
	$z_{pii,\alpha}$ (cm)	0.7					
	f_{awf} (MHz)	6.11	6.10		6.10		#
Other information	p_{rr} (Hz)	3096					
	s_{rr} (Hz)	8.1					
	n_{pps}	14					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	7.5					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	1.1					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	1.6					
	p_r at z_{pii} (MPa)	0.49					
Operating controls	Exam type	Oph	Oph		Oph		
	Mode	CVD	CVD		CVD		
	2D optimization/depth (cm)	Pen/1.9	Pen/5.1		Pen/5.1		
	Color optimization/PRF (Hz)	Low/401	Med/4167		Med/4167		
	Color box position/size	Def/def	Top/short-wide		Top/short-wide		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 32: Transducer model: HSL25x (Ophthalmic use) Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.18	0.12		0.21		(b)
Index component value			0.12	0.08	0.12	0.21	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.44					
	P (mW)		4.0		4.0		#
	P_{1x1} (mW)		4.0		4.0		
	z_s (cm)			0.9			
	z_b (cm)					0.80	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
Other information	f_{awf} (MHz)	6.03	6.03		6.03		#
	pr (Hz)	1953					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	7.4					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	18.4					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	44.9					
p_r at z_{pii} (MPa)	0.56						
Operating controls	Exam type	Oph	Oph		Oph		
	Sample volume size (mm)	1	1		1		
	Sample volume position	Zone 7	Zone 7		Zone 7		
	PRF (Hz)	1953	5208		5208		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 33: Transducer model: HSL25x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.87					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
	f_{awf} (MHz)	6.11	#		#		#
Other information	pr (Hz)	1061					
	srr (Hz)	13.0					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	478					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	12.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	16.4					
	p_r at z_{pii} (MPa)	3.39					
Operating controls	Exam type	Nrv/Msk/ Ven/Art					
	Optimization	Any					
	Depth (cm)	1.9–2.2					
	Mbe	On					

(a) This index is not required for this operating mode; value is < 1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 34: Transducer model: HSL25x Operating mode: Color/CPD

Index label	MI	TIS		TIB		TIC
		At surface	Below surface	At surface	Below surface	At surface
Maximum index value	1.0	(a)		(a)		(b)
Index component value		#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.35				
	P (mW)		#	#		#
	$P_{1 \times 1}$ (mW)		#	#		
	z_s (cm)			—		
	z_b (cm)				—	
	z_{MI} (cm)	0.8				
	$z_{pii,\alpha}$ (cm)	0.8				
	f_{awf} (MHz)	6.11	#		#	
Other information	p_{rr} (Hz)	3079				
	s_{rr} (Hz)	8.0				
	n_{pps}	14				
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	276				
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	47.6				
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	63.9				
	p_r at z_{pii} (MPa)	2.78				
Operating controls	Exam type	Sup				
	Mode	CVD				
	2D optimization/depth (cm)	Pen/3.1				
	Color optimization/PRF (Hz)	Low/401				
	Color box position/size	Def/def				

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 35: Transducer model: HSL25x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		1.5		(b)
Index component value			#	#	0.8	1.5	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		28.1		#
	P_{1x1} (mW)		#		28.1		
	z_s (cm)			#			
	z_b (cm)					0.75	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	#		6.00		#
Other information	pr (Hz)	#					
	srr (Hz)	#					
	n_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#					
p_r at z_{pii} (MPa)	#						
Operating controls	Exam type				Nrv		
	Sample volume size (mm)				8		
	Sample volume position				Zone 7		
	PRF (Hz)				1953		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 36: Transducer model: ICTx Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		1.2		(b)
Index component value			#	#	0.3	1.2	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		16.3		#
	P_{1x1} (mW)		#		16.3		
	z_s (cm)			#			
	z_b (cm)					1.60	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
Other information	f_{awf} (MHz)	#	#		4.36		#
	pr_r (Hz)	#					
	srr (Hz)	#					
	η_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#					
Operating controls	p_r at z_{pii} (MPa)	#					
	Exam type				Any		
	Sample volume size (mm)				3		
	Sample volume position				Zone 1		
	PRF (Hz)				Any		

(a) This index is not required for this operating mode; value is < 1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 37: Transducer model: L25x (Ophthalmic use) Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.02		0.02		(b)
Index component value			0.02	0.02	0.02	0.02	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.47					
	P (mW)		1.62		1.62		#
	P_{1x1} (mW)		0.70		0.70		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
	f_{awf} (MHz)	7.65	6.97		6.97		#
Other information	pr (Hz)	12580					
	srr (Hz)	12.3					
	n_{pps}	4					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	13.4					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	0.6					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	1.0					
p_r at z_{pii} (MPa)	0.58						
Operating controls	Exam type	Oph	Oph		Oph		
	Optimization	Res	Pen		Pen		
	Depth (cm)	1.9	4.3		4.3		
	MB	On	On		On		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 38: Transducer model: L25x (Ophthalmic use) Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.010		0.020		(b)
Index component value			0.010	0.009	0.013	0.020	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.47					
	P (mW)		0.45		0.45		#
	P_{1x1} (mW)		0.45		0.45		
	z_s (cm)			0.9			
	z_b (cm)					0.85	
	z_{MI} (cm)	1.0					
	$z_{pii,\alpha}$ (cm)	1.0					
Other information	f_{awf} (MHz)	7.59	6.25		6.25		#
	pr (Hz)	1600					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	14.9					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	2.3					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	4.0					
Operating controls	p_r at z_{pii} (MPa)	0.61					
	Exam type	Oph	Oph		Oph		
	Optimization	Res	Pen		Pen		
	Depth (cm)	1.9	4.3		4.3		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 39: Transducer model: L25x (Ophthalmic use) Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.06		0.06		(b)
Index component value			0.06	0.06	0.06	0.06	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.42					
	P (mW)		2.9		2.9		#
	$P_{1 \times 1}$ (mW)		1.9		1.9		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.7					
	$z_{pii,\alpha}$ (cm)	0.7					
	f_{awf} (MHz)	6.11	6.10		6.10		#
Other information	pr_r (Hz)	3096					
	srr (Hz)	8.1					
	n_{pps}	14					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	7.5					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	1.1					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	1.6					
	p_r at z_{pii} (MPa)	0.49					
Operating controls	Exam type	Oph	Oph		Oph		
	Mode	CVD	CVD		CVD		
	2D optimization/depth (cm)	Pen/1.9	Pen/5.1		Pen/5.1		
	Color optimization/PRF (Hz)	Low/401	Med/4167		Med/4167		
	Color box position/size	Def/def	Top/short-wide		Top/short-wide		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 40: Transducer model: L25x (Ophthalmic use) Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.18	0.12		0.21		(b)
Index component value			0.12	0.08	0.12	0.21	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.44					
	P (mW)		4.0		4.0		#
	P_{1x1} (mW)		4.0		4.0		
	z_s (cm)			0.9			
	z_b (cm)					0.80	
	z_{MI} (cm)	1.2					
	$z_{pii,\alpha}$ (cm)	1.2					
Other information	f_{awf} (MHz)	6.03	6.03		6.03		#
	pr (Hz)	1953					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	7.4					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	18.4					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	44.9					
Operating controls	p_r at z_{pii} (MPa)	0.56					
	Exam type	Oph	Oph		Oph		
	Sample volume size (mm)	1	1		1		
	Sample volume position	Zone 7	Zone 7		Zone 7		
	PRF (Hz)	1953	5208		5208		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 41: Transducer model: L25x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.87					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
	f_{awf} (MHz)	6.11	#		#		#
Other information	pr (Hz)	1061					
	srr (Hz)	13.0					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	478					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	12.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	16.4					
	p_r at z_{pii} (MPa)	3.39					
Operating controls	Exam type	Nrv/Msk/ Ven/Art					
	Optimization	Any					
	Depth (cm)	1.9–2.2					
	Mbe	On					

(a) This index is not required for this operating mode; value is < 1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 42: Transducer model: L25x Operating mode: Color/CPD

Index label	MI	TIS		TIB		TIC
		At surface	Below surface	At surface	Below surface	At surface
Maximum index value	1.0	(a)		(a)		(b)
Index component value		#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.35				
	P (mW)		#	#	#	#
	$P_{1 \times 1}$ (mW)		#	#	#	
	z_s (cm)			—		
	z_b (cm)				—	
	z_{MI} (cm)	0.8				
	$z_{pii,\alpha}$ (cm)	0.8				
	f_{awf} (MHz)	6.11	#		#	#
Other information	p_{rr} (Hz)	5261				
	s_{rr} (Hz)	13.7				
	n_{pps}	14				
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	276				
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	81.5				
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	109.5				
	p_r at z_{pii} (MPa)	2.78				
Operating controls	Exam type	Ven				
	Mode	CVD				
	2D optimization/depth (cm)	Pen/3.1				
	Color optimization/PRF (Hz)	Low/779				
	Color box position/size	Def/def				

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 43: Transducer model: L25x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		1.7		(b)
Index component value			#	#	0.9	1.7	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		32.1		#
	P_{1x1} (mW)		#		32.1		
	z_s (cm)			#			
	z_b (cm)					0.75	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	#		#		#
Other information	pr (Hz)	#					
	srr (Hz)	#					
	n_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#					
p_r at z_{pii} (MPa)	#						
Operating controls	Exam type				Vas/Ven/Nrv		
	Sample volume size (mm)				8		
	Sample volume position				Zone 7		
	PRF (Hz)				1953		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 44: Transducer model: L38xi Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.3					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
Other information	f_{awf} (MHz)	4.82	#		#		#
	pr (Hz)	1312					
	srr (Hz)	10.3					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	605					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	10.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	13.5					
p_r at z_{pii} (MPa)	3.79						
Operating controls	Exam type	Nrv					
	Optimization	Res					
	Depth (cm)	2.0					
	MB	N/A					
	Needle vision	On					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 45: Transducer model: L38xi Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	(a)		1.2		(b)
Index component value			#	#	0.9	1.2	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.54					
	P (mW)		#		37.1		#
	P_{1x1} (mW)		#		37.1		
	z_s (cm)			#			
	z_b (cm)					0.9	
	z_{MI} (cm)	1.0					
	$z_{pii,\alpha}$ (cm)	1.0					
	f_{awf} (MHz)	5.76		#	5.20		#
Other information	pr (Hz)	1600					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	776					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sij,\alpha}$ (mW/cm^2)	181.8					
	I_{spta} at z_{pii} or z_{sij} (mW/cm^2)	280.5					
	p_r at z_{pii} (MPa)	4.32					
Operating controls	Exam type	Art			Art		
	Optimization	Gen			Pen		
	Depth (cm)	4.7			7.3		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 46: Transducer model: L38xi Operating mode: Color/CPD

Index label	MI	TIS		TIB		TIC	
		At surface	Below surface	At surface	Below surface	At surface	
Maximum index value	1.5	1.1		1.1		(b)	
Index component value		1.1	1.1	1.1	1.1		
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	3.3					
	P (mW)		64.7		64.7		#
	$P_{1 \times 1}$ (mW)		49.0		49.0		
	z_s (cm)			—			
	z_b (cm)				—		
	z_{MI} (cm)	0.8					
	$z_{pii,\alpha}$ (cm)	0.8					
	f_{awf} (MHz)	4.82	4.83		4.83		#
Other information	p_{rr} (Hz)	2190					
	s_{rr} (Hz)	4.5					
	n_{pps}	16					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	605					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	35.6					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	47.4					
	p_r at z_{pii} (MPa)	3.79					
Operating controls	Exam type	Art	Ven		Ven		
	Mode	CVD	CVD		CVD		
	2D optimization/depth (cm)	Pen/2.0	Pen/3.1		Pen/3.1		
	Color optimization/PRF (Hz)	Low/393	Low/2315		Low/2315		
	Color box position/size	Def/def	Bottom/short-narrow		Bottom/short-narrow		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 47: Transducer model: L38xi Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	2.6		3.7		(b)
Index component value			2.6	1.8	2.6	3.7	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.59					
	P (mW)		114.5		114.5		#
	P_{1x1} (mW)		114.5		114.5		
	z_s (cm)			1.2			
	z_b (cm)					1.2	
	z_{MI} (cm)	0.7					
	$z_{pii,\alpha}$ (cm)	0.7					
	f_{awf} (MHz)	4.06	4.78		4.78		#
Other information	pr (Hz)	1008					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	32.3					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	399.8					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	495.1					
p_r at z_{pii} (MPa)	2.86						
Operating controls	Exam type	Art	Nrv		Nrv		
	Sample volume size (mm)	1	1		1		
	Sample volume position	Zone 0	Zone 7		Zone 7		
	PRF (Hz)	1008	10417		10417		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 48: Transducer model: P10x Operating mode: Color

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		(a)		1.1
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		#		42.2
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	#		#		3.89
Other information	pr_r (Hz)	#					
	srr (Hz)	#					
	n_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	#					
	p_r at z_{pii} (MPa)	#					
Operating controls	Exam type						Crd
	Mode						CVD
	2D optimization/depth (cm)/sector width						Pen/8.9/narrow
	Color optimization/PRF (Hz)						Low/2033
	Color box position/size						Top/short-wide

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 49: Transducer model: P10x Operating mode: CW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		1.8		1.7
Index component value			#	#	0.7	1.8	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		34.8		25.7
	P_{1x1} (mW)		#		34.8		
	z_s (cm)			#			
	z_b (cm)					0.70	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	#		4.00		4.00
Other information	pr (Hz)	#					
	srr (Hz)	#					
	n_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sij,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sij} (mW/cm^2)	#					
Operating controls	Exam type				Crd		Crd
	Sample volume position				Zone 3		Zone 0

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 50: Transducer model: P10x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.0	1.1		1.9		1.5
Index component value			1.1	0.6	0.6	1.9	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.92					
	P (mW)		34.4		31.9		26.9
	P_{1x1} (mW)		34.4		31.9		
	z_s (cm)		1.4				
	z_b (cm)					0.90	
	z_{MI} (cm)	2.1					
	$z_{pii,\alpha}$ (cm)	2.1					
Other information	f_{awf} (MHz)	3.87	6.86		3.84		3.86
	pr (Hz)	1562					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	200					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	400.0					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	729.9					
p_r at z_{pii} (MPa)	2.54						
Operating controls	Exam type	Crd	Crd		Abd		Crd
	Sample volume size (mm)	1	7		12		1
	Sample volume position	Zone 2	Zone 6		Zone 1		Zone 0
	PRF (Hz)	1562	1008		1953		15625
	TDI	Off	On		Off		Off

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 51: Transducer model: rC60xi Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	(a)		(a)		(b)
Index component value			#	#	#	#	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.31					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	4.3					
	$z_{pii,\alpha}$ (cm)	4.3					
	f_{awf} (MHz)	2.36	#		#		#
Other information	pr (Hz)	3584					
	srr (Hz)	28.0					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	356					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	24.1					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	44.9					
p_r at z_{pii} (MPa)	3.29						
Operating controls	Exam type	Abd					
	Optimization	Res					
	Depth (cm)	11					
	MB (multi beam)	Off					
	THI	On					

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 52: Transducer model: rC60xi Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	(a)		1.0		(b)
Index component value			#	#	0.36	1.00	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.18					
	P (mW)		#		69.8		#
	P_{1x1} (mW)		#		25.9		
	z_s (cm)			#			
	z_b (cm)					4.2	
	z_{MI} (cm)	4.3					
	$z_{pii,\alpha}$ (cm)	4.3					
Other information	f_{awf} (MHz)	2.66	#		2.89		#
	pr_r (Hz)	800					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	290					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	144.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	328.2					
p_r at z_{pii} (MPa)	3.25						
Operating controls	Exam type	Abd			Msk		
	Optimization	Pen			Pen		
	Depth (cm)	6.6			9.2		
	THI	Off			Off		

(a) This index is not required for this operating mode; value is < 1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 53: Transducer model: rC60xi Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	1.2		1.2		(b)
Index component value			1.2	1.2	1.2	1.2	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.21					
	P (mW)		185.8		185.8		#
	$P_{1 \times 1}$ (mW)		107.5		107.5		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	4.3					
	$z_{pii,\alpha}$ (cm)	4.3					
	f_{awf} (MHz)	2.22	2.21		2.21		#
Other information	pr (Hz)	1265					
	srr (Hz)	9.89					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	342					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	8.9					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	15.8					
	p_r at z_{pii} (MPa)	3.07					
Operating controls	Exam type	Abd	Abd		Abd		
	Mode	CVD	CVD		CVD		
	2D optimization/depth (cm)/THI	Gen/11 / On	Gen/4.7 /Off		Gen/4.7 /Off		
	Color optimization/PRF (Hz)	Low/342	High/3125		High/3125		
	Color box position/size	Bottom/tall-narrow	Bottom/tall-narrow		Bottom/tall-narrow		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 54: Transducer model: rC60xi Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.2	2.0		4.0		(b)
Index component value			0.7	2.0	0.8	4.0	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.73					
	P (mW)		386.5		291.8		#
	P_{1x1} (mW)		67.5		74.2		
	z_s (cm)			4.0			
	z_b (cm)					3.6	
	z_{MI} (cm)	4.5					
	$z_{pii,\alpha}$ (cm)	4.5					
Other information	f_{awf} (MHz)	2.2	2.23		2.23		#
	pr (Hz)	1302					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	267					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	399.7					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	793.3					
Operating controls	p_r at z_{pii} (MPa)	2.43					
	Exam type	Abd	Abd		Abd		
	Sample volume size (mm)	3	7		7		
	Sample volume position	Zone 3	Zone 6		Zone 5		
	PRF (Hz)	1302	2604		2604		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 55: Transducer model: rP19x (Orbital use) Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.03		0.03		0.07
Index component value			0.03	0.03	0.03	0.03	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.25					
	P (mW)		4.4		4.4		4.7
	P_{1x1} (mW)		2.9		2.9		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	3.4					
	$z_{pii,\alpha}$ (cm)	3.4					
	f_{awf} (MHz)	2.06	2.06		2.06		1.90
Other information	pr (Hz)	6413					
	srr (Hz)	15.6					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	4.1					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	0.4					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	0.6					
	p_r at z_{pii} (MPa)	0.31					
Operating controls	Exam type	Orb	Orb		Orb		Orb
	Optimization	Res	Res		Res		Gen
	Depth (cm)	4.7	4.7		4.7		16
	MB	Off	Off		Off		Off

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 56: Transducer model: rP19x (Orbital use) Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.009		0.020		0.021
Index component value			0.006	0.009	0.006	0.020	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.25					
	P (mW)		1.34		1.34		1.34
	P_{1x1} (mW)		0.67		0.67		
	z_s (cm)			2.5			
	z_b (cm)					3.15	
	z_{MI} (cm)	3.4					
	$z_{pii,\alpha}$ (cm)	3.4					
	f_{awf} (MHz)	2.06	1.83		1.83		1.83
Other information	pr (Hz)	800					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	4.05					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	1.7					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	2.7					
	p_r at z_{pii} (MPa)	0.31					
Operating controls	Exam type	Orb	Orb		Orb		Orb
	Optimization	Res	Gen		Gen		Gen
	Depth (cm)	4.7	35		35		35

(a) This index is not required for this operating mode; value is <1 .

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 57: Transducer model: rP19x (Orbital use) Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.17	0.09		0.09		0.23
Index component value			0.09	0.09	0.09	0.09	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.25					
	P (mW)		15.47		15.47		15.50
	P_{1x1} (mW)		9.50		9.50		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	0.7					
	$z_{pii,\alpha}$ (cm)	0.7					
	f_{awf} (MHz)	2.14	2.11		2.11		2.11
Other information	pr (Hz)	5443					
	srr (Hz)	15.9					
	n_{pps}	16					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	1.82					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	3.2					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	3.5					
	p_r at z_{pii} (MPa)	0.26					
Operating controls	Exam type	Orb	Orb		Orb		Orb
	Mode	CVD	CVD		CVD		CVD
	2D optimization/depth (cm)	Gen/4.7	Gen/24		Gen/24		Gen/24
	Color optimization/PRF (Hz)	Low/1157	Low/3125		Low/3125		Low/3125
	Color box position/size	Def/def	Top/short-wide		Top/short-wide		Top/short-wide

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 58: Transducer model: rP19x (Orbital use) Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		0.18	0.27		0.59		0.57
Index component value			0.19	0.27	0.18	0.59	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	0.27					
	P (mW)		37.4		35.3		37.4
	P_{1x1} (mW)		17.5		17.0		
	z_s (cm)			2.5			
	z_b (cm)					3.35	
	z_{MI} (cm)	3.5					
	$z_{pii,\alpha}$ (cm)	3.5					
	f_{awf} (MHz)	2.23	2.23		2.23		2.23
Other information	pr (Hz)	1953					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	2.49					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	28.9					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	69.3					
	p_r at z_{pii} (MPa)	0.36					
Operating controls	Exam type	Orb	Orb		Orb		Orb
	Sample volume size (mm)	5	14		14		14
	Sample volume position	Zone 6	Zone 7		Zone 5		Zone 7
	PRF (Hz)	1953	1953		1953		1953

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 59: Transducer model: rP19x Operating mode: 2D

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	1.0		1.0		2.7
Index component value			1.0	1.0	1.0	1.0	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.1					
	P (mW)		152.6		152.6		177.8
	P_{1x1} (mW)		96.1		96.1		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	4.8					
	$z_{pii,\alpha}$ (cm)	4.8					
	f_{awf} (MHz)	1.99	2.08		2.08		1.53
Other information	pr (Hz)	6186					
	srr (Hz)	48.3					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	184					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	25.4					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	38.6					
	p_r at z_{pii} (MPa)	2.92					
Operating controls	Exam type	Abd	Crd		Crd		Crd
	Optimization	Gen	Res		Res		Pen
	Depth (cm)	10	10		10		4.7
	MB/THI	Off/Off	Off/On		Off/On		Off/On
	Sector width	N/A	Narrow		Narrow		N/A

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 60: Transducer model: rP19x Operating mode: M mode

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	(a)		1.7		1.0
Index component value			#	#	0.2	1.7	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.1					
	P (mW)		#		55.0		62.1
	P_{1x1} (mW)		#		28.5		
	z_s (cm)			#			
	z_b (cm)					4.33	
	z_{MI} (cm)	4.8					
	$z_{pii,\alpha}$ (cm)	4.8					
Other information	f_{awf} (MHz)	1.99	#		1.81		1.77
	pr (Hz)	800					
	srr (Hz)	—					
	η_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	184					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	73.5					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	140.8					
Operating controls	p_r at z_{pii} (MPa)	2.92					
	Exam type	TCD			Abd		Abd
	Optimization	Gen			Res		Res
	Depth (cm)	7.5			10		16
	THI	Off			On		On

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 6 1: Transducer model: rP19x Operating mode: Color/CPD

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.5	1.2		1.2		2.5
Index component value			1.2	1.2	1.2	1.2	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	2.1					
	P (mW)		128.0		128.0		170.5
	P_{1x1} (mW)		115.6		115.6		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	4.8					
	$z_{pii,\alpha}$ (cm)	4.8					
	f_{awf} (MHz)	1.99	2.14		2.14		2.12
Other information	pr (Hz)	505					
	srr (Hz)	7.9					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	184					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	2.1					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	3.2					
	p_r at z_{pii} (MPa)	2.92					
Operating controls	Exam type	Abd	TCD		TCD		Crd
	Mode/THI	CVD/Off	CVD/Off		CVD/Off		CVD/On
	2D optimization/depth (cm)/sector width	Gen/10/N/A	Pen/7.5 /N/A		Pen/7.5 /N/A		Gen/16/narrow
	Color optimization/PRF (Hz)	Low/300	Low/3125		Low/3125		High/5208
	Color box position/size	Def/def	Def/narrow		Def/narrow		Def/def

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 62: Transducer model: rP19x Operating mode: CW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	1.2		4.0		4.0
Index component value			1.2	1.1	1.2	4.0	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		125.4		125.4		125.4
	P_{1x1} (mW)		125.4		125.4		
	z_s (cm)			0.9			
	z_b (cm)					0.9	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	2.00		2.00		2.00
Other information	pr_r (Hz)	#					
	srr (Hz)	#					
	η_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#					
	p_r at z_{pii} (MPa)	#					
Operating controls	Exam type		Crd		Crd		Crd
	Sample volume position		Zone 0		Zone 0		Zone 0

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 63: Transducer model: rP19x Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		1.3	1.8		4.0		3.9
Index component value			1.3	1.8	1.2	4.0	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	1.94					
	P (mW)		253.7		240.2		251.1
	$P_{1 \times 1}$ (mW)		118.6		116.0		
	z_s (cm)			2.5			
	z_b (cm)					3.35	
	z_{MI} (cm)	3.0					
	$z_{pii,\alpha}$ (cm)	3.0					
	f_{awf} (MHz)	2.14	2.23		2.23		2.10
Other information	pr (Hz)	1562					
	srr (Hz)	—					
	n_{pps}	1					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm ²)	180					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm ²)	374.9					
	I_{spta} at z_{pii} or z_{sii} (mW/cm ²)	594.7					
	p_r at z_{pii} (MPa)	2.42					
Operating controls	Exam type	Crd	Crd		Crd		Crd
	Sample volume size (mm)	1	12		1		1
	Sample volume position	Zone 1	Zone 7		Zone 5		Zone 5
	PRF (Hz)	1562	1562		39062		39062
	TDI	Off	Off		Off		Off

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 64: Transducer model: TEExi Operating mode: CW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		1.7		(b)
Index component value			#	#	0.7	1.7	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		34.4		#
	P_{1x1} (mW)		#		34.4		
	z_s (cm)			#			
	z_b (cm)					1.10	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
Other information	f_{awf} (MHz)	#	#		4.00		#
	pr (Hz)	#					
	srr (Hz)	#					
	η_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#					
Operating controls	p_r at z_{pii} (MPa)	#					
	Exam type				Crd		
	Sample volume				Zone 2		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Table 65: Transducer model: TEExi Operating mode: PW Doppler

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	At surface
Maximum index value		(a)	(a)		1.4		(b)
Index component value			#	#	0.7	1.4	
Acoustic parameters	$p_{r,\alpha}$ at z_{MI} (MPa)	#					
	P (mW)		#		35.8		#
	P_{1x1} (mW)		#		35.8		
	z_s (cm)			#			
	z_b (cm)					2.57	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	#		3.81		#
Other information	pr (Hz)	#					
	srr (Hz)	#					
	n_{pps}	#					
	$I_{pa,\alpha}$ at $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ at $z_{pii,\alpha}$ or $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} at z_{pii} or z_{sii} (mW/cm^2)	#					
Operating controls	p_r at z_{pii} (MPa)	#					
	Exam type				Crd		
	Sample volume size (mm)				1		
	Sample volume position				Zone 3		
	PRF (Hz)				2604		

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference global maximum index value line.)

— Data are not applicable for this transducer/mode.

Terms used in the acoustic output tables

Table 66: Terms used in the acoustic output tables

Term	Definition
α	Attenuation coefficient used for derating. Equal to 0.3 dB/cm/MHz ² .
f_{awf}	Acoustic working frequency.
$I_{pa,\alpha}$	Attenuated pulse-average intensity.
I_{spta}	Spatial-peak temporal-average intensity.
$I_{spta,\alpha}$	Attenuated spatial-peak temporal-average intensity.
MI	Mechanical index.
P	Output power.
$P_{1 \times 1}$	Bounded-square output power.
$p_{r,\alpha}$	Attenuated peak-rarefactional acoustic pressure.
p_r	Peak-rarefactional acoustic pressure.
p_{ii}	Pulse-intensity integral.
$p_{ii,\alpha}$	Attenuated pulse-intensity integral.
n_{pps}	Number of pulses per ultrasonic scan line.
prr	Pulse repetition rate.
srr	Scan repetition rate.
TI	Thermal index.
TIB	Bone thermal index.
TIC	Cranial-bone thermal index.
TIS	Soft-tissue thermal index.
z_b	Depth for TIB.
z_{MI}	Depth for mechanical index.
$z_{p_{ii}}$	Depth for peak pulse-intensity integral.
$z_{p_{ii,\alpha}}$	Depth for peak attenuated pulse-intensity integral.

Table 66: Terms used in the acoustic output tables

Term	Definition
z_{sij}	Depth for peak sum of pulse-intensity integrals.
$z_{sij,\alpha}$	Depth for peak sum of attenuated pulse-intensity integrals.
z_s	Depth for TIS.

Glossary (SonoSite Edge II)

The IMT term has been removed from the SonoSite Edge II user guide abbreviation list; the revision will be made in the next update.

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Table 7: Marking symbols for standards (continued)

Symbol	Title	Standard development organization	Reference number	Description
	Defibrillator-safe used parts of type CF	IEC 60601-1 Electromedical equipment – Part 1: General safety requirements and essential performance characteristics	D.1-27	Indicates defibrillator-safe used parts of type BF, which comply with IEC 60601-1
	Units that are sensitive to electrostatic discharge	IEC 60417:2002 Graphic symbols for use on equipment	5134	Indicates packages containing units sensitive to electrostatic discharge, or identifies a unit or part that has not been tested for immunity to electrostatic discharge
	Regulatory compliance mark (RCM)	AS/NZS3820	–	Indicates C-Tick – regulatory compliance mark for Australia and New Zealand The unit is in compliance with the relevant Australian and New Zealand requirements for electronic equipment.
	Batch code, date code or part code, control type number	ISO 15223-1 Medical devices – Symbols to be used on labels for medical devices, marking and information to be provided – Part 1: General requirements	5.1.5	Indicates the manufacturer's batch code, so that the responsible batch or lot can be identified
	Biological hazard	ISO 7010 – Graphic symbols – Safety colors and safety signs	W009	Warns of a biological hazard

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Transducermodel: ICTx Billeddannelsestype: PW-doppler	152
Transducermodel: L25x (oftalmisk anvendelse) Billeddannelsestype: 2D	153
Transducermodel: L25x (oftalmisk anvendelse) Billeddannelsestype: M Mode	154
Transducermodel: L25x (oftalmisk anvendelse) Billeddannelsestype: Color/CPD	155
Transducermodel: L25x (oftalmisk anvendelse) Billeddannelsestype: PW-doppler	156
Transducermodel: L25x Billeddannelsestype: 2D	157
Transducermodel: L25x Billeddannelsestype: Color/CPD	158
Transducermodel: L25x Billeddannelsestype: PW-doppler	159
Transducermodel: L38xi Billeddannelsestype: 2D	160
Transducermodel: L38xi Billeddannelsestype: M Mode	161
Transducermodel: L38xi Billeddannelsestype: Color/CPD	162
Transducermodel: L38xi Billeddannelsestype: PW-doppler	163
Transducermodel: P10x Billeddannelsestype: Color	164
Transducermodel: P10x Billeddannelsestype: CW-doppler	165
Transducermodel: P10x Billeddannelsestype: PW-doppler	166
Transducermodel: rC60xi Billeddannelsestype: 2D	167
Transducermodel: rC60xi Billeddannelsestype: M Mode	168
Transducermodel: rC60xi Billeddannelsestype: Color/CPD	169
Transducermodel: rC60xi Billeddannelsestype: PW-doppler	170
Transducermodel: rP19x (orbital anvendelse) Billeddannelsestype: 2D	171
Transducermodel: rP19x (orbital anvendelse) Billeddannelsestype: M Mode	172
Transducermodel: rP19x (orbital anvendelse) Billeddannelsestype: Color/CPD	173
Transducermodel: rP19x (orbital anvendelse) Billeddannelsestype: PW-doppler	174
Transducermodel: rP19x Billeddannelsestype: 2D	175
Transducermodel: rP19x Billeddannelsestype: M Mode	176
Transducermodel: rP19x Billeddannelsestype: Color/CPD	177
Transducermodel: rP19x Billeddannelsestype: CW-doppler	178
Transducermodel: rP19x Billeddannelsestype: PW-doppler	179
Transducermodel: TEExi Billeddannelsestype: CW-doppler	180
Transducermodel: TEExi Billeddannelsestype: PW-doppler	181

Tabel 10: Transducermodel: C8x Billeddannelsestype: 2D

Indeksbetegnelse		MI	TIS		TIB		TIC
			På overfladen	Under overfladen	På overfladen	Under overfladen	På overfladen
Maksimal indekssværdi		1,1	(a)		(a)		(b)
Indekskomponentværdi			#	#	#	#	
Akustiske parametre	$p_{r,\alpha}$ ved z_{MI} (MPa)	2,48					
	P (mW)		#		#		#
	P_{1x1} (mW)		#		#		
	z_s (cm)			—			
	z_b (cm)					—	
	z_{MI} (cm)	1,2					
	$z_{pii,\alpha}$ (cm)	1,2					
	f_{awf} (MHz)	5,53	#		#		#
Andre oplysninger	p_{rr} (Hz)	9524					
	s_{rr} (Hz)	18,6					
	η_{pps}	1					
	$I_{pa,\alpha}$ ved $z_{pii,\alpha}$ (W/cm^2)	264					
	$I_{spta,\alpha}$ ved $z_{pii,\alpha}$ eller $z_{sii,\alpha}$ (mW/cm^2)	18,3					
	I_{spta} ved z_{pii} eller z_{sii} (mW/cm^2)	25,6					
	p_r ved z_{pii} (MPa)	3,11					
Drifts-kontroller	Undersøgelsestype	Pro					
	Optimering	Pen					
	Dybde (cm)	2,5 - 3,2					
	MB	Fra					

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— Der findes ingen data for denne transducer/tilstand.

Tabel 11: Transducermodel: C8x Billeddannelsestype: M Mode

Indeksbetegnelse		MI	TIS		TIB		TIC
			På overfladen	Under overfladen	På overfladen	Under overfladen	På overfladen
Maksimal indekssværdi		1,3	(a)		(a)		(b)
Indekskomponentværdi			#	#	#	#	
Akustiske parametre	$p_{r,\alpha}$ ved z_{MI} (MPa)	2,91					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			#			
	z_b (cm)					#	
	z_{MI} (cm)	1,1					
	$z_{pii,\alpha}$ (cm)	1,1					
	f_{awf} (MHz)	5,07	#		#		#
Andre oplysninger	p_{rr} (Hz)	800					
	s_{rr} (Hz)	–					
	n_{pps}	1					
	$I_{pa,\alpha}$ ved $z_{pii,\alpha}$ (W/cm^2)	433					
	$I_{spta,\alpha}$ ved $z_{pii,\alpha}$ eller $z_{sii,\alpha}$ (mW/cm^2)	149					
	I_{spta} ved z_{pii} eller z_{sii} (mW/cm^2)	226					
	p_r ved z_{pii} (MPa)	3,57					
Drifts-kontroller	Undersøgelsestype	Pro					
	Optimering	Pen					
	Dybde (cm)	4,2					

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— Der findes ingen data for denne transducer/tilstand.

Tabel 12: Transducermodel: C8x Billeddannelsestype: Color/CPD

Indeksbetegnelse		MI	TIS		TIB		TIC
			På overfladen	Under overfladen	På overfladen	Under overfladen	På overfladen
Maksimal indekssværdi		1,2	(a)		(a)		(b)
Indekskomponentværdi			#	#	#	#	
Akustiske parametre	$p_{r,\alpha}$ ved z_{MI} (MPa)	2,68					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			–			
	z_b (cm)					–	
	z_{MI} (cm)	0,8					
	$z_{pii,\alpha}$ (cm)	0,8					
	f_{awf} (MHz)	4,82	#		#		#
Andre oplysninger	p_{rr} (Hz)	2548					
	s_{rr} (Hz)	26					
	n_{pps}	12					
	$I_{pa,\alpha}$ ved $z_{pii,\alpha}$ (W/cm ²)	381					
	$I_{spta,\alpha}$ ved $z_{pii,\alpha}$ eller $z_{sij,\alpha}$ (mW/cm ²)	132					
	I_{spta} ved z_{pii} eller z_{sij} (mW/cm ²)	176					
	p_r ved z_{pii} (MPa)	3,1					
Drifts-kontroller	Undersøgelsestype	Pro					
	Type	CVD					
	2D-optimering/dybde (cm)	Pen/1,5 – 1,9					
	Farveoptimering/PRF (Hz)	Høj/alle					
	Farveboksposition/-størrelse	Smal/alle					

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— Der findes ingen data for denne transducer/tilstand.

Tabel 13: Transducermodel: C8x Billeddannelsestype: PW-doppler

Indeksbetegnelse		MI	TIS		TIB		TIC
			På overfladen	Under overfladen	På overfladen	Under overfladen	På overfladen
Maksimal indeksværdi		1,0	(a)		1,4		(b)
Indekskomponentværdi			#	#	0,5	1,4	
Akustiske parametre	$p_{r,\alpha}$ ved z_{MI} (MPa)	2,28					
	P (mW)		#		23,1		#
	$P_{1 \times 1}$ (mW)		#		23,1		
	z_s (cm)			#			
	z_b (cm)					1,0	
	z_{MI} (cm)	1,8					
	$z_{pii,\alpha}$ (cm)	1,8					
	f_{awf} (MHz)	4,80	#		4,80		#
Andre oplysninger	pr_r (Hz)	1008					
	srr (Hz)	–					
	n_{pps}	1					
	$I_{pa,\alpha}$ ved $z_{pii,\alpha}$ (W/cm^2)	263					
	$I_{spta,\alpha}$ ved $z_{pii,\alpha}$ eller $z_{sii,\alpha}$ (mW/cm^2)	334					
	I_{spta} ved z_{pii} eller z_{sii} (mW/cm^2)	616					
	p_r ved z_{pii} (MPa)	3,1					
Drifts-kontroller	Undersøgelsestype	Pro			Pro		
	Prøvevolumenstørrelse (mm)	1			1		
	Prøvevolumenposition	Zone 4			Zone 4		
	PRF (Hz)	1008			1008		

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— Der findes ingen data for denne transducer/tilstand.

Tabel 15: Transducermodel: C35x Billeddannelsestype: 2D

Indeksbetegnelse		MI	TIS		TIB		TIC
			På overfladen	Under overfladen	På overfladen	Under overfladen	På overfladen
Maksimal indeksværdi		1,0	(a)		(a)		(b)
Indekskomponentværdi			#	#	#	#	
Akustiske parametre	$p_{r,\alpha}$ ved z_{MI} (MPa)	1,8					
	P (mW)		#		#		#
	$P_{1 \times 1}$ (mW)		#		#		
	z_s (cm)			–			
	z_b (cm)					–	
	z_{MI} (cm)	3,3					
	$z_{pii,\alpha}$ (cm)	3,3					
	f_{awf} (MHz)	3,45	#		#		#
Andre oplysninger	p_{rr} (Hz)	1021					
	s_{rr} (Hz)	7,98					
	η_{pps}	1					
	$I_{pa,\alpha}$ ved $z_{pii,\alpha}$ (W/cm^2)	250					
	$I_{spta,\alpha}$ ved $z_{pii,\alpha}$ eller $z_{sii,\alpha}$ (mW/cm^2)	8,6					
	I_{spta} ved z_{pii} eller z_{sii} (mW/cm^2)	16,5					
	p_r ved z_{pii} (MPa)	2,61					
Drifts-kontroller	Undersøgelsestype	Msk					
	Optimering	Res					
	Dybde (cm)	8,3					
	MB	Ikke relevant					

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— Der findes ingen data for denne transducer/tilstand.

Tabel 16: Transducermodel: C35x Billeddannelsestype: PW-doppler

Indeksbetegnelse		MI	TIS		TIB		TIC
			På overfladen	Under overfladen	På overfladen	Under overfladen	På overfladen
Maksimal indekssværdi		(a)	1,5		2,6		(b)
Indekskomponentværdi			1,5	1,0	1,0	2,6	
Akustiske parametre	$p_{r,\alpha}$ ved z_{MI} (MPa)	#					
	P (mW)		72,8		47,1		#
	$P_{1 \times 1}$ (mW)		71,1		47,1		
	z_s (cm)			1,4			
	z_b (cm)					0,50	
	z_{MI} (cm)	#					
	$z_{pii,\alpha}$ (cm)	#					
	f_{awf} (MHz)	#	4,35		4,37		#
Andre oplysninger	p_{rr} (Hz)	#					
	s_{rr} (Hz)	#					
	η_{pps}	#					
	$I_{pa,\alpha}$ ved $z_{pii,\alpha}$ (W/cm^2)	#					
	$I_{spta,\alpha}$ ved $z_{pii,\alpha}$ eller $z_{sii,\alpha}$ (mW/cm^2)	#					
	I_{spta} ved z_{pii} eller z_{sii} (mW/cm^2)	#					
	p_r ved z_{pii} (MPa)	#					
Drifts-kontroller	Undersøgelsestype		Rygsøjle		Rygsøjle		
	Prøvevolumenstørrelse (mm)		2		1		
	Prøvevolumenposition		Zone 5		Zone 0		
	PRF (Hz)		6250		15.625		

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