SonoSite SII

User Guide Errata

Manufacturer	EC Authorized Representative	Australia Sponsor
FUJIFILM SonoSite, Inc. 21919 30th Drive SE Bothell, WA 98021 USA T: 1-888-482-9449 or 1-425-951-1200 F: 1-425-951-1201	FUJIFILM SonoSite B.V. Joop Geesinkweg 140 1114 AB Amsterdam, The Netherlands	FUJIFILM SonoSite Australasia Pty Ltd 114 Old Pittwater Road BROOKVALE, NSW, 2100 Australia

Caution

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite SII, SonoSite, and the SonoSite logo are registered and unregistered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions.

All other trademarks are the property of their respective owners.

Part Number: P22338-04

Publication Date: November 2017

Copyright © 2017 FUJIFILM SonoSite, Inc. All Rights reserved.

SonoSite SII Ultrasound System User Guide

Errata

Introduction	1
Document conventions	1
Getting Help	2
Saving images and clips	2
Getting Started	3
System Controls	3
Imaging	3
Freezing, viewing frames, and zooming	3
Images and clips	4
Measurements and calculations	4
Calculations menu	4
Cardiac calculations	5
Gynecology (Gyn) calculations	8
OB calculations	9
Cleaning and disinfecting	10
Safety	11
Hazardous materials	11
IT Network	11
Functions	11
Network for connecting the device	11
Specifications for the connection	11

Introduction

Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ♦.

For a description of labeling symbols that appear on the product, see "Labeling Symbols" in the user guide.

Getting Help

For technical support, please contact FUJIFILM SonoSite as follows:

Phone (U.S. or Canada)	877-657-8118
Phone (outside U.S. or Canada)	425-951-1330, or call your local representative
Fax	425-951-6700
Email	ffss-service@fujifilm.com
Web	www.sonosite.com
Europe Service Center	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
Asia Service Center	+65 6380-5581

Printed in the U.S.

The following content was either not included or incorrect in the system user guide; the revisions will be made in the next update.


Saving images and clips

The percentage icon in the system status area shows the percentage of space **available** in internal storage.

Getting Started

System Controls

- | | | |
|---|------------------------|--|
| 1 | Control knobs | Turn to adjust gain, depth, cine buffer, brightness, and more, depending on context. Current functions appear on-screen above the knobs. |
| 2 | Freeze key | Press and hold to freeze or unfreeze the image. |
| 3 | Touchpad | Moves the pointer and other items. |
| 4 | Touchpad key | Works in conjunction with the touchpad. Tap to activate an item on-screen, or to switch between color box functions. (active only when the image is frozen.) |
| 5 | Print key | Available only when a printer is connected to the system. Tap to print from a live or frozen scan. |
| 6 | Save keys | Tap one of these keys to save an image or a clip. |
| 7 | Image mode | Tap one of these keys to change the imaging mode. |
| 8 | System controls | Change system settings, switch transducers, add labels, or see patient information. |
| 9 | Image controls | Use these to adjust the image. |

Imaging

Freezing, viewing frames, and zooming

To freeze or unfreeze an image

❖ Press and hold the Freeze button (▶◀).

When the image is frozen, the button color is blue. When the image is unfrozen, the button color is white.

On a frozen image, the cine icon () and frame number appear above the left knob.

To move forward or backward in the cine buffer

- ❖ On a frozen image, do one of the following:
 - ▶ Turn the Cine (left) knob
 - ▶ Drag your finger left or right on the touchscreen
 - ▶ Drag your finger left or right on the touchpad

The total number of frames appears next to the cine icon. The number changes to the current frame number as you move forward or backward.

WARNING

To avoid loss of data, be careful not to touch the Freeze button () while turning the Cine knob.

Images and clips

Printing, exporting, and deleting images and clips

To manually archive images and clips

You can send patient exams to a DICOM printer or archiver, or to a PC using SonoSite Patient Data Archival Software (PDAS). DICOM and SonoSite PDAS are optional features. For more information about archiving, refer to the SonoSite PDAS and DICOM documentation.

Measurements and calculations

Calculations menu

To select from the calculations menu

- 1 On a frozen image, tap **Calcs**. The calculations menu appears.

The list of calculations or measurements can be too long to fit on a page. To see the next page of calculations or measurements, tap **Next**. To see the previous page, tap **Previous**.

- 2 To start a calculation, tap the desired calculation name.

Note

Only calculations and measurements that are compatible with the current imaging mode are shown.

Many calculations include more than one measurement. The measurements for each calculation appear under the calculation name. You can perform the measurements in any order.

- 3 To perform a measurement within a calculation, tap the measurement name.

- 4 To close the calculations menu, tap **Back**.

Tapping **Back** will not save your calculation.

Cardiac calculations

To measure Ao, LA, AAo, or LVOT D

- 1 On a frozen 2D image or M Mode trace, tap **Calcs**.

- 2 From the calculations menu, tap **Ao/LA**.

- 3 From the **Ao/LA** menu, select the measurement you want to take.

- 4 Position the calipers by dragging.

- 5 Tap **Save Calc.**

- 6 To save a picture of the finished calculation, tap .

- 7 Tap **Back** to exit the calculation.

To measure TAPSE

- 1 On a frozen M Mode trace, tap **Calcs**.

- 2 In the calculations menu, tap **TAPSE**.

- 3 Position the calipers by dragging.

- 4 Tap **Save Calc.**

- 5 To save a picture of the finished calculation, tap .

- 6 Tap **Back** to exit the calculation.

To calculate MV or AV area

- 1 On a frozen 2D image, tap **Calcs**.

- 2 In the calculations menu, tap **Area**.

- 3 In the **Area** menu, select **MV** or **AV**.

- 4** If you are using the touchscreen, lift your finger from the screen momentarily to activate the trace. If you are using the touchpad, tap  to set the position.

The trace feature becomes active.

- 5** Using the touchpad or the touchscreen, trace the desired area.

To make a correction, tap **Undo**, or tap the measurement name to restart the measurement.

WARNING

When using the touchpad to trace a shape, be careful not to touch  until you are finished with the trace. Doing so may complete the trace prematurely, causing an incorrect measurement and delay of care.

- 6** When you are done, tap **Set**. The two ends of the trace are joined automatically.

- 7** Tap **Save Calc.**

- 8** To save a picture of the finished calculation, tap .

- 9** Tap **Back** to exit the calculation.

To calculate LV mass

- 1** On a frozen 2D image, tap **Calcs**.

- 2** In the calculations menu, tap **LV mass**.

- 3** Do the following for these cardiac measurements, **EPI** and **Endo**:

- a** Select the measurement name from the **LV mass** menu.

- b** Using the touchpad or the touchscreen, position the caliper where you want to begin the trace.

- c** If you are using the touchscreen, lift your finger from the screen momentarily to activate the trace. If

you are using the touchpad, tap  to activate the trace.

To make a correction, tap **Undo**, or tap the measurement name to restart the measurement.

- d** Using the touchpad or the touchscreen, trace the desired area, and then tap **Set**.

WARNING

When using the touchpad to trace a shape, be careful not to touch  until you are finished with the trace. Doing so may complete the trace prematurely, causing an incorrect measurement and delay of care.

- e** Tap **Save Calc.**

A check mark appears next to the saved measurement.

4 Select **Apical** from the **LV mass** menu.

5 Positioning the calipers, measure the ventricular length.

6 Tap **Save Calc.**

7 To save a picture of the finished calculation, tap .

8 Tap **Back** to exit the calculation.

To measure LVd and LVs

1 On a frozen 2D image or M Mode trace, tap **Calcs**.

2 Tap **LVd** or **LVs**.

3 Repeat the following for each measurement you want to take:

a On the **LV** calculation list, tap the measurement you want to take.

b Position the calipers by dragging.

c Tap **Save Calc.**

A check mark appears next to the saved measurement.

4 To save a picture of the finished calculation, tap .

5 Tap **Back** to exit the calculation.

To measure IVC Collapse

1 On a frozen 2D or M Mode image, tap **Calcs**.

2 Tap **IVC**.

3 Do the following for both **Max D** and **Min D** measurements.

a On the **IVC** calculations list, tap the measurement you want to take

b Position the calipers by dragging.

c Tap **Save Calc.**

A check mark appears next to the saved measurement.

4 To save a picture of the finished calculation, tap .

5 Tap **Back** to exit the calculation.

To measure EF

1 On a frozen M Mode trace, tap **Calcs**.

2 Tap **EF**.

3 Do the following for both **LVDd** and **LVDs** measurements.

- a** On the **EF** calculations list, tap the measurement that you want to take.
- b** Position the calipers by dragging.
- c** Tap **Save Calc.**

A check mark appears next to the saved measurement.

4 To save a picture of the finished calculation, tap .

5 Tap **Back** to exit the calculation.

Gynecology (Gyn) calculations

To measure uterus or ovary

1 On a frozen 2D image, tap **Calcs.**

2 Tap the name of the structure you want to measure: **Uterus**, **R Ovary**, or **L Ovary**.

3 Do the following for each of the length, height, and width measurements:

- a** Select the measurement name from the calculations menu.
- b** Position the calipers by dragging.
- c** Tap **Save Calc.**

A check mark appears next to the saved measurement.

4 To save a picture of the finished calculation, tap .

5 Tap **Back** to exit the calculation.

To measure follicles

On each side, you can save up to three distance measurements on a follicle, for up to 10 follicles. If you measure a follicle twice, the average appears in the report. If you measure a follicle three times, the average and a volume calculation appear in the report.

1 On a frozen 2D image, tap **Calcs.**

2 Tap the name of the structure you want to measure: **R Follicle** or **L Follicle**.

3 Do the following for each follicle you want to measure:

- a** Tap a measurement number.
- b** Position the calipers by dragging.
- c** Tap **Save Calc.**

A check mark appears next to the saved measurement.

4 To save a picture of the finished calculation, tap .

5 Tap **Back** to exit the calculation.

OB calculations

To measure gestational growth (2D)

For each 2D OB measurement (except CxLen and YS), the system saves up to three individual measurements and their average. If you take more than three measurements, the earliest measurement is deleted.

1 In the patient information form, select **OB** exam type, and enter the **LMP** or **Estab.DD** for the patient, if known.

2 On a frozen 2D image, tap **Calcs**.

3 Do the following for each measurement you want to take:

a From the calculations menu, select the measurement name.

Note

The caliper tool may change depending on the measurement selected.

b Position the calipers by dragging.

c Tap **Save Calc**.

A check mark appears next to the saved measurement.

4 To save a picture of the finished calculation, tap .

5 Tap **Back** to exit the calculation.

To measure fetal heart rate (M Mode)

1 On a frozen M Mode trace, tap **Calcs**.

2 Select **FHR** from the calculations menu. A vertical caliper appears.

3 Position the first vertical caliper at the start of the first heartbeat.

4 Position the second vertical caliper at the start of the next heartbeat.

5 Tap **Save Calc**.

6 To save a picture of the finished calculation, tap .

7 Tap **Back** to exit the calculation.

Cleaning and disinfecting

Table 1: Approved cleaners/disinfectants

Product	Compatible Transducer	Minimum wet contact time ¹
Sani-Cloth AF3 ²	C8x, C11x, C35x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, P10x, P11x ³ , rC60xi, rP19x	3 minutes
Sani-Cloth Plus	HFL50x, L52x ⁴	3 minutes
PI-Spray II	C8x, C11x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, L52x ⁴ , P10x, P11x ³ , rC60xi, rP19x	10 minutes

¹For maximum effectiveness, the component being cleaned must remain wet with disinfectant for a minimum period of time.

²Qualified for use as an intermediate-level disinfectant for mycobacteria.

³The P11x Transducer is not licensed for use in Canada.

⁴The L52x Transducer is for veterinary use only.

Refer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.

Table 2: High-level compatible disinfectants

Disinfectant	Compatible transducers	Temperature	Disinfectant Soak Duration
Cidex	C8x, C11x, C35x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, P10x, P11x ¹ , rC60xi, rP19x	25°C, 77°F	45 minutes
Cidex OPA	C8x, C11x, HFL38xi, HFL50x, HSL25x, ICTx, L25x, L38xi, L52x ² , P10x, rC60xi, rP19x	20°C, 68°F	12 minutes

¹The P11x Transducer is not licensed for use in Canada.

²The L52x Transducer is for veterinary use only.

Refer to the cleaners and disinfection document available at www.sonosite.com for a complete list of the most current cleaners and disinfectants.

Cleaning the ECG cable

This section is no longer applicable and has been removed.

Safety

Hazardous materials

WARNING

Products and accessories may contain hazardous materials. Ensure that products and accessories are disposed of in an environmentally responsible manner and meet federal and local regulations for disposing of hazardous materials.

IT Network

Functions

This device can be connected to an IT network to perform the following functions:

- ▶ Storing the examination data (static images, clips) acquired by this device in Picture Archiving and Communication System (PACS) by DICOM communication.
- ▶ Querying examination orders from the Modality Worklist (MWL) server by DICOM communication and start them.
- ▶ Upload images to the PDAS or SiteLink Image Manager.
- ▶ Setting the time of this device by inquiring the network time service.
- ▶ Communicating procedure status via the Modality Performed Procedure Step (MPPS) service.
- ▶ Requesting transfer of responsibility for image ownership to another system via the Storage Commitment service.

Network for connecting the device

To ensure safety, use an IT network that is isolated from the external environment by a firewall.

Specifications for the connection

Hardware specification

- ▶ 802.11 b/g/n
- ▶ Ethernet 100BASE-TX/10BASE-T using RJ45 port with patch cable

Software specifications

- ▶ This device is connected to PACS and MWL by DICOM standard. Refer to the DICOM Conformance Statement of this device for details.
- ▶ When available, this device connects to the network time server at startup.

- The system conforms to the DICOM standard as specified in the *SonoSite Edge II, SonoSite SII DICOM Conformance Statement*, available at www.sonosite.com. This statement provides information about the purpose, characteristics, configuration, and specifications of the network connections supported by the system.

Security

- The port for DICOM communication (specified by the user in the system settings; typically port 104, 2762, or 11112) is used for outgoing communication to the WLAN.
- Anti-virus software is not installed on this device.
- This device has a single configurable listening port for the purposes of DICOM Echo and Storage Commitment.

Data flow

DICOM

MWL Server-----> SII system-----> PACS

Study order (DICOM MWL)	Study data (DICOM Storage)
----------------------------	-------------------------------

Please refer to the SonoSite SII DICOM Conformance Statement (D18493) for details.

Caution

- 1 Connection of equipment to an IT network that includes other systems, could result in previously unidentified risks to patients, operators or third parties. Before connecting the equipment to an uncontrolled IT Network, make sure that all potential risks resulting from such connections, were identified and evaluated, and suitable countermeasures were put in place. IEC 80001-1:2010 provides guidance for addressing these risks.
 - 2 When a setting of the IT network to which this device is connected has been changed, check that the change does not affect this device and take measures if necessary. Changes to the IT network include:
 - Changes in network configuration (IP address, router etc.)
 - Connection of additional items
 - Disconnection of items
 - Update of equipment
 - Upgrade of equipment
- Any changes to the IT network could introduce new risks requiring additional evaluation to be performed as per item 1 above.

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