---

# SonoSite SII

## Doppler and ECG

*User Guide Supplement*

**Manufacturer**

FUJIFILM SonoSite, Inc.  
21919 30th Drive SE  
Bothell, WA 98021 USA  
T: +1-888-482-9449 or  
+1-425-951-1200  
F: 1-425-951-1201

**EC Authorized Representative**

FUJIFILM SonoSite B.V.  
Joop Geesinkweg 140  
1114 AB Amsterdam,  
The Netherlands

**Australia Sponsor**

FUJIFILM SonoSite Australasia Pty Ltd  
114 Old Pittwater Road  
BROOKVALE, NSW, 2100  
Australia

**Caution**

United States federal law restricts this device to sale by or on the order of a physician.

SonoSite, the SonoSite logo, SonoSite SII are trademarks and registered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P27384-02

Publication date: June 2019

Copyright © 2019 FUJIFILM SonoSite, Inc. All rights reserved.

**CE**  
**2797**

# SonoSite SII Doppler and ECG User Guide Supplement

|   |           |
|---|-----------|
| <b>Introduction</b>                                     | <b>1</b>  |
| Document conventions                                    | 2         |
| Getting Help  | 2         |
| <b>Getting Started</b>                                  | <b>3</b>  |
| Preparing the system                                    | 3         |
| System controls   | 4         |
| Intended uses   | 4         |
| <b>System Setup</b>                                     | <b>5</b>  |
| Cardiac Calculations setup                              | 5         |
| Presets setup   | 5         |
| <b>Imaging</b>  | <b>6</b>  |
| 2D imaging  | 6         |
| PW and CW Doppler imaging                               | 6         |
| Imaging modes and exams available by transducer         | 9         |
| ECG   | 15        |
| <b>Measurements and calculations</b>                    | <b>16</b> |
| Doppler measurements                                    | 16        |
| General calculations                                    | 19        |
| Arterial calculations                                   | 20        |
| Cardiac calculations                                    | 21        |
| <b>Measurement references</b>                           | <b>34</b> |
| Measurement accuracy                                    | 34        |
| Measurement publications and terminology                | 35        |
| <b>Cleaning and disinfecting</b>                        | <b>43</b> |
| Cleaning and disinfecting the ECG cable and slave cable | 43        |
| <b>Safety</b>   | <b>43</b> |
| Electrical safety classification                        | 43        |
| Electrical safety                                       | 43        |
| Compatible accessories and peripherals                  | 44        |
| <b>Acoustic output</b>                                  | <b>45</b> |
| Guidelines for reducing TI                              | 45        |
| Output display  | 46        |
| Acoustic output tables                                  | 48        |

## Introduction

This user guide supplement provides information on PW and CW Doppler modes and the ECG option, now available with the SonoSite SII ultrasound system.

## Document conventions

The document follows these conventions:

- ▶ A **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A **Caution** describes precautions necessary to protect the products.
- ▶ A **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ❖.

For a description of labeling symbols that appear on the product, see "Labeling Symbols" in the ultrasound system user guide.

## Getting Help

For technical support, please contact FUJIFILM SonoSite as follows:

|   |  |
|---|--|
| <b>Phone<br/>(U.S. or Canada)</b>             | 877-657-8118   |
| <b>Phone<br/>(outside U.S. or<br/>Canada)</b> | 425-951-1330, or call your local representative  |
| <b>Fax</b>                                    | 425-951-6700   |
| <b>Email</b>                                  | ffss-service@fujifilm.com  |
| <b>Web</b>                                    | www.sonosite.com   |
| <b>Europe Service Center</b>                  | Main: +31 20 751 2020<br>English support: +44 14 6234 1151<br>French support: +33 1 8288 0702<br>German support: +49 69 8088 4030<br>Italian support: +39 02 9475 3655<br>Spanish support: +34 91 123 8451 |
| <b>Asia Service Center</b>                    | +65 6380-5581  |


Printed in the U.S.

# Getting Started

## Preparing the system


### Components and connectors

You can now connect an ECG cable to the back of the system.

## System controls

- |    |                                  |  |
|----|----------------------------------|--|
| 1  | Control knobs                    | Turn to adjust gain, depth, cine buffer, brightness, and more, depending on context. Current functions appear on-screen above the knobs. |
| 2  | Freeze key                       | Press and hold to freeze or unfreeze the image.  |
| 3  | Touchpad                         | When the touchpad is lit, use it to control items displayed on the screen. Double-tap the touchpad to switch between functions.          |
| 4  | Touchpad key                     | Works in conjunction with the touchpad. Tap to activate an item on-screen, or to switch between functions.                               |
| 5  | Print key                        | Available only when a printer is connected to the system. Tap to print from a live or frozen scan.                                       |
| 6  | Save keys                        | Tap one of these keys to save an image or a clip.  |
| 7  | Image mode                       | Tap one of these keys to change the imaging mode.  |
| 8  | System controls                  | Change system settings, switch transducers, add labels, or see patient information.  |
| 9  | Image, ECG, and Doppler controls | Use these to adjust the image, select the ECG function, or select the Doppler imaging mode.  |
| 10 | Touchscreen                      | Use the touchscreen the same way you would use the touchpad.   |

## Intended uses

### Cardiac Imaging Applications

You can use the licensed FUJIFILM SonoSite ECG function to display the patient's heart rate and to provide a cardiac cycle reference when viewing an ultrasound image.

#### **WARNING**

Do not use the SonoSite ECG to diagnose cardiac arrhythmias or to provide long-term cardiac monitoring.

# System Setup

## Cardiac Calculations setup

On the Cardiac Calculations settings page, you can specify measurement names that appear in the Tissue Doppler Imaging (TDI) calculations menu and on the report page. See “**Cardiac calculations**” on page 21.

### To specify cardiac measurement names

❖ Under **TDI Walls** on the Cardiac Calculations settings page, select a name for each wall.

## Presets setup

The Presets setup page has settings for general preferences.

### Doppler Scale

Select **cm/s** or **kHz**.

### Duplex

Specifies the screen layout for displaying the M Mode trace and the Doppler spectral trace:

- ▶ **1/3 2D, 2/3 Trace**
- ▶ **1/2 2D, 1/2 Trace**
- ▶ **Full 2D, Full Trace**

### Live Trace

Select **Peak** or **Mean** velocity trace.

# Imaging

## 2D imaging

**Table 1: 2D controls**

| Control      | Description  |
|--------------|--|
| <b>Guide</b> | <b>Guide</b> is not available when the ECG cable is connected.                                   |
| <b>ECG</b>   | Displays the ECG signal.<br>This feature is optional and requires a FUJIFILM SonoSite ECG cable. |

## PW and CW Doppler imaging

Pulsed wave (PW) Doppler and continuous wave (CW) Doppler imaging modes are optional features. The default Doppler imaging mode is PW Doppler. In cardiac exams, you can select the CW Doppler or TDI Doppler on-screen control.

PW Doppler is a Doppler recording of blood flow velocities in a range specific area (sample volume) along the length of the beam. CW Doppler is a Doppler recording of blood flow velocities along the length of the beam.

### To display the D-line

1 Tap the **Doppler** control at the bottom of the touchscreen.

**Note**

If the D-line does not appear, make sure that the image isn't frozen.

2 Do any of the following as needed:

- ▶ Adjust controls.
- ▶ Drag your finger on the touchscreen or touchpad to position the D-line and gate where desired. Horizontal movements position the D-line. Vertical movements position the gate.
- ▶ To change the gate size, repeatedly press the right knob or tap the on-screen control above the knob until **Gate** appears, and then turn the knob to the gate size you want. To correct the angle, repeatedly press the right knob or tap the on-screen control above the knob until **Angle** appears, and then turn the knob to the correct angle.

**WARNING**

We do not recommend angle correction for the cardiac exam type.

## To display the spectral trace

### Note

Moving the baseline, scrolling, or inverting the trace while the image is frozen will clear displayed cardiac output results.

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The time scale above the trace has small marks at 200 ms intervals and large marks at one-second intervals.

3 Do any of the following as needed:

- ▶ Adjust the sweep speed (**Med**, **Fast**, **Slow**).
- ▶ Tap **Update** to toggle between the D-line and spectral trace.

## Doppler controls

**Table 2: Doppler on-screen controls**

| Control                        | Description   |
|--------------------------------|---|
| <b>PW Dop, CW Dop, TDI Dop</b> | Toggle between PW Doppler, CW Doppler, and TDI Doppler. The current selection appears in the upper left-hand screen. CW Doppler and TDI Doppler are available only in cardiac exams.  |
| <b>Gate</b>                    | Settings depend on transducer and exam type. Use the right knob to adjust the Doppler gate size. The Doppler gate size indicator is on the upper left-hand screen.  |
| <b>Angle</b>                   | Press the right knob to select <b>Angle</b> , and then turn the knob to choose between: <b>0°</b> , <b>+60°</b> , or <b>-60°</b> . We do not recommend angle correction for the cardiac exam type.  |
| <b>Steering</b>                | Select the desired steering angle setting. Settings available depend on the transducer. The PW Doppler angle correction automatically changes to the optimum setting. <ul style="list-style-type: none"> <li>▶ <b>-15</b> and <b>-20</b> have an angle correction of <math>-60^\circ</math>.</li> <li>▶ <b>0</b> has an angle correction of <math>0^\circ</math>.</li> <li>▶ <b>+15</b> and <b>+20</b> have an angle correction of <math>+60^\circ</math>.</li> </ul> You can manually correct the angle after selecting a steering angle setting. Available on select transducers. |

**Table 2: Doppler on-screen controls (continued)**

| Control  | Description  |
|--|--|
| <b>Volume</b><br> | Increases or decreases Doppler speaker volume ( <b>0-10</b> ). |
| <b>Zoom</b>  | Magnifies the image.   |

**Spectral trace controls****Table 3: Spectral trace on-screen controls**

| Control   | Description  |
|---|--|
| <b>Scale</b>  | Press the right knob to select <b>Scale</b> , and then turn the knob to choose the desired velocity setting [pulse repetition frequency (PRF)] in cm/s or kHz.                       |
| <b>Line</b>   | Press the right knob to select <b>Line</b> , and then turn the knob to set the baseline position. (On a frozen trace, the baseline can be adjusted if <b>Trace</b> is off.)          |
| <b>Invert</b>   | Press the right knob to select <b>Invert</b> , and then turn the knob to vertically flip the spectral trace. (On a frozen trace, <b>Invert</b> is available if <b>Trace</b> is off.) |
| <b>Volume</b><br>        | Increases or decreases Doppler speaker volume ( <b>0-10</b> ).   |
| <b>Wall Filter</b><br>   | Settings include <b>Low, Med, High</b> .   |
| <b>Sweep Speed</b><br> | Settings include <b>Slow, Med, Fast</b> .  |
| <b>Trace</b>  | Displays a live trace of the peak or mean. Specify peak or mean on the Presets setup page. Select <b>Above</b> or <b>Below</b> to position the trace above or below the baseline.    |

## Imaging modes and exams available by transducer

**Table 4: Imaging modes and exams available by transducer**

| Transducer | Exam type <sup>a</sup> | Imaging mode              |                  |                    |                            |               |
|------------|------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|            |                        | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Color <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| C8x        | Pro                    | ✓                         | ✓                | ✓                  | ✓                          |               |
| C11x       | Abd                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Art                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Neo                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Ven                    | ✓                         | ✓                | ✓                  | ✓                          |               |
| C35x       | Abd                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Msk                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | OB                     | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Spn                    | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup>Exam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Ven = Venous.

<sup>b</sup>The optimization settings for 2D are Res, Gen, and Pen.

<sup>c</sup>The optimization settings for CPD and Color are low, medium, and high (flow velocity range) with a range of PRF settings for Color depending on the setting selected.

<sup>d</sup>For the cardiac exam type, PW TDI is also available. See **"Doppler controls"** on page 7.

<sup>e</sup>For more information refer to the *P11x Transducer User Guide*, included with the P11x transducer. The P11x transducer is not licensed for use in Canada.

**Table 4: Imaging modes and exams available by transducer (continued)**

| Transducer                  | Exam type <sup>a</sup> | Imaging mode              |                  |                    |                            |               |
|-----------------------------|------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                             |                        | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Color <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| rC60xi standard/<br>armored | Abd                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Gyn                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Msk                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | OB                     | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup>Exam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Ven = Venous.

<sup>b</sup>The optimization settings for 2D are Res, Gen, and Pen.

<sup>c</sup>The optimization settings for CPD and Color are low, medium, and high (flow velocity range) with a range of PRF settings for Color depending on the setting selected.

<sup>d</sup>For the cardiac exam type, PW TDI is also available. See **“Doppler controls”** on page 7.

<sup>e</sup>For more information refer to the *P11x Transducer User Guide*, included with the P11x transducer. The P11x transducer is not licensed for use in Canada.

**Table 4: Imaging modes and exams available by transducer (continued)**

| Transducer                   | Exam type <sup>a</sup> | Imaging mode              |                  |                    |                            |               |
|------------------------------|------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                              |                        | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Color <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| HFL38xi standard/<br>armored | Art                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Bre                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Lung                   | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Msk                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Oph                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | SmP                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Ven                    | ✓                         | ✓                | ✓                  | ✓                          |               |
| HFL50x                       | Bre                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Msk                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                              | SmP                    | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup>Exam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Ven = Venous.

<sup>b</sup>The optimization settings for 2D are Res, Gen, and Pen.

<sup>c</sup>The optimization settings for CPD and Color are low, medium, and high (flow velocity range) with a range of PRF settings for Color depending on the setting selected.

<sup>d</sup>For the cardiac exam type, PW TDI is also available. See **"Doppler controls"** on page 7.

<sup>e</sup>For more information refer to the *P11x Transducer User Guide*, included with the P11x transducer. The P11x transducer is not licensed for use in Canada.

**Table 4: Imaging modes and exams available by transducer (continued)**

| Transducer | Exam type <sup>a</sup> | Imaging mode              |                  |                    |                            |               |
|------------|------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|            |                        | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Color <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| HSL25x     | Art                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Lung                   | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Msk                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Oph                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Sup                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Ven                    | ✓                         | ✓                | ✓                  | ✓                          |               |
| ICTx       | Gyn                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | OB                     | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup>Exam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Ven = Venous.

<sup>b</sup>The optimization settings for 2D are Res, Gen, and Pen.

<sup>c</sup>The optimization settings for CPD and Color are low, medium, and high (flow velocity range) with a range of PRF settings for Color depending on the setting selected.

<sup>d</sup>For the cardiac exam type, PW TDI is also available. See **“Doppler controls”** on page 7.

<sup>e</sup>For more information refer to the *P11x Transducer User Guide*, included with the P11x transducer. The P11x transducer is not licensed for use in Canada.

**Table 4: Imaging modes and exams available by transducer (continued)**

| Transducer                 | Exam type <sup>a</sup> | Imaging mode              |                  |                    |                            |               |
|----------------------------|------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                            |                        | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Color <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| L25x standard/<br>armored  | Art                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Lung                   | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Msk                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Oph                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Sup                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Ven                    | ✓                         | ✓                | ✓                  | ✓                          |               |
| L38xi standard/<br>armored | Art                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Lung                   | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Nrv                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | SmP                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Ven                    | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup>Exam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Ven = Venous.

<sup>b</sup>The optimization settings for 2D are Res, Gen, and Pen.

<sup>c</sup>The optimization settings for CPD and Color are low, medium, and high (flow velocity range) with a range of PRF settings for Color depending on the setting selected.

<sup>d</sup>For the cardiac exam type, PW TDI is also available. See **"Doppler controls"** on page 7.

<sup>e</sup>For more information refer to the *P11x Transducer User Guide*, included with the P11x transducer. The P11x transducer is not licensed for use in Canada.

**Table 4: Imaging modes and exams available by transducer (continued)**

| Transducer                 | Exam type <sup>a</sup> | Imaging mode              |                  |                    |                            |               |
|----------------------------|------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                            |                        | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Color <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| P10x                       | Abd                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Crd                    | ✓                         |                  | ✓                  | ✓                          | ✓             |
|                            | Neo                    | ✓                         | ✓                | ✓                  | ✓                          |               |
| P11x <sup>e</sup>          | Art                    | ✓                         | ✓                | ✓                  |                            |               |
|                            | Ven                    | ✓                         | ✓                | ✓                  |                            |               |
| rP19x standard/<br>armored | Abd                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Crd                    | ✓                         |                  | ✓                  | ✓                          | ✓             |
|                            | Lung                   | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | OB                     | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Orb                    | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | TCD                    | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup>Exam type abbreviations are as follows: Abd = Abdomen, Art = Arterial, Bre = Breast, Crd = Cardiac, Gyn = Gynecology, Msk = Musculoskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrical, Oph = Ophthalmic, Orb = Orbital, SmP = Small Parts, Sup = Superficial, TCD = Transcranial Doppler, Ven = Venous.

<sup>b</sup>The optimization settings for 2D are Res, Gen, and Pen.

<sup>c</sup>The optimization settings for CPD and Color are low, medium, and high (flow velocity range) with a range of PRF settings for Color depending on the setting selected.

<sup>d</sup>For the cardiac exam type, PW TDI is also available. See **“Doppler controls”** on page 7.

<sup>e</sup>For more information refer to the *P11x Transducer User Guide*, included with the P11x transducer. The P11x transducer is not licensed for use in Canada.

## ECG

ECG is an option and requires a FUJIFILM SonoSite ECG cable.

### WARNINGS

- ▶ Do not use the SonoSite ECG to diagnose cardiac arrhythmias or to provide long-term cardiac monitoring.
- ▶ To avoid electrical interference with aircraft systems, do not use the ECG cable on aircraft. Such interference may have safety consequences.

### Caution

- ▶ Use only accessories recommended by FUJIFILM SonoSite with the system. Connecting an accessory not recommended by FUJIFILM SonoSite can damage the system.

### To use the ECG

- 1 Connect the ECG cable to the ECG connector on the back of the ultrasound system. ECG turns on automatically if the system is in live imaging mode.

### Note

The ECG signal may take up to one minute to restabilize after defibrillator use on the patient.


- 2 Tap the **ECG** control at the bottom of the touchscreen.

The ECG controls appear on the screen.

- 3 Adjust controls as desired.

### ECG controls

**Table 5: ECG on-screen controls**

| Control   | Description   |
|---|---|
| <b>Show/Delay/Hide</b>  | Turns on and off ECG signal with and without the Delay line.  |
| <b>ECG Gain</b>   | Tap the ECG gain control  , and then tap the up or down arrows to increase or decrease the ECG Gain from 0-20. |
| <b>Position</b>   | Press the right knob to select <b>Position</b> , and then turn the knob to set the position of the ECG signal.  |
| <b>Sweep Speed</b><br> | Settings are <b>Slow</b> , <b>Med</b> , and <b>Fast</b> .   |

**Table 5: ECG on-screen controls**

| Control   | Description  |
|---|--|
| <b>Delay</b><br> | Tap <b>Delay</b> , then select the position of the delay line on the ECG signal by tapping one of the icons. The delay line indicates where the clip acquisition is triggered. Select <b>Save</b> to save the current position on the ECG signal. (You can change the position of the delay line temporarily. Starting a new patient information form or cycling system power reverts the delay line to the most recently saved position.) |
| <b>Clips</b>  | Tap <b>Clips</b> , then tap <b>Time</b> to change the clips control to <b>ECG</b> . With <b>ECG</b> , you have the option to capture clips based on the number of heart beats. Tap the <b>beats</b> control, then the up or down arrows, to select the number of beats. If <b>Time</b> is selected, capturing is based on number of seconds. Select the time duration.   |

## Measurements and calculations

You can perform basic measurements in any imaging mode and can save the image with the measurements displayed. Except for the M Mode HR measurement, the results do not automatically save to a calculation and the patient report. To save measurements as part of a calculation, you can first begin a calculation and then measure.

### Doppler measurements

The basic measurements that you can perform in Doppler imaging are:

- ▶ Velocity (cm/s)
- ▶ Pressure Gradient
- ▶ Elapsed Time
- ▶ +/x Ratio
- ▶ Resistive Index (RI)
- ▶ Acceleration

You can also trace manually or automatically. For Doppler measurements, the Doppler scale must be set to cm/s on the Presets setup page.

#### To measure Velocity (cm/s) and Pressure Gradient

This measurement involves a single caliper from the baseline.

**1** On a frozen Doppler spectral trace, tap **Calipers**.

A single caliper appears.

- 2 Drag your finger on either the touchpad or the touchscreen to position the caliper to a peak velocity waveform.

### To measure Velocities, Elapsed Time, Ratio, and Resistive Index (RI) or Acceleration

- 1 On a frozen Doppler spectral trace, tap **Calipers**.

A single vertical caliper appears.

- 2 Using the touchpad or the touchscreen, position the caliper to a peak velocity waveform. Tap  to set the position.

A second vertical caliper appears.

- 3 Drag your finger on either the touchpad or the touchscreen to position the second vertical caliper at the end diastole on the waveform, and then tap .

To make a correction, tap **Delete** above the right knob or press the right knob.

Elapsed time between the times indicated by the two calipers is calculated. Measured velocities are given as results, and a generic ratio between the velocities indicated by the two calipers is calculated.

If the absolute value of the earlier velocity is less than that of the later velocity identified by the calipers, Acceleration is calculated; otherwise, in non-cardiac exams, RI is calculated.

### To measure time duration

- 1 On a frozen Doppler spectral trace, tap **Calipers**.
- 2 Navigate to the second page by tapping the arrow.

- 3 Select **Time** .

A vertical caliper appears.

- 4 Using the touchpad or the touchscreen, position the caliper where desired, and then tap .

A second vertical caliper appears.

- 5 Using the touchpad or the touchscreen, position the second caliper where desired.

### To perform manual trace measurements in Doppler

- 1 On a frozen Doppler spectral trace, tap **Calipers**.
- 2 Navigate to the second page by tapping the arrow.

### 3 Tap **Manual** .

A single caliper appears.

### 4 Using the touchpad or the touchscreen, position the caliper at the beginning of the desired waveform, and then tap to activate the trace.

### 5 Using the touchpad or the touchscreen, trace the waveform, and then tap **Set** or .

To make a correction, tap **Undo** or **Delete**.

#### **WARNING**

When using the touchpad to trace a shape, be careful not to touch  until you are finished with the trace. Doing so may complete the trace prematurely, causing an incorrect measurement and delay of care.

### **To perform automatic trace measurements in Doppler**

#### 1 On a frozen Doppler spectral trace, tap **Calipers**.

#### 2 Navigate to the second page by tapping the arrow.

### 3 Tap **Auto** .

A vertical caliper appears.

### 4 Using the touchpad or the touchscreen, position the caliper at the beginning of the desired waveform, and then tap .

A second vertical caliper appears.

### 5 Using the touchpad or the touchscreen, position the caliper at the end of the desired waveform, and then tap **Set**.

To make a correction, tap **Undo** or **Delete**.

### **Automatic trace results**

Depending on the exam type, the results from automatic tracing include the following:

- ▶ Velocity Time Integral (VTI)
- ▶ Cardiac Output (CO)
- ▶ Peak Velocity (Vmax)
- ▶ Peak Systolic Velocity (PSV)
- ▶ Mean Pressure Gradient (PGmean)
- ▶ Time Average Mean (TAM)
- ▶ Mean Velocity on Peak Trace (Vmean)
- ▶ +/x or Systolic/Diastolic (S/D)

- ▶ Pressure Gradient (PGmax)
- ▶ End Diastolic Velocity (EDV)
- ▶ Acceleration Time (AT)
- ▶ Gate Depth
- ▶ Pulsatility Index (PI)
- ▶ Resistive Index (RI)
- ▶ Time Average Peak (TAP)
- ▶ Minimum Diastolic Velocity (MDV)

## General calculations

### Volume flow calculation

The volume flow calculation is available in the following exam types: Abdomen and Arterial.

Both a 2D and a Doppler measurement are required for the volume flow calculation. For the 2D measurement, you can do either of the following:

- ▶ Measure the diameter of the vessel. This approach is more precise. The measurement overrides the gate size.
- ▶ Use the gate size. If you do not measure the diameter of the vessel, the system automatically uses the gate size and "(gate)" appears in the calculation results. Using this option may result in significant error.

The Doppler sample volume should completely insonate the vessel. You can measure either the time average mean (TAM) or time average peak (TAP).

## Arterial calculations

### WARNINGS

- ▶ To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate.
- ▶ To avoid misdiagnosis or harming the patient outcome, start a new patient form before starting a new patient exam and performing calculations. Starting a new patient form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

In the Arterial exam, you can calculate ICA/CCA ratio, volume, volume flow, and percent reduction. The Arterial calculations that you can perform are listed in the following table.

**Table 6: Arterial calculations**

| Calculation list | Measurement name   | Results                        |
|------------------|--|--------------------------------|
| CCA              | <ul style="list-style-type: none"><li>▶ Prox (Proximal)</li><li>▶ Mid (Middle)</li><li>▶ Dist (Distal)</li><li>▶ Bulb</li></ul>  | s (systolic),<br>d (diastolic) |
| ICA              | <ul style="list-style-type: none"><li>▶ Prox (Proximal)</li><li>▶ Mid (Middle)</li><li>▶ Dist (Distal)</li></ul>                 | s (systolic),<br>d (diastolic) |
| ECA              | <ul style="list-style-type: none"><li>▶ Prox (Proximal)</li><li>▶ Mid (Middle)</li><li>▶ Dist (Distal)</li><li>▶ VArty</li></ul> | s (systolic),<br>d (diastolic) |

### WARNINGS

- ▶ Trace only a single heartbeat. The VTI calculation is not valid if measured with more than one heartbeat.
- ▶ Diagnostic conclusions about blood flow based on VTI alone can lead to improper treatment. Accurate blood flow volume calculations require both the vessel area and velocity of blood flow. In addition, accurate blood flow velocity is dependent on a correct Doppler angle of incidence.

## To perform an Arterial calculation

After you perform arterial measurements, values in the ICA/CCA ratios are selectable on the Arterial page of the patient report.

- 1 On a frozen Doppler spectral trace, tap **Calcs**.
- 2 Do the following for each measurement you want to take:
  - a Under **Left** or **Right**, select the measurement name.
  - b Using the touchpad or touchscreen, position the caliper at the peak systolic waveform, and then tap .
 

A second caliper appears.
  - c Using the touchpad, position the second caliper at the end diastole point on the waveform.
- 3 Tap **Save Calc** to save the calculation.
- 4 To save a picture of the finished calculation, tap .
- 5 Tap **Back** to exit the calculation.

## Cardiac calculations

### WARNINGS

- ▶ To avoid incorrect calculations, verify that the patient information, date, and time settings are accurate.
- ▶ To avoid misdiagnosis or harming the patient outcome, start a new patient form before starting a new patient exam and performing calculations. Starting a new patient form clears the previous patient's data. The previous patient's data will be combined with the current patient if the form is not first cleared.

When performing cardiac calculations, the system uses the heart rate (HR) value present in the patient information form. The HR value can be obtained in any four different ways:

- ▶ Manual entry in the patient information form
- ▶ Doppler measurement
- ▶ M-Mode measurement
- ▶ ECG measurement

The ECG heart rate measurement is only used if the other methods are not available. If the ECG measurement is used, and the HR value in the patient information form is empty, the new HR value is automatically inserted in the patient information form.

The following table shows the measurements required to complete different cardiac calculations.

| Calculation list | Measurement name (imaging mode)   | Results  |
|------------------|---|--|
| EF<br>EF         | <ul style="list-style-type: none"> <li>▶ LVDd (2D or M Mode)</li> <li>▶ LVDs (2D or M Mode)</li> </ul>  | EF<br>LVDFS  |
| LV Vol (EF)      | <ul style="list-style-type: none"> <li>▶ A4Cd (2D)</li> <li>▶ A4Cs (2D)</li> <li>▶ A2Cd (2D)</li> <li>▶ A2Cs (2D)</li> </ul>                  | A4C EF<br>A2C EF<br>LV Vol<br>CO <sup>a</sup><br>SV<br>CI <sup>a</sup><br>SI |
| IVC              | <ul style="list-style-type: none"> <li>▶ Max D (2D or M Mode)</li> <li>▶ Min D (2D or M Mode)</li> </ul>                                      | Collapse ratio   |
| LV<br>LVd        | <ul style="list-style-type: none"> <li>▶ RVW (2D)</li> <li>▶ RVD (2D)</li> <li>▶ IVS (2D)</li> <li>▶ LVD (2D)</li> <li>▶ LVPW (2D)</li> </ul> | EF<br>LVDFS<br>CO <sup>a</sup><br>SV<br>LVESV<br>LVEDV<br>IVSFT<br>LVPWFT    |
| LVs              | <ul style="list-style-type: none"> <li>▶ RVW (2D)</li> <li>▶ RVD (2D)</li> <li>▶ IVS (2D)</li> <li>▶ LVD (2D)</li> <li>▶ LVPW (2D)</li> </ul> | CI <sup>a</sup><br>SI<br>LV Mass (M Mode only)                               |
| HR <sup>a</sup>  | HR (M Mode or Doppler)  | HR   |

<sup>a</sup>HR needed for CO and CI. You can enter the HR measurement on the patient form, or obtain it by measuring in M Mode or Doppler.

<sup>b</sup>dP:dT performed at 100 cm/s and 300 cm/s.

<sup>d</sup>Specified on the cardiac patient report.

<sup>e</sup>Need to measure E (MV measurement) to get E/e' ratio.

| Calculation list | Measurement name (imaging mode)   | Results   |
|------------------|---|---|
| CO               | <ul style="list-style-type: none"> <li>▶ LVOT D (2D)</li> <li>▶ HR (Doppler)</li> <li>▶ LVOT VTI (Doppler)</li> </ul> | CO <sup>a</sup><br>SV<br>CI <sup>a</sup><br>SI<br>VTI<br>HR<br>LVOT D |
| Ao/LA            | ▶ Ao (2D or M Mode)   | Ao<br>LA/Ao   |
|                  | ▶ AAO (2D)  | AAo   |
|                  | ▶ LA (2D or M Mode)   | LA<br>LA/Ao   |
|                  | ▶ LVOT D (2D)   | LVOT D<br>LVOT area   |
|                  | ▶ ACS (M Mode)  | ACS   |
|                  | ▶ LVET (M Mode)   | LVET  |

<sup>a</sup>HR needed for CO and CI. You can enter the HR measurement on the patient form, or obtain it by measuring in M Mode or Doppler.

<sup>b</sup>dP:dT performed at 100 cm/s and 300 cm/s.

<sup>d</sup>Specified on the cardiac patient report.

<sup>e</sup>Need to measure E (MV measurement) to get E/e' ratio.

| Calculation list | Measurement name (imaging mode)   | Results                                 |
|------------------|-----------------------------------|---|
| MV               | ▶ EF: Slope (M Mode)              | EF Slope                                |
|                  | ▶ EPSS (M Mode)                   | EPSS                                    |
|                  | ▶ E (Doppler)                     | E<br>E PG<br>A                          |
|                  | ▶ A (Doppler)                     | A PG<br>E:A                             |
|                  | ▶ PHT (Doppler)                   | PHT<br>MVA<br>Decel time                |
|                  | ▶ VTI (Doppler)                   | VTI<br>Vmax<br>PGmax<br>Vmean<br>PGmean |
|                  | ▶ IVRT (Doppler)                  | time                                    |
|                  | ▶ Adur (Doppler)                  | time                                    |
| MV<br>MR         | ▶ dP:dT <sup>b</sup> (CW Doppler) | dP:dT                                   |
| Area             | ▶ MVA (2D)                        | MV Area                                 |
|                  | ▶ AVA (2D)                        | AV Area                                 |
| Atria            | ▶ LA A4C (2D)<br>▶ LA A2C (2D)    | LA Area<br>LA Volume<br>Biplane         |
|                  | ▶ RA (2D)                         | RA Area<br>RA Volume                    |
| LV mass          | ▶ Epi (2D)                        | LV Mass                                 |
|                  | ▶ Endo (2D)                       | Epi Area                                |
|                  | ▶ Apical (2D)                     | Endo Area<br>D Apical                   |

<sup>a</sup>HR needed for CO and CI. You can enter the HR measurement on the patient form, or obtain it by measuring in M Mode or Doppler.

<sup>b</sup>dP:dT performed at 100 cm/s and 300 cm/s.

<sup>d</sup>Specified on the cardiac patient report.

<sup>e</sup>Need to measure E (MV measurement) to get E/e' ratio.

| Calculation list | Measurement name (imaging mode) | Results                                 |
|------------------|---------------------------------|---|
| AV<br>AV         | ▶ Vmax (Doppler)                | Vmax<br>PGmax                           |
|                  | ▶ VTI (Doppler)                 | VTI<br>Vmax<br>PGmax<br>Vmean<br>PGmean |
| LVOT             | ▶ Vmax (Doppler)                | Vmax<br>PGmax                           |
|                  | ▶ VTI (Doppler)                 | VTI<br>Vmax<br>PGmax<br>Vmean<br>PGmean |
| AI               | ▶ PHT (Doppler)                 | AI PHT<br>AI slope                      |
|                  |                                 |   |
| TV               | ▶ RA pressure <sup>d</sup>      | RVSP                                    |
|                  | ▶ TR Vmax (Doppler)             | Vmax<br>PGmax                           |
|                  | ▶ E (Doppler)<br>▶ A (Doppler)  | E<br>E PG<br>A<br>A PG<br>E:A           |
|                  | ▶ PHT (Doppler)                 | PHT<br>TVA<br>Decel time                |

<sup>a</sup>HR needed for CO and CI. You can enter the HR measurement on the patient form, or obtain it by measuring in M Mode or Doppler.

<sup>b</sup>dP:dT performed at 100 cm/s and 300 cm/s.

<sup>d</sup>Specified on the cardiac patient report.

<sup>e</sup>Need to measure E (MV measurement) to get E/e' ratio.

| Calculation list | Measurement name (imaging mode)  | Results   |
|------------------|--|---|
|                  | ▶ VTI (Doppler)  | VTI<br>Vmax<br>PGmax<br>Vmean<br>PGmean                                   |
| PV               | ▶ Vmax (Doppler)   | Vmax<br>PGmax   |
|                  | ▶ PV VTI (Doppler)<br>▶ AT (Doppler)   | VTI<br>Vmax<br>PGmax<br>Vmean<br>PGmean<br>AT                             |
| P Vein           | ▶ A (Doppler)  | Vmax  |
|                  | ▶ Adur (Doppler)   | time  |
|                  | ▶ S (Doppler)<br>▶ D (Doppler)   | Vmax<br>S/D ratio   |
| PISA             | ▶ Radius (Color)<br>▶ MR VTI (Doppler)<br>▶ Ann D (2D)<br>▶ MV VTI (Doppler)   | PISA Area<br>ERO<br>MV Rate<br>Regurgitant Volume<br>Regurgitant Fraction |
| Qp/Qs            | ▶ LVOT D (2D)<br>▶ RVOT D (2D)<br>▶ LVOT VTI (Doppler)<br>▶ RVOT VTI (Doppler) | D<br>VTI<br>Vmax<br>PGmax<br>Vmean<br>PGmean<br>SV<br>Qp/Qs               |

<sup>a</sup>HR needed for CO and CI. You can enter the HR measurement on the patient form, or obtain it by measuring in M Mode or Doppler.

<sup>b</sup>dP:dT performed at 100 cm/s and 300 cm/s.

<sup>d</sup>Specified on the cardiac patient report.

<sup>e</sup>Need to measure E (MV measurement) to get E/e' ratio.

| Calculation list | Measurement name (imaging mode)  | Results                 |
|------------------|--|-------------------------|
| TDI              | <ul style="list-style-type: none"> <li>▶ Sep e' (Doppler)</li> <li>▶ Sep a' (Doppler)</li> <li>▶ Lat e' (Doppler)</li> <li>▶ Lat a' (Doppler)</li> <li>▶ Inf e' (Doppler)</li> <li>▶ Inf a' (Doppler)</li> <li>▶ Ant e' (Doppler)</li> <li>▶ Ant a' (Doppler)</li> </ul> | E/e' ratio <sup>e</sup> |
| TAPSE            | TAPSE (M Mode)   | TAPSE cm                |

<sup>a</sup>HR needed for CO and CI. You can enter the HR measurement on the patient form, or obtain it by measuring in M Mode or Doppler.

<sup>b</sup>dP:dT performed at 100 cm/s and 300 cm/s.

<sup>d</sup>Specified on the cardiac patient report.

<sup>e</sup>Need to measure E (MV measurement) to get E/e' ratio.

## To measure heart rate in Doppler

### Note

Saving the heart rate to the patient report overwrites any heart rate entered on the patient information form.

- 1 On a frozen Doppler spectral trace, tap **Calcs**.
- 2 From the calculations menu, tap **HR**.  
A vertical caliper appears.
- 3 Drag the first vertical caliper to the peak of the heartbeat, and then tap  to set the caliper position.  
A second vertical caliper appears and is active.
- 4 Drag the second vertical caliper to the peak of the next heartbeat.
- 5 Tap **Save Calc** to save the calculation.
- 6 To save a picture of the finished calculation, tap .
- 7 Tap **Back** to exit the calculation.

## To calculate Proximal Isovelocity Surface Area (PISA)

The PISA calculation requires a measurement in 2D, a measurement in Color, and two measurements in Doppler spectral trace. After all measurements are saved, the result appears in the patient report.

### 1 Measure from Ann D:

- a On a frozen 2D image, tap **Calcs**.
- b On the calculations menu, tap **PISA**.
- c On the **PISA** calculations list, tap **Ann D**.
- d Position the calipers by dragging.
- e Tap **Save Calc** to save the calculation.  
A check mark appears next to the saved measurement.

### 2 Measure from Radius:

- a On a frozen Color image, tap **Calcs**.
- b On the calculations menu, tap **Radius**.
- c Position the calipers by dragging.
- d Tap **Save Calc** to save the calculation.  
A check mark appears next to the saved measurement.

### 3 On a frozen Doppler spectral trace, tap **Calcs**.

### 4 On the calculations menu, tap **PISA**.

### 5 Do the following for both **MR VTI** and **MV VTI**:

- a On the **PISA** calculations list, select the measurement you want to make.
- b Use the automatic trace tool to trace the waveform. See ["To perform automatic trace measurements in Doppler"](#) on page 18.
- c Tap **Save Calc** to save the calculation.

### 6 To save a picture of the finished calculation, tap .

### 7 Tap **Back** to exit the calculation.

### 8 To measure peak velocity

For each cardiac measurement, the system saves up to five individual measurements and calculates their average. If you take more than five measurements, the most recent measurement replaces the oldest measurement. If you delete a saved measurement from the patient report, the next measurement taken replaces the deleted one in the patient report. The most recently saved measurement appears at the bottom of the calculations menu.

- 1 On a frozen Doppler spectral trace, tap **Calcs**.
- 2 On the calculations menu, tap **MV**, **TV**, **TDI**, or **P.Vein**.
- 3 Do the following for each measurement you want to take:
  - a Select the measurement name from the calculations menu.
  - b Position the calipers by dragging.
  - c Tap **Save Calc** to save the calculation.

A check mark appears next to the saved measurement.

### To calculate Velocity Time Integral (VTI)

This calculation computes other results in addition to VTI including Vmax, PGmax, Vmean, and PGmean.

- 1 On a frozen Doppler spectral trace, tap **Calcs**.
- 2 On the calculations menu, tap **VTI** under **MV**, **AV**, **TV**, or **PV**.
- 3 Use the automatic trace tool to trace the waveform. See [“To perform automatic trace measurements in Doppler”](#) on page 18.
- 4 Tap **Save Calc** to save the calculation.
- 5 To save a picture of the finished calculation, tap .
- 6 Tap **Back** to exit the calculation.

### To calculate Right Ventricular Systolic Pressure (RVSP)

- 1 On a frozen Doppler spectral trace, tap **Calcs**.
- 2 On the calculations menu, tap **TV** and then **TRmax**.
- 3 Position the caliper by dragging.
- 4 Tap **Save Calc** to save the calculation.

#### Note:

This calculation requires the RA pressure. If RA pressure has not been adjusted, the default value of 5 mmHg is used. Adjust the RA pressure in the Cardiac patient report.

- 5 To save a picture of the finished calculation, tap .
- 6 Tap **Back** to exit the calculation.

### To calculate Pressure Half Time (PHT) in MV, AV, or TV

- 1 On a frozen Doppler spectral trace, tap **Calcs**.
- 2 On the calculations menu, tap **MV**, **AV**, or **TV**, and then **PHT**.

Position the first caliper at the peak, and then tap . A second caliper appears.

- 3 Position the second caliper:

- ▶ In MV, position the caliper along the EF slope.
- ▶ In AV, position the caliper at the end diastole.

4 Tap **Save Calc** to save the calculation.

5 To save a picture of the finished calculation, tap .

6 Tap **Back** to exit the calculation.

### To calculate Isovolumic Relaxation Time (IVRT)

1 On a frozen Doppler spectral trace, tap **Calcs**.

On the calculations menu, tap **MV**, and then **IVRT**. A vertical caliper appears.

2 Position the caliper at the aortic valve closure.

3 Tap . A second vertical caliper appears.

4 Position the second caliper at onset of mitral inflow.

5 Tap **Save Calc** to save the calculation.

6 To save a picture of the finished calculation, tap .

7 Tap **Back** to exit the calculation.

### To calculate Delta Pressure: Delta Time (dP:dT)

To perform the dP:dT measurements, the CW Doppler scale must include velocities of 300 cm/s or greater on the negative side of the baseline.

1 On a frozen Doppler spectral trace, tap **Calcs**.

2 On the calculations menu, tap **MV**, and then **dP:dT**.

A horizontal dotted line with an active caliper appears at 100 cm/s.

3 Position the first caliper along the waveform at 100 cm/s.

4 Tap .

A second horizontal dotted line with an active caliper appears at 300 cm/s.

5 Position the second caliper along the waveform at 300 cm/s. Tap **Save Calc** to save the calculation.

6 To save a picture of the finished calculation, tap .

7 Tap **Back** to exit the calculation.

### To calculate Aortic Valve Area (AVA)

The AVA calculation requires a measurement in 2D and two measurements in Doppler. After the measurements are saved, the result appears in the patient report.

**1** In 2D:

- a** On a frozen 2D image, tap **Calcs**.
- b** On the calculations menu, tap **Ao/LA**.
- c** From the **Ao/LA** calculation list, select **LVOT D**.
- d** Position the calipers.
- e** Tap **Save Calc** to save the calculation.

**2** In PW Doppler, measure either LVOT Vmax or LVOT VTI.

- ▶ **Vmax** - Tap **AV**, then tap the **Vmax** measurement under **LVOT**. Position the caliper, and then save the measurement.
- ▶ **VTI** - Tap **AV**, then tap the **VTI** measurement under **LVOT**. Use the automatic trace tool to trace the waveform, and then save the measurement.

**Notes**

If **VTI** is chosen, the Vmax value derived from the trace is used as input to the AVA calculation.

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**Notes**

- ▶ If **VTI** is chosen, the Vmax value derived from the trace is used as input to the AVA calculation.
- ▶ If VTI measurements are made for both LVOT and AV, a second AVA result is provided.

**To calculate Qp/Qs**

The Qp/Qs calculation requires two measurements in 2D and two measurements in Doppler. After the measurements are saved, the result appears in the patient report.

- 1** On a frozen 2D image, tap **Calcs**.
- 2** Do the following to measure from LVOT D and again to measure from RVOT D:
  - a** From the **Qp/Qs** calculations list, select **LVOT D** or **RVOT D**.
  - b** Position the calipers.
  - c** Tap **Save Calc** to save the calculation.
- 3** On a frozen Doppler spectral trace, tap **Calcs**.
- 4** Do the following to measure from LVOT VTI and again to measure from RVOT VTI:
  - a** On the calculations menu, tap **Qp/Qs** and then **LVOT VTI** or **RVOT VTI**.

- b** Use the automatic trace tool to trace the waveform. See [“To perform automatic trace measurements in Doppler”](#) on page 18.
- c** Tap **Save Calc** to save the calculation.

### **To calculate Stroke Volume (SV) or Stroke Index (SI)**

The SV and SI calculations require a measurement in 2D and a measurement in Doppler. SI also requires Body Surface Area (BSA). After the measurements are saved, the result appears in the patient report.

- 1** (SI Only) Fill in the **Height** and **Weight** fields on the patient form. The BSA is calculated automatically.
- 2** Measure from LVOT (2D):
  - a** On a frozen 2D image, tap **Calcs**.
  - b** On the calculations menu, tap **Ao/LA** then **LVOT D**.
  - c** Position the calipers.
  - d** Tap **Save Calc** to save the calculation.
- 3** Measure from **LVOT** (Doppler). Refer to [“To calculate Velocity Time Integral \(VTI\)”](#) on page 29. On the calculations menu, tap **AV** and then **LVOT VTI**.

### **To calculate Cardiac Output (CO) or Cardiac Index (CI)**

The CO and CI calculations require Stroke Volume (SV) and Heart Rate (HR) calculations. CI also requires Body Surface Area (BSA). After the measurements are saved, the result appears in the patient report.

- 1** (CI Only) Fill in the **Height** and **Weight** fields on the patient form. The BSA is calculated automatically.
- 2** Calculate SV as described in [“To calculate Stroke Volume \(SV\) or Stroke Index \(SI\)”](#) on page 32.
- 3** Calculate HR as described in [“To measure heart rate in Doppler”](#) on page 27.

## To calculate Cardiac Output (CO) automatically

Make sure that the flow rate is 1 L/min or greater. The system can maintain accuracy of the measurements only if the flow rate is 1 L/min or greater.

### WARNINGS

- ▶ To avoid incorrect calculation results, make sure that the Doppler signal does not alias.
- ▶ To avoid an incorrect diagnosis:
  - ▶ Do not use automatic Cardiac Output calculations as the sole diagnostic criteria. Use them only in conjunction with other clinical information and patient history.
  - ▶ Do not use automatic Cardiac Output calculations in neonatal or pediatric patients.
  - ▶ To avoid inaccurate velocity measurements if you use PW Doppler, make sure that the angle is set to zero

### 1 Measure from LVOT:

- a On a frozen 2D image, tap **Calcs**.
- b On the **CO** calculations menu, tap **LVOT D**.
- c Position the calipers by dragging.
- d Tap **Save Calc** to save the calculation.

### 2 Trace automatically in Doppler. The automatic trace tool always measures the peak regardless of the **Live Trace** setting in Presets setup.

- a Display the live Doppler spectral trace.
- b Tap the arrow to navigate to the next page.
- c Tap **Trace**, and then select **Above** or **Below** to position the automatic trace tool relative to the baseline.
- d Freeze the image, then tap **Calipers**.

### e Tap **Auto** .

A vertical caliper appears.

### f Using the touchpad or the touchscreen, position the caliper at the beginning of the desired waveform,

and then tap .

A second vertical caliper appears

- g Using the touchpad or the touchscreen, position the caliper at the end of the desired waveform, and then tap **Set**.

**Note**

If you invert the frozen image or move the baseline, results are cleared.

- h Tap **Save Calc** to save the calculation.

**To measure a Tissue Doppler Imaging (TDI) waveform**

- 1 Ensure that TDI is on.
- 2 On a frozen Doppler spectral trace, tap **Calcs**.
- 3 Tap **TDI** on the calculations menu, and then do the following for each measurement you want to take:
  - a On the calculations menu, select the measurement name.
  - b Position the calipers.
  - c Tap **Save Calc** to save the calculation.

## Measurement references

### Measurement accuracy

**Table 7: PW Doppler mode measurement and calculation accuracy and range**

| Doppler mode measurement accuracy and range | System tolerance                            | Accuracy by | Test method <sup>a</sup> | Range                  |
|---|---|-------------|--------------------------|------------------------|
| Velocity cursor                             | < +/- 2% plus 1% of full scale <sup>b</sup> | Acquisition | Phantom                  | 0.01 cm/sec-550 cm/sec |
| Frequency cursor                            | < +/- 2% plus 1% of full scale <sup>b</sup> | Acquisition | Phantom                  | 0.01kHz-20.8 kHz       |
| Time  | < +/- 2% plus 1% of full scale <sup>c</sup> | Acquisition | Phantom                  | 0.01-10 sec            |

<sup>a</sup>FUJIFILM SonoSite special test equipment was used.

<sup>b</sup>Full scale for frequency or velocity implies the total frequency or velocity magnitude, displayed on the scrolling graphic image.

<sup>c</sup>Full scale for time implies the total time displayed on the scrolling graphic image.

## Measurement publications and terminology

### Cardiac references

#### Acceleration (ACC) in $\text{cm/s}^2$

Zwiebel, W.J. *Introduction to Vascular Ultrasonography*. 4th ed., W.B. Saunders Company, (2000), p.52.

ACC = abs (delta velocity/delta time)

#### Acceleration Time (AT) in msec

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Lippincott, Williams, and Wilkins, (2007), p.147-148.

[time a - time b]

where: time a = early time;  
time b = later time;

only valid when [a] > [b]

#### Aortic Valve Area (AVA) by Continuity Equation in $\text{cm}^2$

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Lippincott, Williams, and Wilkins, (2007), p.73 and p.191-195.

$$A_2 = A_1 * V_1/V_2$$

where:  $A_2$  =  $A_o$  valve area

$A_1$  = LVOT area;

$V_1$  = Peak LVOT velocity (Vmax) or LVOT VTI

$V_2$  = Peak  $A_o$  valve velocity (Vmax) or  $A_o$  VTI

LVOT = Left Ventricular Outflow Tract

#### Deceleration Time in msec

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Lippincott, Williams, and Wilkins, (2007), p.73-74.

[time a - time b]

where: time a = time associated with Vmax;  
time b = when the line tangent to the envelope and through Vmax crosses the baseline

### **Delta Pressure: Delta Time (dP:dT) in mmHg/s**

Otto, C.M. *Textbook of Clinical Echocardiography*. 2nd ed., W.B. Saunders Company, (2000), p.117-118.

32 mmHg/time interval in seconds

### **E:A Ratio in cm/sec**

E:A = velocity E/velocity A

### **E/Ea Ratio**

Reynolds, Terry. *The Echocardiographer's Pocket Reference*. 2nd ed., School of Cardiac Ultrasound, Arizona Heart Institute, (2000), p.225.

E Velocity/Ea velocity

where: E velocity = Mitral Valve E velocity

Ea = annular E velocity, also known as E prime

### **Effective Regurgitant Orifice (ERO) in mm<sup>2</sup>**

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Lippincott, Williams, and Wilkins, (2007), p.73-76, p.210.

$ERO = MV \text{ Flow Rate} / MR \text{ Vel} * 100$

### **Elapsed Time (ET) in msec**

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Lippincott, Williams, and Wilkins, (2007), p.147, Figure 9-8.

ET = time between velocity cursors in milliseconds

### **Isovolumic Relaxation Time (IVRT) in msec**

Reynolds, Terry, *The Echocardiographer's Pocket Reference*, 2nd Edition, School of Cardiac Ultrasound, Arizona Heart Institute, (2000), p.385.

[time a - time b]

where: time a = mitral valve opening

time b = aortic valve closure

## IVC Percentage Collapse

Lyon, M., N. Verma. "Ultrasound guided volume assessment using inferior vena cava diameter." *The Open Emergency Medicine Journal*. 2010, 3: p.22-24.

$$(IVCd_{exp} - IVCd_{insp}) / IVCd_{exp} \times 100$$

where: expiration (exp) = maximum diameter (Max D)

inspiration (insp) = minimum diameter (Min D)

## LV Ejection Fraction

Schiller, N.B., Shah, P.M., Crawford, M., et al. "Recommendations for Quantification of the Left Ventricle by Two-Dimensional Echocardiography." *Journal of American Society of Echocardiography*. September-October 1989, 2: p.364.

$$EF = ((\text{End Diastolic Volume} - \text{End Systolic Volume}) / \text{End Diastolic Volume}) * 100 (\%)$$

## Mean Velocity (Vmean) in cm/s

Vmean = mean velocity

## Mitral Valve Area (MVA) in cm<sup>2</sup>

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.73-74.

$$MVA = 220 / PHT$$

where: PHT = pressure half time

220 is an empirical derived constant and may not accurately predict mitral valve area in mitral prosthetic heart valves. The mitral valve area continuity equation may be utilized in mitral prosthetic heart valves to predict effective orifice area.

## MV Flow Rate in cc/sec

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.73-76, p.210.

$$\text{Flow} = \text{PISA} * V_a$$

where: PISA = Proximal Isovelocity Surface Area

V<sub>a</sub> = aliasing Velocity

## Pressure Gradient (PGr) in mmHG

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 2nd ed., Philadelphia: Lippincott, Williams, and Wilkins, (1999), p.64.

$$\text{PGr} = 4 * (\text{Velocity})^2$$

Peak E Pressure Gradient (E PG)

$$\text{E PG} = 4 * \text{PE}^2$$

Peak A Pressure Gradient (A PG)

$$\text{A PG} = 4 * \text{PA}^2$$

Peak Pressure Gradient (PGmax)

$$\text{PGmax} = 4 * \text{VMax}^2$$

Mean Pressure Gradient (PGmean)

PGmean = Average pressure gradient during the flow period

Baumgartner, H., Hung, J., Bermejo, J., et al. "Echocardiographic Assessment of Valve Stenosis: EAE/ASE Recommendations for Clinical Practice". *Journal of American Society of Echocardiography*. January 2009, p. 4-5.

$$\text{PG mean} = \text{sum}(4v^2)/N$$

where: v = peak velocity at interval n

N = the number of intervals in the Riemann sum

## Pressure Half Time (PHT) in msec

Reynolds, Terry, *The Echocardiographer's Pocket Reference*, 2nd Edition, School of Cardiac Ultrasound, Arizona Heart Institute, (2000), p. 391.

$$\text{PHT} = \text{DT} * 0.29 \text{ (time required for the pressure gradient to fall half its maximum level)}$$

where: DT = deceleration time

## Proximal Isovelocity Surface Area (PISA) in cm<sup>2</sup>

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.74-76.

$$\text{PISA} = 2 \pi r^2$$

where: r = aliasing radius

**Qp/Qs**

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.70-72.

$$Qp/Qs = SV_{Qp \text{ site}}/SV_{Qs \text{ site}} = RVOT \text{ SV}/LVOT \text{ SV}$$

where:  $RVOT \text{ SV} = RVOT \text{ CSA} * RVOT \text{ VTI} = \pi/4 * RVOT \text{ diameter}^2 * RVOT \text{ VTI}$

$$LVOT \text{ SV} = LVOT \text{ CSA} * LVOT \text{ VTI} = \pi/4 * LVOT \text{ diameter}^2 * LVOT \text{ VTI}$$

**Regurgitant Fraction (RF) in percent**

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.215-217.

$$RF = RV / MV \text{ SV}$$

where:  $RV = \text{Regurgitant Volume}$

$MV \text{ SV} = \text{Mitral Stroke Volume (Mitral CSA} * \text{Mitral VTI)}$

Mitral CSA = cross-sectional area calculated using annulus diameter

**Regurgitant Volume (RV) in cc**

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Lippincott, Williams, and Wilkins, (2007), p.215-217.

$$RV = ERO * MR \text{ VTI}/100$$

**Right Atrial Volume**

Lang, R., M. Bierig, R. Devereux, et al. "Recommendations for chamber quantification: a report from the American Society of Echocardiography's guidelines and standards committee and the chamber quantification writing group, Developed in conjunction with the European Association of Echocardiography, a branch of the European Society of Cardiology." *Journal of the American Society of Echocardiograph*. 2005, 18: p.1440-1463.

$$RA \text{ Vol} = \pi/4 * \sum(a_i) * a_i * L/20 \text{ for } i = 1 \text{ to } 20 \text{ (number of segments)}$$

where:  $RA \text{ Vol} = \text{Right Atrial Volume in ml}$

$a_i = \text{diameter of chamber view slice } i$

$L = \text{length of the chamber view}$

## Right Atrial Volume Index

Wang, Y., J. Gutman, et al. "Atrial volume in a normal adult population by two-dimensional echocardiography." *Chest*. (1984), 86: p.595-601.

$$\text{RA Vol Index} = \text{RA Vol/BSA (ml/L}^2\text{)}$$

## Right Ventricular Systolic Pressure (RVSP) in mmHg

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.66.

$$\text{RVSP} = 4 * (\text{VMax TR})^2 + \text{RAP}$$

where: RAP = Right Atrial Pressure

## S/D

Reynolds, Terry. *The Echocardiographer's Pocket Reference*. 2nd ed., School of Cardiac Ultrasound, Arizona Heart Institute, (2000), p.217.

S velocity/D velocity

where: S velocity = Pulmonary vein S wave

D velocity= Pulmonary vein D wave

## Stroke Volume (SV) Doppler in ml

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.69-71.

$$\text{SV} = (\text{CSA} * \text{VTI})$$

where: CSA = Cross Sectional Area of the orifice (LVOT area)

VTI = Velocity Time Integral of the orifice (LVOT VTI)

## TAPSE

Rudski, L., W. Lai, et al. "Guidelines for the echocardiographic assessment of the right heart in adults: a report from the American Society of Echocardiography." *Journal of the American Society of Echocardiography*. (2010), p.685-713.

M Mode distance measurement of systolic excursion of the right ventricle

## Tricuspid Valve Area (TVA)

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 3rd ed., Philadelphia: Lippincott, Williams, and Wilkins, (2007), p.73-74.

$$TVA = 220 / PHT$$

## Velocity Time Integral (VTI) in cm

Reynolds, Terry. *The Echocardiographer's Pocket Reference*. 2nd ed., School of Cardiac Ultrasound, Arizona Heart Institute, (2000), p.383.

$$VTI = \text{sum of abs (velocities [n])}$$

where: Auto Trace – distance (cm) blood travels with each ejection period. Velocities are absolute values.

## General references

### +/x or S/D Ratio

$$+/x = (\text{Velocity A}/\text{Velocity B})$$

where: A = velocity cursor +  
B = velocity cursor x

## Acceleration Index (AI)

Zwibel, W.J. *Introduction to Vascular Ultrasonography*, 4th ed., W.B. Saunders Company, (2000), p.52.

$$ACC = \text{abs (delta velocity/delta time)}$$

## Elapsed Time (ET)

ET = time between velocity cursors in milliseconds

## Pressure Gradient (PGr) in mmHG

Oh, J.K., J.B. Seward, A.J. Tajik. *The Echo Manual*. 2nd ed., Lippincott, Williams, and Wilkins, (1999), p.64.

$$PG = 4 * (\text{Velocity})^2 \text{ (velocity units must be meters/second)}$$

Peak E Pressure Gradient (E PG)

$$E PG = 4 * PE^2$$

Peak A Pressure Gradient (A PG)

$$A PG = 4 * PA^2$$

Peak Pressure Gradient (PGmax)

$$PG_{max} = 4 * V_{max}^2$$

Mean Pressure Gradient (PGmean)

$$PG_{mean} = 4 * V_{max}^2 \text{ (average pressure gradient during the flow period)}$$

### **Pulsatility Index (PI)**

Zwiebel, W. J., *Introduction to Vascular Ultrasonography*, 4th Edition, W.B. Saunders Company, (2000).

$$PI = (PSV - MDV)/V \text{ (no units)}$$

where: PSV = peak systolic velocity

EDV = minimum diastolic velocity

V = TAP (Time Averaged Peak) flow velocity throughout the cardiac cycle

### **Resistive Index (RI)**

Kurtz, A.B., W.D. Middleton. *Ultrasound-the Requisites*. Mosby Year Book, Inc., (1996), p.467.

$$RI = ((\text{Velocity A} - \text{Velocity B})/\text{Velocity A}) \text{ in measurements}$$

where: A = velocity cursor +

B = velocity cursor x

### **Time Averaged Mean (TAM) in cm/s**

TAM = mean (mean Trace)

### **Time Averaged Peak (TAP) in cm/s**

TAP = mean (peak Trace)

### **Volume Flow (VF) in ml/m**

Robert J. Daigle, BA, RVT: *Techniques in Noninvasive Vascular Diagnosis; Second Edition*, p.210

One of the following, depending on the Live Trace setting:

$$VF = CSA * TAM * 60$$

$$VF = CSA * TAP * 60$$

$$VF = CSA * TAV * 60 \text{ (When manual trace is used)}$$

# Cleaning and disinfecting

## Cleaning and disinfecting the ECG cable and slave cable

**Caution** | To avoid damaging the ECG cable, do not sterilize.

### To clean and disinfect the ECG cables (wipe method)

- 1 Remove the cable from the system.
- 2 Examine the ECG cable for damage such as cracks or splitting.
- 3 Clean the surface using a soft cloth lightly dampened in a mild soap, cleaning solution, or pre-moistened wipe. Apply the solution to the cloth rather than the surface.
- 4 Wipe the surfaces with a FUJIFILM SonoSite Approved cleaner or disinfectant. Refer to the cleaners and disinfection tool available at [www.sonosite.com/support/cleaners-disinfectants](http://www.sonosite.com/support/cleaners-disinfectants).
- 5 Air dry or dry with a clean cloth.

For more information on the ECG slave cable, refer to the *ECG Slave Cable User Guide*.

## Safety

### Electrical safety classification

**Type CF applied parts**                      ECG module/ECG leads

### Electrical safety

**WARNING** | To avoid the risk of electrical shock:

- ▶ Do not allow any part of the system (including the bar code scanner, external mouse, power supply, power supply connector, external keyboard, and so on), except for the transducer or ECG leads, to touch the patient.

## Compatible accessories and peripherals

**Table 8: Accessories and peripherals**

| Description     | Maximum cable length |
|-----------------|----------------------|
| ECG lead wires  | 24 in/0.6 m          |
| ECG module      | 5.8 ft/1.8 m         |
| ECG slave cable | 8 ft/2.4 m           |

# Acoustic output

## Guidelines for reducing TI

**Table 9: Guidelines for reducing TI**

| Transducer   | CPD settings |            |           |     |       |             | PW settings                          |
|--|--------------|------------|-----------|-----|-------|-------------|--------------------------------------|
|  | Box width    | Box height | Box depth | PRF | Depth | Optimize    |                                      |
| C8x  | ↓            |            |           |     | ↑     |             | ↓ (Depth)                            |
| C11x   |              |            | ↑         | ↓   | ↑     |             | ↓ (Depth)                            |
| C35x   | ↑            |            |           | ↓   | ↑     |             | ↓ (Depth)                            |
| rC60xi standard/<br>armored  | ↓            |            |           | ↓   | ↑     |             | ↓ (PRF)                              |
| HFL38xi standard/<br>armored   |              |            | ↑         | ↑   | ↑     |             | ↓ (Depth)                            |
| HFL50x   |              |            | ↑         | ↑   | ↑     |             | ↓ (Depth)                            |
| HSL25x   | ↓            |            |           |     | ↑     |             | ↓ (PRF)                              |
| ICTx   |              | ↑          | ↑         | ↓   |       | Exam<br>Gyn | ↓ (PRF)                              |
| L25x standard/<br>armored  | ↓            |            |           |     | ↑     |             | ↓ (PRF)                              |
| L38xi standard/<br>armored   | ↑            | ↑          |           |     |       |             | ↓ (Sample<br>volume zone<br>or size) |
| P10x   |              |            | ↑         | ↓   |       |             | ↓ (PRF)                              |
| rP19x standard/<br>armored   |              |            |           | ↓   | ↑     |             | ↓ (Depth)                            |
| ↓ Decrease or lower setting of parameter to reduce MI.<br>↑ Increase or raise setting of parameter to reduce MI. |              |            |           |     |       |             |                                      |

## Output display

**Table 10: TI or MI  $\geq$  1.0**

| Transducer                   | Index            | 2D/M Mode | CPD/Color | PW Doppler | CW Doppler |
|------------------------------|------------------|-----------|-----------|------------|------------|
| C8x                          | MI               | Yes       | Yes       | Yes        | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |
| C11x                         | MI               | No        | No        | No         | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |
| C35x                         | MI               | Yes       | No        | No         | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |
| rC60xi standard/<br>armored  | MI               | Yes       | Yes       | Yes        | —          |
|                              | TIC, TIB, or TIS | Yes       | Yes       | Yes        | —          |
| HFL38xi standard/<br>armored | MI               | Yes       | Yes       | Yes        | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |
| HFL50x                       | MI               | Yes       | Yes       | Yes        | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |
| HSL25x                       | MI               | Yes       | Yes       | No         | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |
| ICTx                         | MI               | No        | No        | No         | —          |
|                              | TIC, TIB, or TIS | No        | No        | Yes        | —          |

Even when MI is less than 1.0, the system provides a continuous real-time display of MI in all imaging modes, in increments of 0.1.

The system meets the output display standard for TI and provides a continuous real-time display of TI in all imaging modes, in increments of 0.1.

The TI consists of three user-selectable indices, and only one of these is displayed at any one time. In order to display TI properly and meet the ALARA principle, the user selects an appropriate TI based on the specific exam being performed. FUJIFILM SonoSite provides a copy of *AIUM Medical Ultrasound Safety*, which contains guidance on determining which TI is appropriate.

**Table 10: TI or MI  $\geq$  1.0 (continued)**

| Transducer                 | Index               | 2D/M Mode | CPD/Color | PW Doppler | CW Doppler |
|----------------------------|---------------------|-----------|-----------|------------|------------|
| L25x standard/<br>armored  | MI                  | Yes       | Yes       | No         | —          |
|                            | TIC, TIB, or<br>TIS | No        | No        | Yes        | —          |
| L38xi standard/<br>armored | MI                  | Yes       | Yes       | Yes        | —          |
|                            | TIC, TIB, or<br>TIS | Yes       | Yes       | Yes        | —          |
| P10x                       | MI                  | No        | No        | Yes        | No         |
|                            | TIC, TIB, or<br>TIS | No        | Yes       | Yes        | Yes        |
| rP19x standard/<br>armored | MI                  | Yes       | Yes       | Yes        | No         |
|                            | TIC, TIB, or<br>TIS | Yes       | Yes       | Yes        | Yes        |

Even when MI is less than 1.0, the system provides a continuous real-time display of MI in all imaging modes, in increments of 0.1.

The system meets the output display standard for TI and provides a continuous real-time display of TI in all imaging modes, in increments of 0.1.

The TI consists of three user-selectable indices, and only one of these is displayed at any one time. In order to display TI properly and meet the ALARA principle, the user selects an appropriate TI based on the specific exam being performed. FUJIFILM SonoSite provides a copy of *AIUM Medical Ultrasound Safety*, which contains guidance on determining which TI is appropriate.

## Acoustic output tables

|   |    |
|---|----|
| Transducer model: C8x Operating mode: PW Doppler .....                    | 49 |
| Transducer model: C11x Operating mode: PW Doppler .....                   | 50 |
| Transducer model: C35x Operating mode: PW Doppler .....                   | 51 |
| Transducer model: rC60xi Operating mode: PW Doppler .....                 | 52 |
| Transducer model: HFL38xi Operating mode: PW Doppler .....                | 53 |
| Transducer Model: HFL38xi Ophthalmic Use Operating Mode: PW Doppler ..... | 54 |
| Transducer model: HFL50x Operating mode: PW Doppler .....                 | 55 |
| Transducer model: HSL25x Operating mode: PW Doppler .....                 | 56 |
| Transducer model: HSL25x Ophthalmic Use Operating mode: PW Doppler .....  | 57 |
| Transducer model: ICTx Operating mode: PW Doppler .....                   | 58 |
| Transducer model: L25x Operating mode: PW Doppler .....                   | 59 |
| Transducer model: L25x Ophthalmic Use Operating mode: PW Doppler .....    | 60 |
| Transducer model: L38xi Operating mode: PW Doppler .....                  | 61 |
| Transducer model: P10x Operating mode: PW Doppler .....                   | 62 |
| Transducer Model: P10x Operating Mode: CW Doppler .....                   | 63 |
| Transducer model: rP19x Operating mode: PW Doppler .....                  | 64 |
| Transducer model: rP19x Orbital Use Operating mode: PW Doppler .....      | 65 |
| Transducer model: rP19x Operating mode: CW Doppler .....                  | 66 |

**Table 11: Transducer model: C8x Operating mode: PW Doppler**

| Index label                   |   |              | M.I.   | TIS  |                |                | TIB      | TIC |
|-------------------------------|---|--------------|--------|------|----------------|----------------|----------|-----|
|                               |   |              |        | Scan | Non-scan       |                | Non-scan |     |
|                               |   |              |        |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |          |     |
| Global maximum index value    |   |              | 1.2    | —    | (a)            | —              | 2.0      | (b) |
| Associated acoustic parameter | $P_{r,3}$                               | (MPa)        | 2.59   |      |                |                |          |     |
|                               | $W_0$                                   | (mW)         |        | —    | #              |                | 36.0     | #   |
|                               | min of [ $W_{,3}(z_1), I_{TA,3}(z_1)$ ] | (mW)         |        |      |                | —              |          |     |
|                               | $z_1$                                   | (cm)         |        |      |                | —              |          |     |
|                               | $Z_{bp}$                                | (cm)         |        |      |                | —              |          |     |
|                               | $Z_{sp}$                                | (cm)         | 1.1    |      |                |                | 1.10     |     |
|                               | $d_{eq}(Z_{sp})$                        | (cm)         |        |      |                |                | 0.28     |     |
|                               | $F_c$                                   | (MHz)        | 4.79   | —    | #              | —              | 4.79     | #   |
|                               | Dim of $A_{aprt}$                       | X (cm)       |        | —    | #              | —              | 1.12     | #   |
| Y (cm)                        |   |              | —      | #    | —              | 0.40           | #        |     |
| Other information             | PD                                      | ( $\mu$ sec) | 1.131  |      |                |                |          |     |
|                               | PRF                                     | (Hz)         | 1008   |      |                |                |          |     |
|                               | $p_r@PII_{max}$                         | (MPa)        | 3.10   |      |                |                |          |     |
|                               | $d_{eq}@PII_{max}$                      | (cm)         |        |      |                |                | 0.28     |     |
|                               | Focal Length                            | $FL_x$ (cm)  |        | —    | #              | —              |          | #   |
|                               |   | $FL_y$ (cm)  |        | —    | #              | —              |          | #   |
| $I_{PA,3}@MI_{max}$           | (W/cm <sup>2</sup> )                    | 296          |        |      |                |                |          |     |
| Operating control conditions  | Control 1: Exam type                    |              | Pro    |      |                |                | Pro      |     |
|                               | Control 2: Sample volume size           |              | 1 mm   |      |                |                | 1 mm     |     |
|                               | Control 3: Sample volume position       |              | Zone 5 |      |                |                | Zone 5   |     |
|                               | Control 4: PRF                          |              | 1008   |      |                |                | 3125     |     |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 12: Transducer model: C11x Operating mode: PW Doppler**

| Index label                   |                                       |              | M.I. | TIS  |                |                | TIB      | TIC    |
|-------------------------------|---------------------------------------|--------------|------|------|----------------|----------------|----------|--------|
|                               |                                       |              |      | Scan | Non-scan       |                | Non-scan |        |
|                               |                                       |              |      |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |          |        |
| Global maximum index value    |                                       |              | (a)  | —    | (a)            | —              | 1.5      | 1.1    |
| Associated acoustic parameter | $P_{r.3}$                             | (MPa)        | #    |      |                |                |          |        |
|                               | $W_0$                                 | (mW)         |      | —    | #              |                | 24.6     | 21.7   |
|                               | min of $[W_{.3}(z_1), I_{TA.3}(z_1)]$ | (mW)         |      |      |                | —              |          |        |
|                               | $z_1$                                 | (cm)         |      |      |                | —              |          |        |
|                               | $Z_{bp}$                              | (cm)         |      |      |                | —              |          |        |
|                               | $Z_{sp}$                              | (cm)         |      |      |                |                | 1.70     |        |
|                               | $z@PII_{.3max}$                       | (cm)         | #    |      |                |                |          |        |
|                               | $d_{eq}(Z_{sp})$                      | (cm)         |      |      |                |                | 0.23     |        |
|                               | $F_c$                                 | (MHz)        | #    | —    | #              | —              | 4.37     | 4.36   |
|                               | Dim of $A_{aprt}$                     | X (cm)       |      | —    | #              | —              | 0.64     | 0.40   |
|                               | Y (cm)                                |              | —    | #    | —              | 0.50           | 0.50     |        |
| Other information             | PD                                    | ( $\mu$ sec) | #    |      |                |                |          |        |
|                               | PRF                                   | (Hz)         | #    |      |                |                |          |        |
|                               | $p_r@PII_{max}$                       | (MPa)        | #    |      |                |                |          |        |
|                               | $d_{eq}@PII_{max}$                    | (cm)         |      |      |                |                | 0.22     |        |
|                               | Focal Length                          | $FL_x$ (cm)  |      | —    | #              | —              |          | 1.52   |
|                               |                                       | $FL_y$ (cm)  |      | —    | #              | —              |          | 4.40   |
| $I_{PA.3}@MI_{max}$           | (W/cm <sup>2</sup> )                  | #            |      |      |                |                |          |        |
| Operating control conditions  | Control 1: Exam type                  |              |      |      |                |                | Nrv      | Nrv    |
|                               | Control 2: Sample volume size         |              |      |      |                |                | 1 mm     | 7 mm   |
|                               | Control 3: Sample volume position     |              |      |      |                |                | Zone 1   | Zone 0 |
|                               | Control 4: PRF                        |              |      |      |                |                | 10417    | 6250   |

(a) This index is not required for this operating mode; value is <1.  
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.  
 # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
 — Data are not applicable for this transducer/mode.

**Table 13: Transducer model: C35x Operating mode: PW Doppler**

| Index label                   |                                       |                      | M.I. | TIS  |                | TIB            | TIC    |          |
|-------------------------------|---------------------------------------|----------------------|------|------|----------------|----------------|--------|----------|
|                               |                                       |                      |      | Scan | Non-scan       |                |        | Non-scan |
|                               |                                       |                      |      |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |        |          |
| Global maximum index value    |                                       |                      | (a)  | —    | 1.5            | —              | 2.6    | (b)      |
| Associated acoustic parameter | Pr.3                                  | (MPa)                | #    |      |                |                |        |          |
|                               | $W_0$                                 | (mW)                 |      | —    | 71.1           |                | 47.1   | #        |
|                               | min of $[W_{.3}(z_1), I_{TA.3}(z_1)]$ | (mW)                 |      |      |                | —              |        |          |
|                               | $z_1$                                 | (cm)                 |      |      |                | —              |        |          |
|                               | $Z_{bp}$                              | (cm)                 |      |      |                | —              |        |          |
|                               | $Z_{sp}$                              | (cm)                 |      |      |                |                | 0.50   |          |
|                               | $z@PII_{.3max}$                       | (cm)                 | #    |      |                |                |        |          |
|                               | $d_{eq}(Z_{sp})$                      | (cm)                 |      |      |                |                | 0.36   |          |
|                               | $F_c$                                 | (MHz)                | #    | —    | 4.35           | —              | 4.37   | #        |
|                               | Dim of $A_{aprt}$                     | X (cm)               |      | —    | 1.28           | —              | 0.26   | #        |
|                               | Y (cm)                                |                      | —    | 0.80 | —              | 0.80           | #      |          |
| Other information             | PD                                    | ( $\mu$ sec)         | #    |      |                |                |        |          |
|                               | PRF                                   | (Hz)                 | #    |      |                |                |        |          |
|                               | $p_r@PII_{max}$                       | (MPa)                | #    |      |                |                |        |          |
|                               | $d_{eq}@PII_{max}$                    | (cm)                 |      |      |                |                | 0.28   |          |
|                               | Focal Length                          | $FL_x$ (cm)          |      | —    | 8.42           | —              |        | #        |
|                               |                                       | $FL_y$ (cm)          |      | —    | 5.00           | —              |        | #        |
|                               | $I_{PA.3}@MI_{max}$                   | (W/cm <sup>2</sup> ) | #    |      |                |                |        |          |
| Operating control conditions  | Control 1: Exam type                  |                      |      |      | Spine          |                | Spine  |          |
|                               | Control 2: Sample volume size         |                      |      |      | 2 mm           |                | 1 mm   |          |
|                               | Control 3: Sample volume position     |                      |      |      | Zone 5         |                | Zone 0 |          |
|                               | Control 4: PRF                        |                      |      |      | 6250           |                | 15625  |          |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 14: Transducer model: rC60xi Operating mode: PW Doppler**

| Index label                   |                                       |              | M.I.   | TIS  |                | TIB            | TIC   |          |
|-------------------------------|---------------------------------------|--------------|--------|------|----------------|----------------|-------|----------|
|                               |                                       |              |        | Scan | Non-scan       |                |       | Non-scan |
|                               |                                       |              |        |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |       |          |
| Global maximum index value    |                                       |              | 1.2    | —    | —              | 2.0            | 4.0   | (b)      |
| Associated acoustic parameter | $P_{r.3}$                             | (MPa)        | 1.73   |      |                |                |       |          |
|                               | $W_0$                                 | (mW)         |        | —    | —              |                | 291.8 | #        |
|                               | min of $[W_{.3}(z_1), I_{TA.3}(z_1)]$ | (mW)         |        |      |                | 187.5          |       |          |
|                               | $z_1$                                 | (cm)         |        |      |                | 4.0            |       |          |
|                               | $Z_{bp}$                              | (cm)         |        |      |                | 4.0            |       |          |
|                               | $Z_{sp}$                              | (cm)         |        |      |                |                | 3.60  |          |
|                               | $z@PII_{.3max}$                       | (cm)         | 4.5    |      |                |                |       |          |
|                               | $d_{eq}(Z_{sp})$                      | (cm)         |        |      |                |                | 0.94  |          |
|                               | $F_c$                                 | (MHz)        | 2.20   | —    | —              | 2.23           | 2.23  | #        |
|                               | Dim of $A_{aprt}$                     | X (cm)       |        | —    | —              | 4.77           | 3.28  | #        |
| Y (cm)                        |                                       |              | —      | —    | 1.20           | 1.20           | #     |          |
| Other information             | PD                                    | ( $\mu$ sec) | 1.153  |      |                |                |       |          |
|                               | PRF                                   | (Hz)         | 1302   |      |                |                |       |          |
|                               | $p_r@PII_{max}$                       | (MPa)        | 2.43   |      |                |                |       |          |
|                               | $d_{eq}@PII_{max}$                    | (cm)         |        |      |                |                | 0.54  |          |
|                               | Focal Length                          | $FL_x$ (cm)  |        | —    | —              | 17.97          |       | #        |
|                               |                                       | $FL_y$ (cm)  |        | —    | —              | 6.50           |       | #        |
| $I_{PA.3}@MI_{max}$           | (W/cm <sup>2</sup> )                  | 267          |        |      |                |                |       |          |
| Operating control conditions  | Control 1: Exam type                  |              | Abd    |      |                | Abd            |       | Abd      |
|                               | Control 2: Sample volume size         |              | 3 mm   |      |                | 7 mm           |       | 7 mm     |
|                               | Control 3: Sample volume position     |              | Zone 3 |      |                | Zone 6         |       | Zone 5   |
|                               | Control 4: PRF                        |              | 1302   |      |                | 2604           |       | 2604     |

(a) This index is not required for this operating mode; value is <1.  
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.  
 # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
 — Data are not applicable for this transducer/mode.

**Table 15: Transducer model: HFL38xi Operating mode: PW Doppler**

| Index label                   |                                       |              | M.I.   | TIS  |                |                | TIB      | TIC |
|-------------------------------|---------------------------------------|--------------|--------|------|----------------|----------------|----------|-----|
|                               |                                       |              |        | Scan | Non-scan       |                | Non-scan |     |
|                               |                                       |              |        |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |          |     |
| Global maximum index value    |                                       |              | 1.2    | —    | 1.1            | —              | 2.2      | (b) |
| Associated acoustic parameter | Pr.3                                  | (MPa)        | 2.69   |      |                |                |          |     |
|                               | $W_0$                                 | (mW)         |        | —    | 47.7           |                | 47.7     | #   |
|                               | min of $[W_{.3}(z_1), I_{TA.3}(z_1)]$ | (mW)         |        |      |                | —              |          |     |
|                               | $z_1$                                 | (cm)         |        |      |                | —              |          |     |
|                               | $Z_{bp}$                              | (cm)         |        |      |                | —              |          |     |
|                               | $Z_{sp}$                              | (cm)         |        |      |                |                | 1.10     |     |
|                               | $z@PII_{.3max}$                       | (cm)         | 1.0    |      |                |                |          |     |
|                               | $d_{eq}(Z_{sp})$                      | (cm)         |        |      |                |                | 0.31     |     |
|                               | $F_c$                                 | (MHz)        | 5.34   | —    | 4.86           | —              | 4.86     | #   |
|                               | Dim of $A_{aprt}$                     | X (cm)       |        | —    | 1.08           | —              | 1.08     | #   |
|                               | Y (cm)                                |              | —      | 0.40 | —              | 0.40           | #        |     |
| Other information             | PD                                    | ( $\mu$ sec) | 1.288  |      |                |                |          |     |
|                               | PRF                                   | (Hz)         | 1008   |      |                |                |          |     |
|                               | $p_r@PII_{max}$                       | (MPa)        | 3.23   |      |                |                |          |     |
|                               | $d_{eq}@PII_{max}$                    | (cm)         |        |      |                |                | 0.25     |     |
|                               | Focal Length                          | $FL_x$ (cm)  |        | —    | 3.72           | —              |          | #   |
|                               |                                       | $FL_y$ (cm)  |        | —    | 2.44           | —              |          | #   |
|                               | $I_{PA.3}@MI_{max}$                   | ( $W/cm^2$ ) | 308    |      |                |                |          |     |
| Operating control conditions  | Control 1: Exam type                  |              | Nrv    |      | Art            |                | Art      |     |
|                               | Control 2: Sample volume size         |              | 1 mm   |      | 1 mm           |                | 1 mm     |     |
|                               | Control 3: Sample volume position     |              | Zone 3 |      | Zone 7         |                | Zone 7   |     |
|                               | Control 4: PRF                        |              | 1008   |      | 3125           |                | 3125     |     |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 16: Transducer Model: HFL38xi Ophthalmic Use Operating Mode: PW Doppler**

| Index Label                          |   |                      | M.I.   | TIS  |                | TIB            | TIC    |          |   |
|--------------------------------------|---|----------------------|--------|------|----------------|----------------|--------|----------|---|
|                                      |   |                      |        | Scan | Non-scan       |                |        | Non-scan |   |
|                                      |   |                      |        |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |        |          |   |
| Global Maximum Index Value           |   |                      | 0.18   | —    | 0.09           | —              | 0.17   | (b)      |   |
| Associated Acoustic Parameter        | P <sub>r,3</sub>  | (MPa)                | 0.41   |      |                |                |        |          |   |
|                                      | W <sub>0</sub>  | (mW)                 |        | —    | 3.56           |                | 3.56   | #        |   |
|                                      | min of [W <sub>.3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] |                      | (mW)   |      |                | —              |        |          |   |
|                                      | z <sub>1</sub>  | (cm)                 |        |      |                | —              |        |          |   |
|                                      | Z <sub>bp</sub>   | (cm)                 |        |      |                | —              |        |          |   |
|                                      | Z <sub>sp</sub>   | (cm)                 |        |      |                |                | 1.64   |          |   |
|                                      | z@PII <sub>.3max</sub>  | (cm)                 | 0.9    |      |                |                |        |          |   |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |        |      |                |                | 0.31   |          |   |
|                                      | F <sub>c</sub>  | (MHz)                | 5.34   | —    | 5.33           | —              | 5.33   | #        |   |
|                                      | Dim of A <sub>aprt</sub>  | X (cm)               |        |      | —              | 1.08           | —      | 1.08     | # |
| Y (cm)                               |   |                      | —      | 0.40 | —              | 0.40           | #      |          |   |
| Other Information                    | PD  |                      | (μsec) | 1.28 |                |                |        |          |   |
|                                      | PRF   |                      | (Hz)   | 1302 |                |                |        |          |   |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  |                      | (MPa)  | 0.48 |                |                |        |          |   |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   |                      | (cm)   |      |                |                | 0.19   |          |   |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |        |      | —              | 3.72           | —      |          | # |
|                                      |   | FL <sub>y</sub> (cm) |        |      | —              | 2.44           | —      |          | # |
| I <sub>PA,3</sub> @MI <sub>max</sub> |   | (W/cm <sup>2</sup> ) | 6.6    |      |                |                |        |          |   |
| Operating Control Conditions         | Control 1: Exam Type  |                      | Oph    |      | Oph            |                | Oph    |          |   |
|                                      | Control 2: Sample Volume Size   |                      | 1 mm   |      | 10 mm          |                | 10 mm  |          |   |
|                                      | Control 3: Sample Volume Position   |                      | Zone 1 |      | Zone 7         |                | Zone 7 |          |   |
|                                      | Control 4: PRF  |                      | 1302   |      | 10417          |                | 10417  |          |   |

(a) This index is not required for this operating mode; value is <1.  
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.  
 # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
 — Data are not applicable for this transducer/mode.

**Table 17: Transducer model: HFL50x Operating mode: PW Doppler**

| Index label                   |                                       |                      | M.I.   | TIS  |                |                | TIB         | TIC |
|-------------------------------|---------------------------------------|----------------------|--------|------|----------------|----------------|-------------|-----|
|                               |                                       |                      |        | Scan | Non-scan       |                | Non-scan    |     |
|                               |                                       |                      |        |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |             |     |
| Global maximum index value    |                                       |                      | 1.2    | —    | 1.1            | —              | 1.9         | (b) |
| Associated acoustic parameter | Pr.3                                  | (MPa)                | 2.69   |      |                |                |             |     |
|                               | $W_0$                                 | (mW)                 |        | —    | 42.6           |                | 42.6        | #   |
|                               | min of $[W_{.3}(z_1), I_{TA.3}(z_1)]$ | (mW)                 |        |      |                | —              |             |     |
|                               | $z_1$                                 | (cm)                 |        |      |                | —              |             |     |
|                               | $Z_{bp}$                              | (cm)                 |        |      |                | —              |             |     |
|                               | $Z_{sp}$                              | (cm)                 | 1.0    |      |                |                | 1.1         |     |
|                               | $d_{eq}(Z_{sp})$                      | (cm)                 |        |      |                |                | 0.33        |     |
|                               | $F_c$                                 | (MHz)                | 5.34   | —    | 5.34           | —              | 5.34        | #   |
|                               | Dim of $A_{aprt}$                     | X (cm)               |        | —    | 1.08           | —              | 1.08        | #   |
|                               | Y (cm)                                |                      | —      | 0.40 | —              | 0.40           | #           |     |
| Other information             | PD                                    | ( $\mu$ sec)         | 1.29   |      |                |                |             |     |
|                               | PRF                                   | (Hz)                 | 1008   |      |                |                |             |     |
|                               | $p_r@PII_{max}$                       | (MPa)                | 3.23   |      |                |                |             |     |
|                               | $d_{eq}@PII_{max}$                    | (cm)                 |        |      |                |                | 0.22        |     |
|                               | Focal Length                          | $FL_x$ (cm)          |        | —    | 3.72           | —              |             | #   |
|                               |                                       | $FL_y$ (cm)          |        | —    | 2.44           | —              |             | #   |
|                               | $I_{PA.3}@MI_{max}$                   | (W/cm <sup>2</sup> ) | 308    |      |                |                |             |     |
| Operating control conditions  | Control 1: Exam type                  |                      | Any    | —    | Any            | —              | Any         | —   |
|                               | Control 2: Sample volume size         |                      | 1 mm   | —    | 1 mm           | —              | 1 mm        | —   |
|                               | Control 3: Sample volume position     |                      | Zone 3 | —    | Zone 7         | —              | Zone 7      | —   |
|                               | Control 4: PRF                        |                      | 1008   | —    | 1563 - 3125    | —              | 1563 - 3125 | —   |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 18: Transducer model: HSL25x Operating mode: PW Doppler**

| Index label                   |   |              | M.I. | TIS  |                | TIB            | TIC    |          |
|-------------------------------|---|--------------|------|------|----------------|----------------|--------|----------|
|                               |   |              |      | Scan | Non-scan       |                |        | Non-scan |
|                               |   |              |      |      | $A_{aprt} < 1$ | $A_{aprt} > 1$ |        |          |
| Global maximum index value    |   |              | (a)  | —    | (a)            | —              | 1.5    | (b)      |
| Associated acoustic parameter | $P_{r,3}$                               | (MPa)        | #    |      |                |                |        |          |
|                               | $W_0$                                   | (mW)         |      | —    | #              |                | 28.1   | #        |
|                               | min of [ $W_{,3}(z_1), I_{TA,3}(z_1)$ ] | (mW)         |      |      |                | —              |        |          |
|                               | $z_1$                                   | (cm)         |      |      |                | —              |        |          |
|                               | $Z_{bp}$                                | (cm)         |      |      |                | —              |        |          |
|                               | $Z_{sp}$                                | (cm)         |      |      |                |                | 0.75   |          |
|                               | $z@PII_{,3max}$                         | (cm)         | #    |      |                |                |        |          |
|                               | $d_{eq}(Z_{sp})$                        | (cm)         |      |      |                |                | 0.30   |          |
|                               | $F_c$                                   | (MHz)        | #    | —    | #              | —              | 6.00   | #        |
|                               | Dim of $A_{aprt}$                       | X (cm)       |      | —    | #              | —              | 0.76   | #        |
|                               | Y (cm)                                  |              | —    | #    | —              | 0.30           | #      |          |
| Other Information             | PD                                      | ( $\mu$ sec) | #    |      |                |                |        |          |
|                               | PRF                                     | (Hz)         | #    |      |                |                |        |          |
|                               | $p_r@PII_{max}$                         | (MPa)        | #    |      |                |                |        |          |
|                               | $d_{eq}@PII_{max}$                      | (cm)         |      |      |                |                | 0.21   |          |
|                               | Focal Length                            | $FL_x$ (cm)  |      | —    | #              | —              |        | #        |
|                               |   | $FL_y$ (cm)  |      | —    | #              | —              |        | #        |
| Operating control conditions  | Control 1: Exam type                    |              |      |      |                |                | Nrv    |          |
|                               | Control 2: Sample volume size           |              |      |      |                |                | 8 mm   |          |
|                               | Control 3: Sample volume position       |              |      |      |                |                | Zone 7 |          |
|                               | Control 4: PRF                          |              |      |      |                |                | 1953   |          |

(a) This index is not required for this operating mode; value is <1.  
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.  
 # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
 — Data are not applicable for this transducer/mode.

**Table 19: Transducer model: HSL25x Ophthalmic Use Operating mode: PW Doppler**

| Index label                          |   |                      | M.I.   | TIS  |                          | TIB                   | TIC    |          |
|--------------------------------------|---|----------------------|--------|------|--------------------------|-----------------------|--------|----------|
|                                      |   |                      |        | Scan | Non-scan                 |                       |        | Non-scan |
|                                      |   |                      |        |      | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ |        |          |
| Global maximum index value           |   |                      | 0.18   | —    | 0.12                     | —                     | 0.21   | (b)      |
| Associated acoustic parameter        | Pr,3  | (MPa)                | 0.44   |      |                          |                       |        |          |
|                                      | W <sub>0</sub>  | (mW)                 |        | —    | 4.0                      |                       | 4.0    | #        |
|                                      | min of [W <sub>,3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] | (mW)                 |        |      |                          | —                     |        |          |
|                                      | z <sub>1</sub>  | (cm)                 |        |      |                          | —                     |        |          |
|                                      | Z <sub>bp</sub>   | (cm)                 |        |      |                          | —                     |        |          |
|                                      | Z <sub>sp</sub>   | (cm)                 |        |      |                          |                       | 0.80   |          |
|                                      | z@PII <sub>,3max</sub>  | (cm)                 | 1.2    |      |                          |                       |        |          |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |        |      |                          |                       | 0.32   |          |
|                                      | F <sub>c</sub>  | (MHz)                | 6.03   | —    | 6.03                     | —                     | 6.03   | #        |
|                                      | Dim of A <sub>aprt</sub>  | X (cm)               |        | —    | 0.76                     | —                     | 0.76   | #        |
| Y (cm)                               |   |                      | —      | 0.30 | —                        | 0.30                  | #      |          |
| Other information                    | PD  | (μsec)               | 1.275  |      |                          |                       |        |          |
|                                      | PRF   | (Hz)                 | 1953   |      |                          |                       |        |          |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  | (MPa)                | 0.56   |      |                          |                       |        |          |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   | (cm)                 |        |      |                          |                       | 0.23   |          |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |        | —    | 3.80                     | —                     |        | #        |
|                                      |   | FL <sub>y</sub> (cm) |        | —    | 2.70                     | —                     |        | #        |
| I <sub>PA,3</sub> @MI <sub>max</sub> | (W/cm <sup>2</sup> )  | 7.4                  |        |      |                          |                       |        |          |
| Operating control conditions         | Control 1: Exam type  |                      | Oph    |      | Oph                      |                       | Oph    |          |
|                                      | Control 2: Sample volume size   |                      | 1 mm   |      | 1 mm                     |                       | 1 mm   |          |
|                                      | Control 3: Sample Vvolume position  |                      | Zone 7 |      | Zone 7                   |                       | Zone 7 |          |
|                                      | Control 4: PRF  |                      | 1953   |      | 5208                     |                       | 5208   |          |

(a) This index is not required for this operating mode; value is <1.  
(b) This transducer is not intended for transcranial or neonatal cephalic uses.  
# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
— Data are not applicable for this transducer/mode.

**Table 20: Transducer model: ICTx Operating mode: PW Doppler**

| Index label                   |   |              | M.I. | TIS  |                   | TIB            | TIC    |          |
|-------------------------------|---|--------------|------|------|-------------------|----------------|--------|----------|
|                               |   |              |      | Scan | Non-scan          |                |        | Non-scan |
|                               |   |              |      |      | $A_{aprt} \leq 1$ | $A_{aprt} > 1$ |        |          |
| Global maximum index value    |   |              | (a)  | —    | (a)               | —              | 1.2    | (a)      |
| Associated acoustic parameter | $P_{r,3}$                               | (MPa)        | #    | —    | —                 | —              | —      | —        |
|                               | $W_0$                                   | (mW)         | —    | —    | #                 | —              | 16.348 | #        |
|                               | min of [ $W_{,3}(z_1), I_{TA,3}(z_1)$ ] | (mW)         | —    | —    | —                 | —              | —      | —        |
|                               | $z_1$                                   | (cm)         | —    | —    | —                 | —              | —      | —        |
|                               | $Z_{bp}$                                | (cm)         | —    | —    | —                 | —              | —      | —        |
|                               | $Z_{sp}$                                | (cm)         | —    | —    | —                 | —              | 1.6    | —        |
|                               | $d_{eq}(Z_{sp})$                        | (cm)         | —    | —    | —                 | —              | 0.192  | —        |
|                               | $F_c$                                   | (MHz)        | #    | —    | #                 | —              | 4.36   | #        |
|                               | Dim of $A_{aprt}$                       | X (cm)       | —    | —    | #                 | —              | 0.6    | #        |
|                               | Y (cm)                                  | —            | —    | #    | —                 | 0.5            | #      |          |
| Other information             | PD                                      | ( $\mu$ sec) | #    | —    | —                 | —              | —      | —        |
|                               | PRF                                     | (Hz)         | #    | —    | —                 | —              | —      | —        |
|                               | $p_r@PII_{max}$                         | (MPa)        | #    | —    | —                 | —              | —      | —        |
|                               | $d_{eq}@PII_{max}$                      | (cm)         | —    | —    | —                 | —              | 0.187  | —        |
|                               | Focal Length                            | $FL_x$ (cm)  | —    | —    | #                 | —              | —      | #        |
|                               |   | $FL_y$ (cm)  | —    | —    | #                 | —              | —      | #        |
| Operating control conditions  | Control 1: Exam type                    |              | —    | —    | —                 | —              | Any    | —        |
|                               | Control 2: Sample volume size           |              | —    | —    | —                 | —              | 3 mm   | —        |
|                               | Control 3: Sample volume position       |              | —    | —    | —                 | —              | Zone 1 | —        |
|                               | Control 4: PRF                          |              | —    | —    | —                 | —              | Any    | —        |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 21: Transducer model: L25x Operating mode: PW Doppler**

| Index label                          |   |                      | M.I. | TIS  |                          |                       | TIB         | TIC |   |
|--------------------------------------|---|----------------------|------|------|--------------------------|-----------------------|-------------|-----|---|
|                                      |   |                      |      | Scan | Non-scan                 |                       | Non-scan    |     |   |
|                                      |   |                      |      |      | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ |             |     |   |
| Global maximum index value           |   |                      | (a)  | —    | (a)                      | —                     | 1.7         | (b) |   |
| Associated acoustic parameter        | Pr.3  | (MPa)                | #    |      |                          |                       |             |     |   |
|                                      | W <sub>0</sub>  | (mW)                 |      | —    | #                        |                       | 32.1        | #   |   |
|                                      | min of [W <sub>.3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] | (mW)                 |      |      |                          | —                     |             |     |   |
|                                      | z <sub>1</sub>  | (cm)                 |      |      |                          | —                     |             |     |   |
|                                      | Z <sub>bp</sub>   | (cm)                 |      |      |                          | —                     |             |     |   |
|                                      | Z <sub>sp</sub>   | (cm)                 |      |      |                          |                       | 0.75        |     |   |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |      |      |                          |                       | 0.30        |     |   |
|                                      | F <sub>c</sub>  | (MHz)                | #    | —    | #                        | —                     | 6.00        | #   |   |
|                                      | Dim of A <sub>aprt</sub>  | X (cm)               |      | —    | #                        | —                     | 0.76        | #   |   |
|                                      | Y (cm)  |                      | —    | #    | —                        | 0.30                  | #           |     |   |
| Other Information                    | PD  | (μsec)               | #    |      |                          |                       |             |     |   |
|                                      | PRF   | (Hz)                 | #    |      |                          |                       |             |     |   |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  | (MPa)                | #    |      |                          |                       |             |     |   |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   | (cm)                 |      |      |                          |                       | 0.21        |     |   |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |      | —    | #                        | —                     |             |     | # |
|                                      |   | FL <sub>y</sub> (cm) |      | —    | #                        | —                     |             |     | # |
| I <sub>PA,3</sub> @MI <sub>max</sub> | (W/cm <sup>2</sup> )  | #                    |      |      |                          |                       |             |     |   |
| Operating control conditions         | Control 1: Exam type  |                      | —    | —    | —                        | —                     | Vas/Ven/Nrv | —   |   |
|                                      | Control 2: Sample volume size   |                      | —    | —    | —                        | —                     | 8 mm        | —   |   |
|                                      | Control 3: Sample volume position   |                      | —    | —    | —                        | —                     | Zone 7      | —   |   |
|                                      | Control 4: PRF  |                      | —    | —    | —                        | —                     | 1953        | —   |   |

(a) This index is not required for this operating mode; value is <1.  
(b) This transducer is not intended for transcranial or neonatal cephalic uses.  
# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
— Data are not applicable for this transducer/mode.

**Table 22: Transducer model: L25x Ophthalmic Use Operating mode: PW Doppler**

| Index label                   |                                       |              | M.I.   | TIS  |                   | TIB            | TIC    |          |
|-------------------------------|---------------------------------------|--------------|--------|------|-------------------|----------------|--------|----------|
|                               |                                       |              |        | Scan | Non-scan          |                |        | Non-scan |
|                               |                                       |              |        |      | $A_{aprt} \leq 1$ | $A_{aprt} > 1$ |        |          |
| Global maximum index value    |                                       |              | 0.18   | —    | 0.12              | —              | 0.21   | (b)      |
| Associated acoustic parameter | $P_{r,3}$                             | (MPa)        | 0.44   |      |                   |                |        |          |
|                               | $W_0$                                 | (mW)         |        | —    | 4.0               |                | 4.0    | #        |
|                               | min of $[W_{,3}(z_1), I_{TA,3}(z_1)]$ | (mW)         |        |      |                   | —              |        |          |
|                               | $z_1$                                 | (cm)         |        |      |                   | —              |        |          |
|                               | $Z_{bp}$                              | (cm)         |        |      |                   | —              |        |          |
|                               | $Z_{sp}$                              | (cm)         |        |      |                   |                | 0.80   |          |
|                               | $z@PII_{,3max}$                       | (cm)         | 1.2    |      |                   |                |        |          |
|                               | $d_{eq}(Z_{sp})$                      | (cm)         |        |      |                   |                | 0.32   |          |
|                               | $F_c$                                 | (MHz)        | 6.03   | —    | 6.03              | —              | 6.03   | #        |
|                               | Dim of $A_{aprt}$                     | X (cm)       |        | —    | 0.76              | —              | 0.76   | #        |
| Y (cm)                        |                                       |              | —      | 0.30 | —                 | 0.30           | #      |          |
| Other information             | PD                                    | ( $\mu$ sec) | 1.275  |      |                   |                |        |          |
|                               | PRF                                   | (Hz)         | 1953   |      |                   |                |        |          |
|                               | $p_r@PII_{max}$                       | (MPa)        | 0.56   |      |                   |                |        |          |
|                               | $d_{eq}@PII_{max}$                    | (cm)         |        |      |                   |                | 0.23   |          |
|                               | Focal Length                          | $FL_x$ (cm)  |        | —    | 3.80              | —              |        | #        |
|                               |                                       | $FL_y$ (cm)  |        | —    | 2.70              | —              |        | #        |
| $I_{PA,3}@MI_{max}$           | (W/cm <sup>2</sup> )                  | 7.4          |        |      |                   |                |        |          |
| Operating control conditions  | Control 1: Exam type                  |              | Oph    |      | Oph               |                | Oph    |          |
|                               | Control 2: Sample volume size         |              | 1 mm   |      | 1 mm              |                | 1 mm   |          |
|                               | Control 3: Sample volume position     |              | Zone 7 |      | Zone 7            |                | Zone 7 |          |
|                               | Control 4: PRF                        |              | 1953   |      | 5208              |                | 5208   |          |

(a) This index is not required for this operating mode; value is <1.  
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.  
 # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
 — Data are not applicable for this transducer/mode.

**Table 23: Transducer model: L38xi Operating mode: PW Doppler**

| Index label                          |   |                      | M.I.   | TIS  |                          | TIB                   | TIC    |          |  |
|--------------------------------------|---|----------------------|--------|------|--------------------------|-----------------------|--------|----------|--|
|                                      |   |                      |        | Scan | Non-scan                 |                       |        | Non-scan |  |
|                                      |   |                      |        |      | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ |        |          |  |
| Global maximum index value           |   |                      | 1.3    | —    | 2.6                      | —                     | 3.7    | (b)      |  |
| Associated acoustic parameter        | Pr.3  | (MPa)                | 2.59   |      |                          |                       |        |          |  |
|                                      | W <sub>0</sub>  | (mW)                 |        | —    | 114.5                    |                       | 114.5  | #        |  |
|                                      | min of [W <sub>.3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] |                      | (mW)   |      |                          |                       | —      |          |  |
|                                      | z <sub>1</sub>  | (cm)                 |        |      |                          | —                     |        |          |  |
|                                      | Z <sub>bp</sub>   | (cm)                 |        |      |                          | —                     |        |          |  |
|                                      | Z <sub>sp</sub>   | (cm)                 |        |      |                          |                       | 1.20   |          |  |
|                                      | z@PII <sub>.3max</sub>  | (cm)                 | 0.7    |      |                          |                       |        |          |  |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |        |      |                          |                       | 0.32   |          |  |
|                                      | F <sub>c</sub>  | (MHz)                | 4.06   | —    | 4.78                     | —                     | 4.78   | #        |  |
|                                      | Dim of A <sub>aprt</sub>  | X (cm)               |        | —    | 1.86                     | —                     | 1.86   | #        |  |
| Y (cm)                               |   |                      | —      | 0.40 | —                        | 0.40                  | #      |          |  |
| Other information                    | PD  | (μsec)               | 1.230  |      |                          |                       |        |          |  |
|                                      | PRF   | (Hz)                 | 1008   |      |                          |                       |        |          |  |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  | (MPa)                | 2.86   |      |                          |                       |        |          |  |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   | (cm)                 |        |      |                          |                       | 0.46   |          |  |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |        | —    | 5.54                     | —                     |        | #        |  |
|                                      |   | FL <sub>y</sub> (cm) |        | —    | 1.50                     | —                     |        | #        |  |
| I <sub>PA,3</sub> @MI <sub>max</sub> | (W/cm <sup>2</sup> )  | 323                  |        |      |                          |                       |        |          |  |
| Operating control conditions         | Control 1: Exam type  |                      | Art    |      | Nrv                      |                       | Nrv    |          |  |
|                                      | Control 2: Sample volume size   |                      | 1 mm   |      | 1 mm                     |                       | 1 mm   |          |  |
|                                      | Control 3: Sample volume position   |                      | Zone 0 |      | Zone 7                   |                       | Zone 7 |          |  |
|                                      | Control 4: PRF  |                      | 1008   |      | 10417                    |                       | 10417  |          |  |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 24: Transducer model: P10x Operating mode: PW Doppler**

| Index label   |  |                      | M.I. | TIS    |                   | TIB            | TIC    |          |
|---|--|----------------------|------|--------|-------------------|----------------|--------|----------|
|   |  |                      |      | Scan   | Non-scan          |                |        | Non-scan |
|   |  |                      |      |        | $A_{aprt} \leq 1$ | $A_{aprt} > 1$ |        |          |
| Global maximum index value                                |  |                      | 1.0  | —      | 1.1               | —              | 1.9    | 1.5      |
| Associated acoustic parameter                             | Pr,3 (MPa)   |                      | 1.92 |        |                   |                |        |          |
|   | W <sub>0</sub> (mW)  |                      |      | —      | 34.4              |                | 31.9   | 26.9     |
|   | min of [W <sub>,3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] (mW) |                      |      |        |                   | —              |        |          |
|   | z <sub>1</sub> (cm)  |                      |      |        |                   | —              |        |          |
|   | Z <sub>bp</sub> (cm)   |                      |      |        |                   | —              |        |          |
|   | Z <sub>sp</sub> (cm)   |                      |      |        |                   |                | 0.80   |          |
|   | z@PII <sub>,3max</sub> (cm)  | 2.1                  |      |        |                   |                |        |          |
|   | d <sub>eq</sub> (Z <sub>sp</sub> ) (cm)  |                      |      |        |                   |                | 0.31   |          |
|   | F <sub>c</sub> (MHz)   | 3.87                 | —    | 6.86   | —                 | 3.84           | 6.86   |          |
|   | Dim of A <sub>aprt</sub>   | X (cm)               |      | —      | 0.99              | —              | 0.42   | 0.22     |
| Y (cm)  |  |                      | —    | 0.70   | —                 | 0.70           | 0.70   |          |
| Other information   | PD (μsec)  | 1.277                |      |        |                   |                |        |          |
|   | PRF (Hz)   | 1562                 |      |        |                   |                |        |          |
|   | p <sub>r</sub> @PII <sub>max</sub> (MPa)   | 2.54                 |      |        |                   |                |        |          |
|   | d <sub>eq</sub> @PII <sub>max</sub> (cm)   |                      |      |        |                   | 0.24           |        |          |
|   | Focal Length   | FL <sub>x</sub> (cm) |      | —      | 6.74              | —              |        | 0.92     |
|   |  | FL <sub>y</sub> (cm) |      | —      | 5.00              | —              |        | 5.00     |
| I <sub>PA,3</sub> @MI <sub>max</sub> (W/cm <sup>2</sup> ) | 200  |                      |      |        |                   |                |        |          |
| Operating control conditions                              | Control 1: Exam type   | Crd                  |      | Crd    |                   | Abd            | Crd    |          |
|   | Control 2: Sample volume size  | 1 mm                 |      | 7 mm   |                   | 12 mm          | 1 mm   |          |
|   | Control 3: Sample volume position  | Zone 2               |      | Zone 6 |                   | Zone 1         | Zone 0 |          |
|   | Control 4: PRF   | 1562                 |      | 1008   |                   | 1953           | 15625  |          |
|   | Control 5: TDI   | Off                  |      | On     |                   | Off            | Off    |          |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 25: Transducer Model: P10x Operating Mode: CW Doppler**

| Index label                          |   |                      | M.I. | TIS  |                          |                       | TIB      | TIC    |
|--------------------------------------|---|----------------------|------|------|--------------------------|-----------------------|----------|--------|
|                                      |   |                      |      | Scan | Non-scan                 |                       | Non-scan |        |
|                                      |   |                      |      |      | $A_{\text{aprt}} \leq 1$ | $A_{\text{aprt}} > 1$ |          |        |
| Global maximum index value           |   |                      | (a)  | —    | (a)                      | —                     | 1.8      | 1.7    |
| Associated acoustic parameter        | Pr.3  | (MPa)                | 2.59 |      |                          |                       |          |        |
|                                      | W <sub>0</sub>  | (mW)                 |      | —    | #                        |                       | 34.8     | 25.7   |
|                                      | min of [W <sub>.3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] | (mW)                 |      |      |                          | —                     |          |        |
|                                      | z <sub>1</sub>  | (cm)                 |      |      |                          | —                     |          |        |
|                                      | Z <sub>bp</sub>   | (cm)                 |      |      |                          | —                     |          |        |
|                                      | Z <sub>sp</sub>   | (cm)                 |      |      |                          |                       | 0.70     |        |
|                                      | z@PII <sub>.3max</sub>  | (cm)                 | #    |      |                          |                       |          |        |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |      |      |                          |                       | 0.36     |        |
|                                      | F <sub>c</sub>  | (MHz)                | #    | —    | #                        | —                     | 4.00     | 4.00   |
| Dim of A <sub>aprt</sub>             | X (cm)  |                      | —    | #    | —                        | 0.32                  | 0.16     |        |
|                                      | Y (cm)  |                      | —    | #    | —                        | 0.70                  | 0.70     |        |
| Other information                    | PD  | (μsec)               | #    |      |                          |                       |          |        |
|                                      | PRF   | (Hz)                 | #    |      |                          |                       |          |        |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  | (MPa)                | #    |      |                          |                       |          |        |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   | (cm)                 |      |      |                          |                       | 0.27     |        |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |      | —    | #                        | —                     |          | 0.92   |
|                                      |   | FL <sub>y</sub> (cm) |      | —    | #                        | —                     |          | 5.00   |
| I <sub>PA,3</sub> @MI <sub>max</sub> | (W/cm <sup>2</sup> )  | #                    |      |      |                          |                       |          |        |
| Operating control conditions         | Control 1: Exam type  |                      |      |      |                          |                       | Crd      | Crd    |
|                                      | Control 2: Sample volume position   |                      |      |      |                          |                       |          | Zone 0 |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 26: Transducer model: rP19x Operating mode: PW Doppler**

| Index label                          |   |                      | M.I.    | TIS  |                   | TIB            | TIC     |          |  |
|--------------------------------------|---|----------------------|---------|------|-------------------|----------------|---------|----------|--|
|                                      |   |                      |         | Scan | Non-scan          |                |         | Non-scan |  |
|                                      |   |                      |         |      | $A_{aprt} \leq 1$ | $A_{aprt} > 1$ |         |          |  |
| Global maximum index value           |   |                      | 1.3     | —    | —                 | 1.8            | 4.0     | 3.9      |  |
| Associated acoustic parameter        | P <sub>r,3</sub>  | (MPa)                | 1.94    |      |                   |                |         |          |  |
|                                      | W <sub>0</sub>  | (mW)                 |         | —    | —                 |                | 240.2   | 251.1    |  |
|                                      | min of [W <sub>,3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] |                      | (mW)    |      |                   |                | 173.7   |          |  |
|                                      | z <sub>1</sub>  | (cm)                 |         |      |                   | 2.5            |         |          |  |
|                                      | Z <sub>bp</sub>   | (cm)                 |         |      |                   | 2.5            |         |          |  |
|                                      | Z <sub>sp</sub>   | (cm)                 |         |      |                   |                | 3.35    |          |  |
|                                      | z@PII <sub>,3max</sub>  | (cm)                 | 3.0     |      |                   |                |         |          |  |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |         |      |                   |                | 0.80    |          |  |
|                                      | F <sub>c</sub>  | (MHz)                | 2.14    | —    | —                 | 2.23           | 2.23    | 2.10     |  |
|                                      | Dim of A <sub>aprt</sub>  | X (cm)               |         | —    | —                 | 1.86           | 1.80    | 1.80     |  |
| Y (cm)                               |   |                      | —       | —    | 1.15              | 1.15           | 1.15    |          |  |
| Other information                    | PD  | (µsec)               | 1.334   |      |                   |                |         |          |  |
|                                      | PRF   | (Hz)                 | 1562    |      |                   |                |         |          |  |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  | (MPa)                | 2.42    |      |                   |                |         |          |  |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   | (cm)                 |         |      |                   |                | 0.62    |          |  |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |         | —    | —                 | 29.82          |         | 18.46    |  |
|                                      |   | FL <sub>y</sub> (cm) |         | —    | —                 | 9.00           |         | 9.00     |  |
| I <sub>PA,3</sub> @MI <sub>max</sub> | (W/cm <sup>2</sup> )  | 180                  |         |      |                   |                |         |          |  |
| Operating control conditions         | Control 1: Exam type  |                      | Crd     |      |                   | Crd            | Crd     | Crd      |  |
|                                      | Control 2: Sample volume size   |                      | 1 mm    |      |                   | 12 mm          | 12 mm   | 1 mm     |  |
|                                      | Control 3: Sample volume position   |                      | Zone 1  |      |                   | Zone 7         | Zone 5  | Zone 5   |  |
|                                      | Control 4: PRF  |                      | 1562 Hz |      |                   | 1562 Hz        | 1562 Hz | 3906 Hz  |  |
|                                      | Control 5: TDI  |                      | Off     |      |                   | Off            | Off     | Off      |  |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 27: Transducer model: rP19x Orbital Use Operating mode: PW Doppler**

| Index label                   |   |              | M.I.   | TIS  |                   |                | TIB      | TIC    |
|-------------------------------|---|--------------|--------|------|-------------------|----------------|----------|--------|
|                               |   |              |        | Scan | Non-scan          |                | Non-scan |        |
|                               |   |              |        |      | $A_{aprt} \leq 1$ | $A_{aprt} > 1$ |          |        |
| Global maximum index value    |   |              | 0.18   | —    | —                 | 0.27           | 0.59     | 0.57   |
| Associated acoustic parameter | $P_{r,3}$                               | (MPa)        | 0.27   |      |                   |                |          |        |
|                               | $W_0$                                   | (mW)         |        | —    | —                 |                | 35.3     | 37.4   |
|                               | min of [ $W_{.3}(z_1), I_{TA,3}(z_1)$ ] | (mW)         |        |      |                   | 25.3           |          |        |
|                               | $z_1$                                   | (cm)         |        |      |                   | 2.5            |          |        |
|                               | $Z_{bp}$                                | (cm)         |        |      |                   | 2.5            |          |        |
|                               | $Z_{sp}$                                | (cm)         |        |      |                   |                | 3.35     |        |
|                               | $z@PII_{.3max}$                         | (cm)         | 3.5    |      |                   |                |          |        |
|                               | $d_{eq}(Z_{sp})$                        | (cm)         |        |      |                   |                | 0.80     |        |
|                               | $F_c$                                   | (MHz)        | 2.23   | —    | —                 | 2.23           | 2.23     | 2.23   |
|                               | Dim of $A_{aprt}$                       | X (cm)       |        | —    | —                 | 1.86           | 1.80     | 1.86   |
| Y (cm)                        |   |              | —      | —    | 1.15              | 1.15           | 1.15     |        |
| Other information             | PD                                      | ( $\mu$ sec) | 6.557  |      |                   |                |          |        |
|                               | PRF                                     | (Hz)         | 1953   |      |                   |                |          |        |
|                               | $p_r@PII_{max}$                         | (MPa)        | 0.36   |      |                   |                |          |        |
|                               | $d_{eq}@PII_{max}$                      | (cm)         |        |      |                   |                | 0.64     |        |
|                               | Focal Length                            | $FL_x$ (cm)  |        | —    | —                 | 29.82          |          | 29.82  |
|                               |   | $FL_y$ (cm)  |        | —    | —                 | 9.00           |          | 9.00   |
| $I_{PA,3}@MI_{max}$           | ( $W/cm^2$ )                            | 2.49         |        |      |                   |                |          |        |
| Operating control conditions  | Control 1: Exam type                    |              | Orb    |      |                   | Orb            | Orb      | Orb    |
|                               | Control 2: Sample volume size           |              | 5 mm   |      |                   | 14 mm          | 14 mm    | 14 mm  |
|                               | Control 3: Sample volume position       |              | Zone 6 |      |                   | Zone 7         | Zone 5   | Zone 7 |
|                               | Control 4: PRF                          |              | 1953   |      |                   | 1953           | 1953     | 1953   |

(a) This index is not required for this operating mode; value is <1.

(b) This transducer is not intended for transcranial or neonatal cephalic uses.

# No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)

— Data are not applicable for this transducer/mode.

**Table 28: Transducer model: rP19x Operating mode: CW Doppler**

| Index label                          |   |                      | M.I. | TIS  |                   | TIB            | TIC    |          |
|--------------------------------------|---|----------------------|------|------|-------------------|----------------|--------|----------|
|                                      |   |                      |      | Scan | Non-scan          |                |        | Non-scan |
|                                      |   |                      |      |      | $A_{aprt} \leq 1$ | $A_{aprt} > 1$ |        |          |
| Global maximum index value           |   |                      | (a)  | —    | 1.2               | —              | 4.0    | 4.0      |
| Associated acoustic parameter        | Pr,3  | (MPa)                | #    |      |                   |                |        |          |
|                                      | W <sub>0</sub>  | (mW)                 |      | —    | 125.4             |                | 125.4  | 125.4    |
|                                      | min of [W <sub>,3</sub> (z <sub>1</sub> ), I <sub>TA,3</sub> (z <sub>1</sub> )] | (mW)                 |      |      |                   | —              |        |          |
|                                      | z <sub>1</sub>  | (cm)                 |      |      |                   | —              |        |          |
|                                      | Z <sub>bp</sub>   | (cm)                 |      |      |                   | —              |        |          |
|                                      | Z <sub>sp</sub>   | (cm)                 |      |      |                   |                | 0.90   |          |
|                                      | z@PII <sub>,3max</sub>  | (cm)                 | #    |      |                   |                |        |          |
|                                      | d <sub>eq</sub> (Z <sub>sp</sub> )  | (cm)                 |      |      |                   |                | 0.64   |          |
|                                      | F <sub>c</sub>  | (MHz)                | #    | —    | 2.00              | —              | 2.00   | 2.00     |
|                                      | Dim of A <sub>aprt</sub>  | X (cm)               |      | —    | 0.42              | —              | 0.42   | 0.42     |
| Y (cm)                               |   |                      | —    | 1.15 | —                 | 1.15           | 1.15   |          |
| Other information                    | PD  | (µsec)               | #    |      |                   |                |        |          |
|                                      | PRF   | (Hz)                 | #    |      |                   |                |        |          |
|                                      | p <sub>r</sub> @PII <sub>max</sub>  | (MPa)                | #    |      |                   |                |        |          |
|                                      | d <sub>eq</sub> @PII <sub>max</sub>   | (cm)                 |      |      |                   |                | 0.61   |          |
|                                      | Focal Length  | FL <sub>x</sub> (cm) |      | —    | 1.55              | —              |        | 1.55     |
|                                      |   | FL <sub>y</sub> (cm) |      | —    | 9.00              | —              |        | 9.00     |
| I <sub>PA,3</sub> @MI <sub>max</sub> | (W/cm <sup>2</sup> )  | #                    |      |      |                   |                |        |          |
| Operating control conditions         | Control 1: Exam type  |                      |      |      | Crd               |                | Crd    | Crd      |
|                                      | Control 2: Sample volume position   |                      |      |      | Zone 0            |                | Zone 0 | Zone 0   |

(a) This index is not required for this operating mode; value is <1.  
 (b) This transducer is not intended for transcranial or neonatal cephalic uses.  
 # No data are reported for this operating condition since the global maximum index value is not reported for the reason listed. (Reference Global Maximum Index Value line.)  
 — Data are not applicable for this transducer/mode.

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| Transducer | Undersøgel-<br>sestype <sup>a</sup> | Billeddannelsestype       |                  |                    |                            |               |
|------------|-------------------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|            |                                     | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farve <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| C8x        | Pro                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
| C11x       | Abd                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Art                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Neo                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Nrv                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Ven                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
| C35x       | Abd                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Msk                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Nrv                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | OB                                  | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Spn                                 | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup> Forkortelser for undersøgelsestyper er som følger: Abd = Abdomen, Art = Arteriel, Bre = Bryst, Crd = Hjerte, Gyn = Gynækologi, Msk = Muskuloskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetriske, Oph = Oftalmisk, Orb = Orbital, SmP = Mindre organer, Sup = Superficiel, TCD = Transkraniel doppler, Ven = Venøs.

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| Transducer                  | Undersøgelsestype <sup>a</sup> | Billeddannelsestype       |                  |                    |                            |               |
|-----------------------------|--------------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                             |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farve <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| rC60xi standard/<br>armeret | Abd                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Gyn                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Msk                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Nrv                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | OB                             | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup> Forkortelser for undersøgelsestyper er som følger: Abd = Abdomen, Art = Arteriel, Bre = Bryst, Crd = Hjerte, Gyn = Gynækologi, Msk = Muskuloskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrisk, Oph = Oftalmisk, Orb = Orbital, SmP = Mindre organer, Sup = Superficiel, TCD = Transkraniel doppler, Ven = Venøs.

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| Transducer                  | Undersøgel-<br>sestype <sup>a</sup> | Billeddannelsestype       |                  |                    |                            |               |
|-----------------------------|-------------------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                             |                                     | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farve <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| HFL38xi<br>standard/armeret | Art                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Bre                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Lung                                | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Msk                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Nrv                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Oph                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | SmP                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Ven                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
| HFL50x                      | Bre                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Msk                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | Nrv                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                             | SmP                                 | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup> Forkortelser for undersøgelsestyper er som følger: Abd = Abdomen, Art = Arteriel, Bre = Bryst, Crd = Hjerte, Gyn = Gynækologi, Msk = Muskuloskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetriske, Oph = Oftalmisk, Orb = Orbital, SmP = Mindre organer, Sup = Superficiel, TCD = Transkraniel doppler, Ven = Venøs.

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| Transducer | Undersøgelsestype <sup>a</sup> | Billeddannelsestype       |                  |                    |                            |               |
|------------|--------------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|            |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farve <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| HSL25x     | Art                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Lung                           | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Msk                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Nrv                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Oph                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Sup                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | Ven                            | ✓                         | ✓                | ✓                  | ✓                          |               |
| ICTx       | Gyn                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|            | OB                             | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup> Forkortelser for undersøgelsestyper er som følger: Abd = Abdomen, Art = Arteriel, Bre = Bryst, Crd = Hjerte, Gyn = Gynækologi, Msk = Muskuloskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrisk, Oph = Oftalmisk, Orb = Orbital, SmP = Mindre organer, Sup = Superficiel, TCD = Transkraniel doppler, Ven = Venøs.

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| Transducer                 | Undersøgel-<br>sestype <sup>a</sup> | Billeddannelsestype       |                  |                    |                            |               |
|----------------------------|-------------------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                            |                                     | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farve <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| L25x standard/<br>armeret  | Art                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Lung                                | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Msk                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Nrv                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Oph                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Sup                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Ven                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
| L38xi standard/<br>armeret | Art                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Lung                                | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Nrv                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | SmP                                 | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Ven                                 | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup> Forkortelser for undersøgelsestyper er som følger: Abd = Abdomen, Art = Arteriel, Bre = Bryst, Crd = Hjerte, Gyn = Gynækologi, Msk = Muskuloskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetriske, Oph = Oftalmisk, Orb = Orbital, SmP = Mindre organer, Sup = Superficiel, TCD = Transkraniel doppler, Ven = Venøs.

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| Transducer                 | Undersøgelsestype <sup>a</sup> | Billeddannelsestype       |                  |                    |                            |               |
|----------------------------|--------------------------------|---------------------------|------------------|--------------------|----------------------------|---------------|
|                            |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farve <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| P10x                       | Abd                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Crđ                            | ✓                         |                  | ✓                  | ✓                          | ✓             |
|                            | Neo                            | ✓                         | ✓                | ✓                  | ✓                          |               |
| P11x <sup>e</sup>          | Art                            | ✓                         | ✓                | ✓                  |                            |               |
|                            | Ven                            | ✓                         | ✓                | ✓                  |                            |               |
| rP19x standard/<br>armeret | Abd                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Crđ                            | ✓                         |                  | ✓                  | ✓                          | ✓             |
|                            | Lung                           | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | OB                             | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | Orb                            | ✓                         | ✓                | ✓                  | ✓                          |               |
|                            | TCD                            | ✓                         | ✓                | ✓                  | ✓                          |               |

<sup>a</sup> Forkortelser for undersøgelsestyper er som følger: Abd = Abdomen, Art = Arteriel, Bre = Bryst, Crđ = Hjerte, Gyn = Gynækologi, Msk = Muskuloskeletal, Neo = Neonatal, Nrv = Nerve, OB = Obstetrisk, Oph = Oftalmisk, Orb = Orbital, SmP = Mindre organer, Sup = Superficiel, TCD = Transkraniel doppler, Ven = Venøs.

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| Transduser | Undersøkelsestype <sup>a</sup> | Avbildningsmodus          |                  |                    |                         |            |
|------------|--------------------------------|---------------------------|------------------|--------------------|-------------------------|------------|
|            |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-doppler <sup>d</sup> | CW-doppler |
| C8x        | Pro                            | ✓                         | ✓                | ✓                  | ✓                       |            |
| C11x       | Abd                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|            | Art                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|            | Neo                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|            | Nrv                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|            | Ven                            | ✓                         | ✓                | ✓                  | ✓                       |            |

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| Transduser                 | Undersøkelsestype <sup>a</sup> | Avbildningsmodus          |                  |                    |                         |            |
|----------------------------|--------------------------------|---------------------------|------------------|--------------------|-------------------------|------------|
|                            |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-doppler <sup>d</sup> | CW-doppler |
| C35x                       | Abd                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | Msk                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | Nrv                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | OB                             | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | Spn                            | ✓                         | ✓                | ✓                  | ✓                       |            |
| rC60xi standard/<br>armert | Abd                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | Gyn                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | Msk                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | Nrv                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                            | OB                             | ✓                         | ✓                | ✓                  | ✓                       |            |

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| Transduser                  | Under-<br>søkelses-<br>type <sup>a</sup> | Avbildningsmodus          |                  |                    |                             |                |
|-----------------------------|--|---------------------------|------------------|--------------------|-----------------------------|----------------|
|                             |  | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-<br>doppler <sup>d</sup> | CW-<br>doppler |
| HFL38xi standard/<br>armert | Art                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | Bre                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | Lung                                     | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | Msk                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | Nrv                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | Oph                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | SmP                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|                             | Ven                                      | ✓                         | ✓                | ✓                  | ✓                           |                |

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| Transduser | Under-<br>søkelses-<br>type <sup>a</sup> | Avbildningsmodus          |                  |                    |                             |                |
|------------|--|---------------------------|------------------|--------------------|-----------------------------|----------------|
|            |  | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-<br>doppler <sup>d</sup> | CW-<br>doppler |
| HFL50x     | Bre                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Msk                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Nrv                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | SmP                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
| HSL25x     | Art                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Lung                                     | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Msk                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Nrv                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Oph                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Sup                                      | ✓                         | ✓                | ✓                  | ✓                           |                |
|            | Ven                                      | ✓                         | ✓                | ✓                  | ✓                           |                |

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| Transduser               | Undersøkelsestype <sup>a</sup> | Avbildningsmodus          |                  |                    |                         |            |
|--------------------------|--------------------------------|---------------------------|------------------|--------------------|-------------------------|------------|
|                          |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-doppler <sup>d</sup> | CW-doppler |
| ICTx                     | Gyn                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | OB                             | ✓                         | ✓                | ✓                  | ✓                       |            |
| L25x standard/<br>armert | Art                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | Lung                           | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | Msk                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | Nrv                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | Oph                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | Sup                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                          | Ven                            | ✓                         | ✓                | ✓                  | ✓                       |            |

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| Transduser                | Undersøkelsestype <sup>a</sup> | Avbildningsmodus          |                  |                    |                         |            |
|---------------------------|--------------------------------|---------------------------|------------------|--------------------|-------------------------|------------|
|                           |                                | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-doppler <sup>d</sup> | CW-doppler |
| L38xi standard/<br>armert | Art                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | Lung                           | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | Nrv                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | SmP                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | Ven                            | ✓                         | ✓                | ✓                  | ✓                       |            |
| P10x                      | Abd                            | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | Crd                            | ✓                         |                  | ✓                  | ✓                       | ✓          |
|                           | Neo                            | ✓                         | ✓                | ✓                  | ✓                       |            |
| P11x <sup>e</sup>         | Art                            | ✓                         | ✓                | ✓                  |                         |            |
|                           | Ven                            | ✓                         | ✓                | ✓                  |                         |            |

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| Transduser                | Under-søkelses-type <sup>a</sup> | Avbildningsmodus          |                  |                    |                         |            |
|---------------------------|----------------------------------|---------------------------|------------------|--------------------|-------------------------|------------|
|                           |                                  | 2D <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Farge <sup>c</sup> | PW-doppler <sup>d</sup> | CW-doppler |
| rP19x standard/<br>armert | Abd                              | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | Crd                              | ✓                         |                  | ✓                  | ✓                       | ✓          |
|                           | Lung                             | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | OB                               | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | Orb                              | ✓                         | ✓                | ✓                  | ✓                       |            |
|                           | TCD                              | ✓                         | ✓                | ✓                  | ✓                       |            |

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$$RV = ERO * MR \text{ VTI}/100$$

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Lang, R., M. Bierig, R. Devereux, et al. "Recommendations for chamber quantification: a report from the American Society of Echocardiography's guidelines and standards committee and the chamber quantification writing group, Developed in conjunction with the European Association of Echocardiography, a © of the European Society of ©." *Journal of the American Society of Echocardiograph*. 2005, 18: p.1440-1463.

$$RA \text{ Vol} = \pi/4 * \sum(ai) * ai * L/20 \text{ for } i = 1 \text{ til } 20 \text{ (antall segmenter)}$$

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## TAPSE

Rudski, L., W. Lai, et al. "Guidelines for the echocardiographic assessment of the right heart in adults: a report from the American Society of Echocardiography." *Journal of the American Society of Echocardiograph*. (2010), p.685-713.

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$$RV = ERO * MR VTI / 100$$

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Lang, R., M. Bierig, R. Devereux, et al. "Recommendations for chamber quantification: a report from the American Society of Echocardiography's guidelines and standards committee and the chamber quantification writing group, Developed in conjunction with the European Association of Echocardiography, a branch of the European Society of Cardiology." *Journal of the American Society of Echocardiograph*. 2005, 18: p. 1440–1463.

$$RA \text{ Vol} = \pi/4 * \sum(ai) * ai * L/20 \text{ för } i = 1 \text{ till } 20 \text{ (antal segment)}$$

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Wang, Y., J. Gutman, et al. "Atrial volume in a normal adult population by two-dimensional echocardiography". *Chest*. (1984), 86: p. 595–601.

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Lang, R., M. Bierig, R. Devereux, et al. "Recommendations for chamber quantification: a report from the American Society of Echocardiography's guidelines and standards committee and the chamber quantification writing group, Developed in conjunction with the European Association of Echocardiography, a branch of the European Society of Cardiology." *Journal of the American Society of Echocardiography*. 2005, 18: p.1440-1463.

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Wang Y., Gutman J., et al. "Atrial volume in a normal adult population by two-dimensional echocardiography." *Chest*. (1984), 86: p.595-601.

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## TAPSE

Rudski L., Lai W., et al. "Guidelines for the echocardiographic assessment of the right heart in adults: a report from the American Society of Echocardiography." *Journal of the American Society of Echocardiography*. (2010), p.685-713.

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| Dönüştürücü                   | Muayene türü <sup>a</sup> | Görüntüleme modu          |                  |                     |                         |            |
|-------------------------------|---------------------------|---------------------------|------------------|---------------------|-------------------------|------------|
|                               |                           | 2B <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Renkli <sup>c</sup> | PW Doppler <sup>d</sup> | CW Doppler |
| HFL38xi standart/<br>korumalı | Art                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Bre                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Lung                      | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Msk                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Nrv                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Oph                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | SmP                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Ven                       | ✓                         | ✓                | ✓                   | ✓                       |            |
| HFL50x                        | Bre                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Msk                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | Nrv                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                               | SmP                       | ✓                         | ✓                | ✓                   | ✓                       |            |

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| Dönüştürücü | Muayene türü <sup>a</sup> | Görüntüleme modu          |                  |                     |                            |               |
|-------------|---------------------------|---------------------------|------------------|---------------------|----------------------------|---------------|
|             |                           | 2B <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Renkli <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| HSL25x      | Art                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | Lung                      | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | Msk                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | Nrv                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | Oph                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | Sup                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | Ven                       | ✓                         | ✓                | ✓                   | ✓                          |               |
| ICTx        | Gyn                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|             | OB                        | ✓                         | ✓                | ✓                   | ✓                          |               |

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| Dönüştürücü                 | Muayene türü <sup>a</sup> | Görüntüleme modu          |                  |                     |                         |            |
|-----------------------------|---------------------------|---------------------------|------------------|---------------------|-------------------------|------------|
|                             |                           | 2B <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Renkli <sup>c</sup> | PW Doppler <sup>d</sup> | CW Doppler |
| L25x standart/<br>korumalı  | Art                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Lung                      | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Msk                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Nrv                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Oph                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Sup                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Ven                       | ✓                         | ✓                | ✓                   | ✓                       |            |
| L38xi standart/<br>korumalı | Art                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Lung                      | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Nrv                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | SmP                       | ✓                         | ✓                | ✓                   | ✓                       |            |
|                             | Ven                       | ✓                         | ✓                | ✓                   | ✓                       |            |

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| Dönüştürücü                 | Muayene türü <sup>a</sup> | Görüntüleme modu          |                  |                     |                            |               |
|-----------------------------|---------------------------|---------------------------|------------------|---------------------|----------------------------|---------------|
|                             |                           | 2B <sup>b</sup><br>M Mode | CPD <sup>c</sup> | Renkli <sup>c</sup> | PW<br>Doppler <sup>d</sup> | CW<br>Doppler |
| P10x                        | Abd                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|                             | Crd                       | ✓                         |                  | ✓                   | ✓                          | ✓             |
|                             | Neo                       | ✓                         | ✓                | ✓                   | ✓                          |               |
| P11x <sup>e</sup>           | Art                       | ✓                         | ✓                | ✓                   |                            |               |
|                             | Ven                       | ✓                         | ✓                | ✓                   |                            |               |
| rP19x standart/<br>korumalı | Abd                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|                             | Crd                       | ✓                         |                  | ✓                   | ✓                          | ✓             |
|                             | Lung                      | ✓                         | ✓                | ✓                   | ✓                          |               |
|                             | OB                        | ✓                         | ✓                | ✓                   | ✓                          |               |
|                             | Orb                       | ✓                         | ✓                | ✓                   | ✓                          |               |
|                             | TCD                       | ✓                         | ✓                | ✓                   | ✓                          |               |

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