

NANOMAXX®

TRANSDUCERS



SonoSite




READY FOR A COMPLETE RANGE OF APPLICATIONS

The NanoMaxx® ultrasound system offers a range of transducers with broadband signal processing capabilities, that produce high-resolution diagnostic imaging across a range of clinical applications to support clinical diagnosis and guide interventional procedures. With a splash resistant membrane for easy cleaning and the ability to withstand the impact of a 3-foot drop, you can depend on the NanoMaxx system to perform in the most demanding of situations and environments. The NanoMaxx system and its transducers are covered by SonoSite's industry-leading five-year warranty.



<p>L38n</p>  <p>10-5MHz Linear</p> <p>Applications: CIMT musculoskeletal nerve small parts vascular venous</p> <p>Scan depth: 9 cm</p>	<p>L25n</p>  <p>13-6MHz Linear</p> <p>Applications: CIMT musculoskeletal nerve superficial vascular venous</p> <p>Scan depth: 6 cm</p>	<p>C11n</p>  <p>8-5MHz Curved</p> <p>Applications: abdominal vascular</p> <p>Scan depth: 10 cm</p>	<p>C60n</p>  <p>5-2MHz Curved</p> <p>Applications: abdominal nerve OB GYN</p> <p>Scan depth: 30 cm</p>	<p>P21n</p>  <p>5-1MHz Phased</p> <p>Applications: abdominal cardiology OB</p> <p>Scan depth: 35 cm</p>	<p>ICTn</p>  <p>8-5MHz Curved</p> <p>Applications: OB GYN</p> <p>Scan depth: 13 cm</p>
--	---	---	--	--	---

L52n (Vet)



10-5MHz Linear

Applications:
musculoskeletal
OB
vascular

Scan depth:
15 cm

Needle guide kit is available with the L25n transducer.

For detailed information, please visit www.sonosite.com/products/nanomaxx or speak with your SonoSite customer representative.



5
YEAR

TECHNOLOGY DRIVEN
5-YEAR STANDARD WARRANTY

MADE IN THE USA



SonoSite, Inc.
Worldwide Headquarters
21919 30th Drive SE, Bothell, WA 98021-3904
Tel: +1 (425) 951-1200 or +1 (877) 657-8050
Fax: +1 (425) 951-6800 E-mail: nanomaxx@sonosite.com

www.sonosite.com/products/nanomaxx

SonoSite Worldwide Offices
SonoSite Australasia Pty Ltd – Australia 1300-663-516
– New Zealand 0800-888-204
SonoSite – Brazil +55 61 8118 7100
SonoSite Canada Inc. – Canada +1 888-554-5502
SonoSite GmbH – Germany +49 91-31-97-79-29-0
SonoSite Iberica SL – Spain +34 91-640-49-11
SonoSite India Pvt Ltd – India +91 124-288-1100
SonoSite Italy S.r.l. – Italy +39 02-98490885
SonoSite KK – Japan +81 3-5304-5337
SonoSite Ltd – China – Shanghai +86 21-5239-6693
SonoSite Ltd – United Kingdom +44 1462-444800
SonoSite SARL – France +33 1-69-18-96-30