

Disinfectants for VividTrac™ Device Accessories

This memo summarizes the disinfectant testing that was completed for the following items:

- Nexus 7 Tablet (32G Android™ OS 4.2) with VividVision™ Software (P18235)
- 7-inch USB OTG Adapter Cable (P18390)
- 6-foot USB Extension Cable (P18475)

The tablet, USB extension cable, and OTG cable was tested with Super Sanicloth Wipes, Sanicloth HB Wipes, and Sanicloth Plus wipes. They showed no degradation to any of the materials and the touchscreen on the tablet was still functional after exposure.

For additional questions, please do not hesitate to contact our Service Department directly at 1-877-657-8118. As always, SonoSite is here to serve you.

© SonoSite, Inc. 2014. All rights reserved.

FUJIFILM SonoSite, Inc. the SonoSite logo and other trademarks not owned by third parties are registered and unregistered trademarks of FUJIFILM SonoSite, Inc. in various jurisdictions. All other trademarks are the property of their respective owners. VividTrac device is a trademark of Vivid Medical.