

# Barcode Scanner Expressions Supplement




## Introduction

This supplement provides information on the custom programming feature available with the JADAK Flexpoint HS-1M barcode scanner for use with FUJIFILM SonoSite ultrasound systems.

For information from the manufacturer, see the *JADAK Flexpoint HS-1M Quick Start Guide*.

## Document conventions

The document follows these conventions:

- ▶ A  **WARNING** describes precautions necessary to prevent injury or loss of life.
- ▶ A  **Caution** describes precautions necessary to protect the products.
- ▶ A  **Note** provides supplemental information.
- ▶ Numbered and lettered steps must be performed in a specific order.
- ▶ Bulleted lists present information in list format but do not imply a sequence.
- ▶ Single-step procedures begin with ❖.

For a description of labeling symbols that appear on the product, see "Labeling Symbols" in the system user guide.

## Getting help

For technical support, please contact FUJIFILM SonoSite as follows:

<b>Phone (U.S. or Canada)</b>	877-657-8118
<b>Phone (outside U.S. or Canada)</b>	425-951-1330, or call your local representative
<b>Fax</b>	425-951-6700
<b>Email</b>	ffss-service@fujifilm.com
<b>Web</b>	www.sonosite.com
<b>Europe Service Center</b>	Main: +31 20 751 2020 English support: +44 14 6234 1151 French support: +33 1 8288 0702 German support: +49 69 8088 4030 Italian support: +39 02 9475 3655 Spanish support: +34 91 123 8451
<b>Asia Service Center</b>	+61 2 9938 8700 ffss-apacme-service@fujifilm.com

Printed in the U.S.

## Customizing the barcode scanner

Custom programming is available for the JADAK Flexpoint HS-1M barcode scanner in English only. Contact your customer representative for further information.



**Note** Before using the barcode scanner with patients, verify that your custom program is working as expected.

# Specifications

The ultrasound system with the JADAK Flexpoint HS-1M barcode scanner supports the following scanning formats:

## **1D scanning formats**

Code 39

Code 128

Codabar

UPC-A

UPC-E

EAN-8

EAN-13

Interleaved 2 of 5

Code 93

## **2D scanning formats**

PDF417, except micro PDF417

Data Matrix

QR Code, except micro QR Code

**Manufacturer**

FUJIFILM SonoSite, Inc.  
21919 30th Drive SE  
Bothell, WA 98021 USA  
T: +1-888-482-9449 or +1-425-951-1200  
F: 1-425-951-1201

**EC Authorized Representative**

FUJIFILM SonoSite B.V.  
Joop Geesinkweg 140  
1114 AB Amsterdam,  
The Netherlands

**Australia Sponsor**

FUJIFILM SonoSite Australasia Pty Ltd  
114 Old Pittwater Road  
BROOKVALE, NSW, 2100  
Australia



**Caution** United States federal law restricts this device to sale by or on the order of a physician.

SONOSITE is a registered trademark of FUJIFILM SonoSite, Inc. in various jurisdictions. FUJIFILM is a registered trademark of FUJIFILM Corporation in various jurisdictions. Value from Innovation is a trademark of FUJIFILM Holdings America Corporation.

All other trademarks are the property of their respective owners.

Part number: P22992-02

Publication date: August 2020

Copyright © 2020 FUJIFILM SonoSite, Inc. All rights reserved.



P22992-02

